
baboo®

Duo Sense infrared forehead and
ear thermometer
Operating Manual

Model AOJ-20A

LANGUAGES

ENGLISH	3	LIETUVIŠKAI	59	ქართული	115
ESPAÑOL	18	LATVIEŠU	73	ROMÂNĂ	133
FRANÇAIS	32	EESTI	87	AZƏRBAYCAN	147
SUOMI	46	УКРАЇНСЬКА	100	РУССКИЙ	161
				العربية	175

ENGLISH

INTRODUCTION

Thank you for purchasing this Baboo® Forehead & Ear Infrared Thermometer. Please read the Operating Manual carefully before using the thermometer and put it in a safe and secure place for reference.

FOREWORD

The infrared thermometer can be used to read the body temperature by measuring the ear and forehead temperature, which is suitable for adult and children (ear test mode only for the child above 3 months).

PACKAGE CONTENTS

No.	Name	Quantity
1	Infrared Thermometer	1 pcs.
2	Pouch	1 pcs.
3	Operating Manual	1 pcs.
4	AAA Battery (Optional)	2 pcs.

1. WARNINGS AND PRECAUTIONS	4
2. Product Description	5
2.1 Overview	5
2.2 Structure	5
2.3 Operating principle	5
2.4 Indications for use	6
2.5 Contraindications	6
3. Features	6
4. Product Structure	6
5. Display description	7
6. How to use your thermometer	7
6.1 In order to avoid the inaccuracy	7
6.2 Take your forehead temperature	7
6.3 Take your ear temperature	8
6.4 After a measurement	9
6.5 Read your temperature	9
6.6 Switching between mute and un-mute	9
6.7 Checking 40 sets of memory data	9
6.8 °C/°F conversion	10
6.9 To turn off	10
6.10 Replace the battery	10
7. Temperature taking tips	10
8. Care and cleaning	11
9. Error and Troubleshooting	11
10. Technical Specifications	13
11. Symbols	13
12. Appendix: EMC information-guidance and manufacturer's declarations	14
13. Warranty and after-sale service	17

1. WARNINGS AND PRECAUTIONS

1. Keep out of reach of children under 12 years.
2. Never immerse the thermometer into water or other liquids (it is not waterproof). For cleaning and disinfecting, please follow the instructions in the "Care and cleaning" section.
3. Never use the thermometer for purposes other than those it has been intended for. Please follow the general safety precautions when using on children.
4. Keep the thermometer away from direct exposure to the sun and keep it in a dust- free, dry area, well-ventilated place at a temperature between 5°C (41°F) - 40°C (104°F). Do not use the thermometer in high-humidity environments. (>95% RH). Do not use the thermometer if there are signs of damage.
5. Damage on the measuring sensor or on the instrument itself. If damaged, do not attempt to repair the damage instrument! Please contact the dealer.
6. This thermometer contains of high-quality precision parts. Do not drop the instrument. Protect it from severe impact and shock. Do not twist the instrument or the measuring sensor.
7. Please consult your doctor if you see symptoms such as unexplained irritability, vomiting, diarrhea, dehydration, changes in appetite or activity, seizure, muscle pain, shivering, stiff neck, pain when urinating, etc., even in the absence of fever.
8. Even in the absence of fever, those who exhibit a normal temperature may still need to receive medical attention. People who are on antibiotics, analgesics, or antipyretics should not be assessed solely on temperature readings to determine the severity of their illness.
9. Temperature elevation may signal a serious illness, especially in adults who are old, frail, have a weakened immune system, or neonates and infants. Please seek professional advice immediately when there is a temperature elevation, especially if you are taking the temperature of someone who is:
 - Over 60 years of age (Fever may be blunted or even absent in elderly patients)
 - Having diabetes mellitus or a weakened immune system (e.g., HIV positive, cancer, chemotherapy, chronic steroid treatment, splenectomy)
 - Bedridden (e.g., nursing home patient, stroke, chronic illness)
 - A transplant patient (e.g., liver, heart, lung, kidney)
10. This thermometer is not intended for preterm babies or small-for- gestational age babies. This thermometer is not intended to interpret hypothermic temperatures. Do not allow children to take

their temperatures unattended.

11. Use of this thermometer is not intended as a substitute for consultation with your physician or pediatrician.
12. Clean the thermometer probe after each use.
13. Do not use the thermometer on newborns or for continuous temperature monitoring purposes.
14. Do not take a measurement while or immediately after nursing a baby.
15. Patients should not drink, eat, or be physically active before or during the measurement.

2. PRODUCT DESCRIPTION

1) Overview

Infrared thermometer measures the body temperature based on the infrared energy emitted from the eardrum or the forehead. Users can quickly get measurement results after positioning properly the temperature probe in the ear canal or forehead.

Normal body temperature is a range. The following table shows that this normal range also varies by site. Therefore, readings from different site should not be directly compared. Tell your doctor what type of thermometer you used to take your temperature and on what part of the body. Also bear this in mind if you are diagnosing yourself.

Measuring part	Normal temperature range
Forehead temperature	36.1°C - 37.5°C (97°F - 99.5°F)
Ear temperature	35.8°C - 37.5°C (96.44°F - 99.5°F)
Oral temperature	35.5°C - 37.5°C (95.9°F - 99.5°F)
Rectal temperature	36.6°C - 38°C (97.9°F - 100.4°F)
Axillary temperature	34.7°C - 37.3°C (94.5°F - 99.1°F)

2) Structure

The device consists of an enclosure, a button, a temperature sensor, Infrared temperature-measuring element, microcomputer controlled circuit, LCD, backlight and a buzzer.

3) Operating principle

The infrared thermometer's working principle is based on the fourth-power law: infer an object's radiation temperature by measuring its infrared energy. Thermopile sensor can convert the infrared energy into thermoelectricity, and output as a detected signal after signal processing.

4) Indications for use

The infrared thermometer takes human body temperature via the eardrum or forehead. It applies to all age groups except for babies under three months. Both devices apply to both professional use and home use.

5) Contraindications


Do not use the thermometer if the ear is infected with otitis or suppuration.

3. FEATURES

- Magnetic probe cover and automatic mode switch between forehead & ear mode (OPTIONAL)
- Quick measurement, less than 1 second
- Small body with super large font
- Multi-functional, can measure the ear and forehead
- Easy operation, one button design, to measure ear, forehead
- 40 sets of memories, easy to recall
- Switching between mute and unmute mode
- Audio and visual fever warning
- Switching between °C and °F
- Auto shut-down and power-saving
- Backlight color temperature indication
- Patent technology to realize separate measurement of adults and children


4. PRODUCT STRUCTURE

1. LCD display screen
2. Mode/Memory button
3. Measure/Mute-unmute button
4. Probe
5. Battery cover
6. Probe cover (remove it when measuring ear temperature)

5. DISPLAY DESCRIPTION

1. Adult forehead mode
2. Child forehead mode (under 12 years old)
3. Ear mode
4. Adjust mode
5. Fahrenheit/Celsius degrees
6. Previous temperature value
7. Low battery indication
8. Mute/Unmute icon
9. Current temperature value

6. HOW TO USE YOUR THERMETER

1) In order to avoid inaccuracies:

1. Please ensure that the device is only used in the room and that there is no strong wind or loud noises.
2. Please make sure that there is no intense emotion and movement before measuring.
3. If the device is transferred from one condition to another, which has different ambient temperature, it is suggested to rest for more than 30 minutes.
4. If the tester is transferred from one condition to another, which has different ambient temperature, it is suggested to have a rest for more than 10 minutes.
5. Do not hold the device for a long time, as it is highly sensitive to heat. The device has undergone a clinical test; it is safe and accurate when using in accordance with the operation manual.

2) Take your forehead temperature

When using the thermometer for the first time, please load the batteries.

Press and release the **Measure/Mute-unmute** button for one second to power on. Press the **Mode/Memory** button to choose  or . Then point the temperature probe at the center of the forehead, with a distance of 0-3 cm (0-1.18 inches).

Press and release the **Measure/Mute-unmute** button for 1 second. The beep is heard; you can now read the value.

NOTE: The forehead measurement is an indication. The measured forehead temperature can fluctuate up to 1° F/0.5°C from your actual body temperature. Please be aware of the factors that influence the accuracy as described in the section “Temperature taking tips” and “WARNINGS AND PRECAUTIONS” .

- ⚠ If the eyebrow area is covered with hair, sweat or dirt, please clean the area beforehand to improve the reading accuracy.
- ⚠ Always check if the lens is clean.
- ⚠ Always make sure the user and the thermometer will have been in the same room for at least 30 minutes prior to the measurement.

3) Take your ear temperature

1. Press the **Measure/Mute-unmute** button to power on.
2. Press the **Mode/Memory** button to choose  and take off the probe cover (if with magnetic probe cover, it is automatically switched between forehead & ear mode), fit the probe snugly into the ear canal.
3. Press and release the **Measure/Mute-un-mute** button for 1 second, the beep is heard, you can now read the value.

- ⚠ Please make sure the ear is clean and free of earwax or obstructions.
- ⚠ The right ear reading may differ from the reading taken at the left ear. Therefore, always take the temperature in the same ear.

NOTE: Children under 1 year: Pull the ear straight back. For children aged 1 year and older and for adults: Pull the ear up and back.

- ⚠ Do not force the thermometer into the ear canal.
- ⚠ Otherwise, the ear canal may get injured.
- ⚠ When taking the temperature on an adult, gently pull the ear up and back to make sure the ear canal is straight, so that the temperature probe can receive an infrared ray from the eardrum.
- ⚠ Be careful when taking the temperature of a child, whose ear canal is small.

4) After a measurement

1. Once the reading has been completed, remove the thermometer from the forehead/ear and observe the temperature.
2. After each measurement, you can enter the recall mode and query earlier temperature readings.

- ⚠ Do not hold the thermometer for a long time, because it is sensitive to the ambient temperature.
- ⚠ After each measurement, clean the temperature probe with a soft cloth, and put the thermometer in a dry and well-ventilated place.
- ⚠ You should wait at least 10 seconds between each measurement.
- ⚠ It is dangerous to make a self-diagnosis or self-treatment based on the obtained measurement results.

For such purposes, please consult a doctor.

5) Reading your temperature: T indicates a temperature reading.

1. In forehead or ear mode.
 - If $32^{\circ}\text{C} \leq T \leq 37.5^{\circ}\text{C}$ ($89.6^{\circ}\text{F} \leq T \leq 99.5^{\circ}\text{F}$), the green light will be displayed, with one long beep.
 - If $37.6^{\circ}\text{C} \leq T \leq 38.0^{\circ}\text{C}$ ($99.7^{\circ}\text{F} \leq T \leq 100.4^{\circ}\text{F}$), the orange light will be displayed with 6 short beeps, indicating that you may have a mild fever.
 - If $38.1^{\circ}\text{C} \leq T \leq 42.9^{\circ}\text{C}$ ($100.6^{\circ}\text{F} \leq T \leq 109.2^{\circ}\text{F}$), the red light will be displayed with 6 short beeps, indicating that you may have a high fever.
2. In Adjust mode.
 - If $0^{\circ}\text{C} \leq T \leq 100^{\circ}\text{C}$ ($32.0^{\circ}\text{F} \leq T \leq 199.0^{\circ}\text{F}$), the white light will be displayed with one long beep. This mode is not for human temperature measurement, for maintenance use only.

6) Switching between mute and unmute.

When the thermometer is turned on, keep pressing the **Measure/Mute-unmute** for 5 seconds, to switch from unmute to mute.

7) Checking 40 sets of memory data.

When the thermometer is turned on or off, keep pressing the **Mode/Memory** button for 5 seconds to go to the memory mode, press this button again to check the 40 sets of memory one by one. If no value, it will display “---”.

8) °C/°F conversion.

Open the battery cover, use the toggle switch to change the °C/°F.

9) Automatically turn off.

The unit will shut down automatically after 13 seconds of inactivity.


CAUTION:

⚠ All memory records will be lost when uninstalling or reinstalling the battery.

10) Replace the battery.

Put two AAA batteries correctly into the compartment.

⚠ Remove the batteries if the thermometer will not be used for more than two months.

7. TEMPERATURE TAKING TIPS

- 1) It is important to know everyone's normal temperature when they are well. This is the only way to accurately diagnose a fever. Record readings twice a day (early morning and late afternoon). Take the average of the two temperatures to calculate normal oral equivalent temperature. Always take the temperature in the same location, since the temperature readings may vary across different areas of the forehead.
- 2) A child's normal temperature can be as high as 99.9°F (37.7°C) or as low as 97.0°F (36.1°C). Please note that this unit reads 0.5°C (0.9°F) lower than a rectal digital thermometer.
- 3) External factors may influence ear temperatures, including when an individual has:
 - been lying on one ear or the other
 - had their ears covered
 - been exposed to very hot or very cold temperatures
 - been recently swimming or bathing

In these cases, remove the individual from the situation and wait 20 minutes before taking their temperature.

Use the untreated ear if prescription ear drops or other ear medications have been placed in the ear canal.

- 4) Holding the thermometer for too long in the hand before taking a measurement can cause the device to warm up. This means the measurement could be incorrect.
- 5) Patients and the thermometer should stay in steady-state room condition for at least 30 minutes.
- 6) Before placing the thermometer sensor onto the forehead, remove dirt, hair, or sweat from the forehead area. Wait 10 minutes after cleaning before taking a measurement.
- 7) Use an alcohol swab to carefully clean the sensor and wait for 5 minutes before taking a measurement on another patient. Wiping the forehead with a warm or cool cloth may change your reading. It is advised to wait 10 minutes before taking a reading.
- 8) In the following situations, it is recommended to take 3-5 measurements in the same location and use the highest as the reading:
 - Newborn infants in the first 100 days.
 - Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
 - When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the instrument and obtains consistent readings.

8. CARE AND CLEANING

Use an alcohol swab or cotton swab moistened with 70% alcohol to clean the thermometer casing and the measuring probe. After the alcohol has completely dried out, you can take a new measurement.

Ensure that no liquid enters the interior of the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the instrument in water or other cleaning liquids. Take care not to scratch the surface of the LCD screen.

9. ERROR AND TROUBLESHOOTING

Symptom	Possible Cause	Description & Solution
Failed to power on	The battery level is too low.	Replace with a new battery.
	The polarities of the batteries are reversed.	Ensure the batteries are in the right position.
	The polarities thermometer is damaged.	Contact the dealer.

The reading is too low	The lens of the probe is dirty.	Clean the lens with a cotton swab.
	The distance between the item and the target is too far.	Keep in contact with the forehead or put the probe into the ear canal.
	You have just come from a cold environment.	Stay in a warmer room for at least 30 minutes before taking a reading.
The reading is too high	You have just come from a hot environment.	Stay in an adequately cool room for at least 30 minutes before taking a reading.
	The ambient temperature is out of range.	One long beep and red backlight for 3 seconds. Take a measurement when the ambient temperature is between 5°C (41.0°F) and 40°C (104°F).
	Memory error or calibration is not finished. The unit will turn off automatically after 3 seconds.	One long beep; the screen will be displayed in red. Contact the dealer.

Symptom	Possible Cause	Description & Solution
	In ear or forehead mode, $T > 42.9^{\circ}\text{C}$ (109.2°F) In adjust mode, $T > 100^{\circ}\text{C}$ (199°F)	One long beep, the screen will be displayed in red.
	In ear or forehead mode, $T < 32^{\circ}\text{C}$ (89.6°F) In adjust mode, $T < 0^{\circ}\text{C}$ (32°F)	One long beep, the screen will be displayed in red.
	The power voltage is between 2.4V-2.7V	The battery icon will flicker, which means the battery is low, but you can still use the thermometer.
	The power voltage is lower than 2.4V	The battery icon will flicker, and the unit will turn off automatically after 3 seconds. Please replace with a new battery.

10. TECHNICAL SPECIFICATIONS

Product name	Forehead & Ear infrared thermometer	
Applicable regulations and laws	ASTM E 1965-98	
	ISO80601-2-56	
Display	Segment LCD, four color LED backlight (white, green, orange, red)	
Temperature units	°C / °F, switchable	
Power supply	2x1.5V AAA Batteries	
Measurement range	Measurement: 32.0°C–42.9°C (89.6°F–109.2°F)	
	Direct mode: 0°C–100°C (32°F–199°F)	
Accuracy (Laboratory)	Measurement mode	±0.2°C/±0.4°F
	Direct mode	±1.0°C /±2.0°F
Display resolution	0.1°C / ° F	
Automatic	13s±1s	
Shutdown memory	40 groups of measured temperature	
Operational conditions	Temperature: 5°C-40°C(41°F-104°F)	
	Humidity: 15-95%RH, non-condensing Atmospheric pressure: 70-106kPa	
Storage condition	Ambient temperature: -20°C-55°C(-4°F-131°F)	
	Relative humidity: 0-95% RH, non-condensing Atmospheric pressure: 50kPa to 106kPa	
Battery	2*AAA, can be used for more than 3000 times	
Weight & Dimensions	60g (without battery), 143x35x41mm	

11. SYMBOLS

Symbol	Description
	Type BF applied part.
	Information about the manufacturer.
	Please read the instructions carefully.
	Waste electrical equipment should be sent to a dedicated collection point for recycling.
	Production Date
LOT	Batch number
	IMPORTANT Inaccurate readings or thermometer damage may occur if the thermometer is not used correctly.

IP22	2 Protected against solid foreign objects of 12,5 mm Ø and greater. 2. If the thermometer is kept at a 15-degree angle, it can still prevent water drops.
 0123	CE marking
	Authorized Representative in the European Community
	Keep away from the sunlight
	Fragile, handle with care
	Keep dry
	This way up

12. EMC INFORMATION-GUIDANCE AND MANUFACTURER'S DECLARATIONS

CAUTION:

- The infrared thermometer needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided for in the ACCOMPANYING DOCUMENTS.
- Portable and mobile RF communications equipment can affect the infrared thermometer.
- The infrared thermometer should not be used adjacent to or stacked with other equipment. Guidance and manufacturer's declaration – electromagnetic emission – for all EQUIPMENT AND SYSTEMS.

ELECTROMAGNETIC EMISSIONS

Guidance and manufacturer's declaration - Electromagnetic emission	
The infrared thermometer is intended for use in the electromagnetic environment specified below. The customer or user of the infrared thermometer should ensure that it is used in such an environment.	
Emission test	Compliance
RF emissions CIS PR 11	Group 1

RF emissions CIS PR 11	Class B
Harmonic emissions IEC61000-3-2	Not applicable
Voltage fluctuations/flicker emissions IEC61000-3-3	Not applicable

ELECTROMAGNETIC IMMUNITY

Guidance and manufacturer's declaration - Electromagnetic Immunity		
The infrared thermometer is intended for use in the electromagnetic environment specified below. The customer or user of the infrared thermometer should ensure that it is used in such an environment.		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kv contact ±2 kv, ±4 kv, ±8kv, ±15kv air	±8 kv contact ±2 kv, ±4 kv, ±8kv, ±15kv air
Electrical fast transient/burst IEC 61000-4-4	Not applicable	Not applicable
Surge IEC61000-4-5	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable
Power frequency magnetic field IEC 61000-4-8	30A/m 50HZ/60HZ	30A/m 50HZ/60HZ
Conducted RF IEC61000-4-6	Not applicable	Not applicable
Radiated RF IEC61000-4-3	10v/m 80 MHZ- 2,7 GHZ 80%AM at 1 KHZ	10v/m 80 MHZ- 2,7GHZ 80%AM at 1 KHZ
NOTE: UT is the AC mains voltage prior to application of the test level.		

ELECTRO MAGNETIC IMMUNITY

Guidance and manufacturer's declaration - Electromagnetic Immunity

The infrared thermometer is intended for use in the electromagnetic environment specified below. The customer or user of the infrared thermometer should ensure that it is used in such an environment.

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (w)	Distance (m)	Test Level (w/m)	Compliance level (w/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY for RF wireless communications equipment)	385	380-390	TETRA 400	Pulse modulation 18 HZ	1.8	0,3	27	27
	400	430-470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1kHz sine	2	0,3	28	28
	710	704-787	LTE Band 13,17	Pulse modulation 217 HZ	0.2	0,3	9	9
	745							
	780							
	810	800-960	GSM 800-900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 HZ	2	0,3	28	28
	870							
	930							
	810	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 HZ	2	0,3	28	28
	870							
	930							
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450,	Pulse modulation 217 HZ	2	0,3	28	28
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450,	Pulse modulation 217 HZ	2	0,3	28	28
			LTE Band 7					
	810	5100-5800	WLAN 802.11a/n	Pulse modulation 217 HZ	0,2	0,3	9	9
	870							
	930							

13. WARRANTY AND AFTER-SALE SERVICE

This product is covered by the Global Standard Warranty, which is valid from the date of purchase. The warranty period is two years in European countries and one year in all other regions worldwide. The batteries, the packaging, and any damage caused by improper use are not covered by the warranty.

The warranty excludes the following user-caused failures: failures resulting from unauthorized disassembly or modification, from an unexpected drop during application or transportation, and from not following the instructions in the operating manual.

Discover more Baboo® products at www.baboo.baby. Baboo® is a registered trademark of Baby's Kingdom LTD, 10 John Street, London, WC1N 2EB, UK. Made in China.

 Importer in EU: UAB „Sentija“, Vilnius, Tuskulenu str. 13A-22, LT-09217, Lithuania, baboo@sentija.com.

 Shenzhen AOJ Medical Technology Co., Ltd.
Add: Room 301&4F, Block A, Building A, Jingfa Intelligent Manufacturing Park, Xiawei Yuan, Gushu Community, Xixiang Street, Bao'an District, 518126 Shenzhen, CHINA

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Discover more Baboo® products at www.baboo.baby. Baboo® is a registered trademark of Baby's Kingdom LTD, 10 John Street, London, WC1N 2EB, UK. Made in China. Importer in EU: UAB „Sentija”, Vilnius, Tuskulenu str. 13A-22, LT-09217, Lithuania. baboo@sentija.com.