NebuFlow

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Firma: Dongguan Aidsy Machinery & Electronic Equipment Co.,Ltd
Adres: Part B, 3rd F, Block A, Wentang Industrial Park, Longhua Road, Zhouwu
Dongcheng District Dongguan City, Guangdong Province, 523118 P.R. China

2024/05/01

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01. INTRODUCTION

Thank you for choosing our Mesh Nebulizers. Our Mesh Nebulizer is a handheld device, designed to aerosolize medication for respiratory therapy purpose. Considering about the difference of your needs between traveling and indoor use.

This device can operate on internal lithium titanate rechargeable battery. With battery and compact size, you are able to continue your treatment while you go on traveling. To help save the planet.

As this device is a medical instrument. Please read instruction. Manual before use, and be sure to follow the instructions of a doctor and use the device correctly.

Product Features

1. Microporous (Mesh) and oscillator combines the technique of atomization
2. USB cable charging. Easy to obtain re-charge process
3. Compact body design, handy to carry around
4. Palm size device easy to carry
5. Easy to operate
6. Easy to clean & maintain

Medical purpose	This device is designed to nebulize the solution into an aerosol for respiratory therapy purpose.	
Intended user	Legally certified medical experts, such as doctors, nurses, therapists, or patients under the guidance of qualified medical experts. The user should also be capable of understanding general operation of this device and the content of instruction manual.	
Intended patients	This device design for all ages of patients Except for the patients who are unconscious, not able to breathing spontaneously or having pulmonary edema.	
Recommended Operation Environment	This device is intended for using in the medical facility, such as hospital, clinic and doctor's office, a room of general household and open-air environment with a roof. • Temperature range: 10 °C to 40 °C • Humidity: 30 to 85 % RH	
Durable period	Durable periods as follows. The product is used to nebulize 0.9% Saline 2 times a day, 10 minutes each time at room temperature (23°C). Durable period may vary depending on environment and solution.	
	Main unit:	36 months
	Medication cup:	12 months
Precautions for use	Warnings and cautions described in the manual should be observed.	

02. SAFETY PRECAUTIONS

⚠ READ ALL WARNINGS AND INSTRUCTIONS BEFORE YOU USE THIS DEVICE.

As with any medical device, this product may become unusable due to an electrical outage, battery depletion, or mechanical impaction. We recommend having a backup device available to you.

When using electrical products, always follow basic safety precautions. Take extra care around children.

⚠ WARNINGS

- Only use this device for medications prescribed by your doctor.
- The nebulizer is only intended for respiratory therapy and any other application of this device is improper and dangerous. The manufacturer cannot be held liable for any damage caused by improper or incorrect use.
- Do not share your nebulizer with others. It is intended to be used by a single user. If more than one person uses it, there is a risk of spreading infectious illnesses.
- Clean all parts of your nebulizer before use, after each use, and after extended storage.
- Adult supervision is required when this device is used by children and individuals who require special assistance.
- Never operate this device if any of the parts are not working properly or have been damaged.
- Be sure the device has been properly cleaned before use to avoid possible contamination.
- Do not plug in or unplug the adapter from the electrical outlet with wet hands.
- Do not drop medication on the main unit or USB port. If you drop medication on either area, immediately wipe it off with gauze. If you use the device while it is still wet, it may cause trouble or injury.
- Do not leave the device in a car if it be subjected to significant heat or cold. The battery is not replaceable. Do not attempt to replace the battery yourself, as you may damage it, cause overheating, and injury. Could result, NEVER incinerate the battery, NEVER immerse the


battery in water, as it may cause the battery to rupture or injury.

- All battery deteriorate over time if they are not useful or charged, do not store the device for long period of time without charging it periodically. If the device is not useful or charged for long period of time, the battery may create a hazardous condition.
- Never use the device while it is charging.
- If you are not using the device for a long period of time, disconnect the USB charger.
- The device contains sensitive components, including a stainless steel disk with carefully measured 3–5 micron holes. Do not drop, crush, puncture, bend, heat, incinerate, or apply strong stock to the device or its parts.
- Do not wash the medication cup and mesh under running water or any other liquid.

CAUTIONS

- Do not attempt to clean mesh with any foreign objects; it may damage the mesh.
- Avoid dropping the device, otherwise, it may not function normally.
- Keep the device away from direct sunlight, excessive heat, or cold to avoid damaging the batteries.
- Do not attempt to open, repair, or modify this device.
- Do not drop any liquids on the main unit or on the adapter. If you drop liquids on them, wipe it off immediately.
- Follow local laws and recycling plans regarding disposal or recycling of components, batteries, and packaging.
- If the USB wall charger and USB port do not join with reasonable ease, they probably do not match. Make sure that you are using the USB wall charger provided and that you have positioned it correctly.

03. CLASSIFICATION AND EXPLANATION OF SYMBOLS

	Warning / Caution / Note
	Class II equipment per IEC 60601-1
	Type BF equipment per IEC 60601-1
IP22	IP22: Protected against foreign objects equal to or greater than 12.5mm in diameter and against drops of water falling at up to 15° from vertical.
	Consult instructions for use
	ON / OFF Button
	Keep dry
	Operating temperature limits: 10°C to 40°C
	Storage and delivery temperature limits: -20°C to 60°C
	Operating humidity limits: 30% to 85% R,H
	Storage and delivery humidity limits: 20% to 75% R,H
	Manufacturer
	Date of manufacture
	Serial number

LOT	Batch code
	Disposal of Electrical & Electronic Equipment (WEEE): Do not treat this product as household waste.
CE 1639	This device complies with the requirements of the Medical Devices Directive (93/42/EEC)
EC REP	Authorized representative in the European Community

04. PACKAGE CONTENTS

Check before use


The following items are contained in the package. Please check all parts for visible damage. Replace any damaged parts before you use this device. In the case of missing parts, malfunction, or damage, please contact the store where you purchased it or the nearest dealer.

Main unit

Medication cup

Instruction Manual

Mask


Mouth piece

USB cable




Nosepiece
(Optional)


USB charger
(Optional)

05. SYSTEM OVERVIEW

Press Release Button to remove the medication cup

Bottom view

06. HOW TO USE THE BATTERY


Battery information:

The device has a built-in, replaceable lithium titanate battery that cannot be replaced. NEVER use the device while it is charging. The device is charged via USB, like many mobile phones and portable electronic devices. The battery can only be charged when connected to a wall outlet using the supplied USB charger.

NOTE: The USB charger is an optional accessory. Only use an external charger that complies with IEC60601-1:2012, with an output voltage of 5,0V 1-2A.

How to charge the battery from a wall outlet:

1. Open the USB port cover to expose the USB port.
2. Gently insert the USB cable into the USB port.
3. Gently connect the USB cable to the USB wall charger port.
4. Plug the USB charger into the wall outlet.

Battery Status:

When the device is not connected to the power source, the ON-OFF indicator can indicate three states:

- no light: off
- blue light: in use
- orange light: low battery

When the device is connected to the power source, the ON-OFF indicator can indicate two states:

- flashing orange light: battery is charging
- flashing blue light: battery is fully charged

NOTE: Batteries do not have an unlimited life, Battery life and number of charging cycles will vary depending on use.

⚠ CLEAN ALL PARTS OF THE PORTABLE NEBULIZER BEFORE USE, AFTER EACH USE, AND AFTER LONG-TERM STORAGE.

1. Preparing the medication cup

Step 1: Remove the medication cup and the main unit.

Step 2: Connect the medication cup to the main unit.

Step 3: Make sure the medication cup is properly connected to the main unit.

2. Filling with medication

Make sure the power indicator is off before adding medication.

- Open the medicine container lid.
- Fill the medicine container WITHOUT EXCEEDING THE LIMIT (max. 6 ml).
- Close the medicine container lid properly.

3. Install the mask or mouthpiece into the mask adapter

The device is now ready for use. Refer to the next section on how to inhale.

⚠ Caution: For type, dose, and regime of medication, follow the instructions of your doctor or pharmacist.

1. Slightly tilt as illustrated in the figure.

In this position, the vibrating element is immersed in the medication, and nebulization will start after power-on. After the vibrating element has been immersed in the medication, the unit can be used under any angle.

Note: In some positions (e.g., upright), nebulization may stop after a short while. In that case, briefly tilt the unit again to re-immers the vibrating element in the medication.

2. Place the mouthpiece into your mouth or place the mask to cover your mouth and nose.

3. If you push the ON/OFF () button, nebulization starts.

If you push the ON/OFF () button again, nebulization speed is lower.

Press the ON/OFF button again, and it will turn off.

Long press the ON/OFF button for 3 seconds to start the cleaning mode.

The blue power indicator lights up during nebulization,

Always be calm and relaxed when you inhale. Breathe in slowly and deeply so that the medication can reach deep into the lungs. Hold your breath briefly, then breathe out slowly while removing the mouthpiece from your mouth. Do not breathe too quickly. Pause when you feel like having a rest.

4. After inhalation, always switch off the nebulizer.

5. Clean the unit after each inhalation.

⚠ 1. Tilt the device slightly when the solution is almost exhausted to ensure that the residual solution contacts the mesh and is sprayed completely.

2. If a high-viscosity solution is used, nebulization may be reduced.

3. If excessive solution is accumulated on the mesh, nebulization may stop. In this case, turn off the power and absorb the solution with gauze or a lint-free towel.

4. The device will auto power-off after 10 minutes.

2. Do not drop medication on the main unit or into its USB port. If you use the device when it is still wet, it may cause trouble or injury.

07. CLEANING AND DISINFECTION

⚠ If the device is not cleaned and disinfected correctly and frequently as indicated, microorganisms may remain in the unit and cause risk of infection.

Caution: **NEVER** use the device when it is charging.

NEVER pour liquid into it. Liquid may **DAMAGE** the main unit.

If not cleaned immediately after use, the atomized mesh will be blocked, which may lead to the absence of fog in the next use.

1. After every use:

Clean the medication cup with mineral water after each inhalation.

1. Pour out the residual solutions in the medication cup.
2. Pour some mineral water into the medication cup.
3. Turn on the device to nebulize the mineral water for 2 to 3 minutes to clean the mesh.
4. During use, keep the main unit in a dry state.
5. Remove the medication cup from the main unit.
6. Wash and rinse the medication cup with mineral water.
7. Shake off excess water and allow parts to fully air dry in a clean place.
8. Use gauze or a clean towel to wipe off stains on the main unit if necessary.
9. Make sure that all cleaned parts are completely dry before storing or using them next time.

2. If the atomized mesh is blocked, please refer to the following operations:

Clean the medication cup with approximately 60% edible white vinegar after using it 9 times.

1. Pour out the residual medication in the medication cup.
2. Fill some mineral water into the medication cup.
3. Pour out the mineral water from the medication cup.
4. Fill 3ml of approximately 60% edible white vinegar into the medication cup.
5. Turn on the device to nebulize the vinegar for 5 to 10 minutes to clean the mesh.
6. Disconnect the adapter and remove the medication cup from the main unit.
7. Wash and rinse the medication cup with mineral water.
8. Shake off excess water and allow parts to fully air dry in a clean place.
9. Use gauze or a clean towel to wipe off the main unit if necessary.
10. Make sure that the medication cup and other cleaned parts are completely dry before storing or using them next time.

- ⚠**
1. Keep the battery compartment dry at all times.
 2. Do NOT poke the mesh with fingers, cotton swabs, or any objects.
 3. Do NOT clean parts in a dishwasher.
 4. Do NOT use a microwave to dry any parts.
 5. Do NOT wash the medication cup and mesh under running water or any other liquid.

Daily disinfections:

It is important to disinfect the medication cup on a daily basis.

1. Disinfection by boiled water:

- Rinse the medication cup with distilled water.
- Bring a saucepan of distilled water to the boil.
- **NOTE: Do NOT boil the medication cup directly.**
- Carefully immerse the medication cup in the boiled water for a maximum of 15 minutes.
- Carefully remove the medication cup from the boiled water and shake off excess water.
- Allow parts to cool down and fully air dry on a clean, dry towel, out of reach of children.
- Make sure all cleaned parts are completely dry before storing or using them next time.

2. Disinfection by alcohol:

- Rinse the medication cup with distilled water.
- Immerse the medication cup in 75% ethyl alcohol for 1 minute.
- Rinse the medication cup with distilled water again, shake off excess water, and allow parts to fully air dry on a clean, dry towel.
- Make sure all cleaned parts are completely dry before storing or using them next time.

3. Cleaning

Mouthpiece & Mask Cleaning

Cleaning & Disinfection for the Nebulizer's Accessories (including: nebulizer cup, mouthpiece, child mask).

After Use the Nebulizer Accessories:

- Clean the Nebulizer's Accessories with distilled water.
- Completely immerse the Nebulizer's Accessories into approx. 75% medicinal alcohol.
- Soak for approx. 30 minutes.
- Wash and rinse the Nebulizer's Accessories with about 50°C warm distilled water.
- Shake out water, dry them with a clean, soft towel, or leave them to air-dry in a clean place.
- Make sure that the Nebulizer's Accessories and other cleaned parts are completely dry, then put all these parts in a dry, sealed container.

- ⚠ 1. Alcohol is highly flammable. Do NOT use alcohol in close vicinity of open fire or smoke.
The following item is able to be disinfected by boiled water or alcohol

The below parts are NOT able to be disinfected by boiled water or alcohol

Main unit

USB cable

- ⚠ Do NOT rinse or immerse the main unit in any liquid.

08. CARRYING AND STORING

Close the mouthpiece's cover, store the device and medication cup in a dry and clean environment.

- ⚠ 1. Do NOT leave or carry the device containing residual liquid in the medication cup.
2. Do NOT leave the device under direct sunlight, in high humidity, extreme heat, or cold environment.
3. Keep this device away from fire, high electromagnetic fields, and out of the reach of children.

09. TROUBLESHOOTING

If any trouble occurs while you are using the device, please check the following list first.

Problem	Possible Cause	Action
Low atomization	Batteries are low (orange light)	Please charge the battery
	The stains on the electrode cause a fault connection	Use rubbing alcohol to clean electrodes
	The mesh holes have been clogged or stained	Refer to Cleaning and Disinfection procedure to clean the medication cup page
	The mesh is broken	Replace the medication cup
After turning on the power, power indicator is not on and no mist comes out	Batteries are too low	Please charge the battery
	Fault connection between adapter and main unit	Check and reconnect the USB charger to the main unit
After turning on the power, the power indicator is on but no atomization	The medication cup is not installed properly	Refer to Set up the medication cup procedure to re-install the medication cup
	The mesh holes have been clogged or stained	Refer to Cleaning and Disinfection procedure to clean the medication cup page
	The mesh is broken	Replace the medication cup
Power indicator turns into orange light	Low battery power	Please charge the battery

- ⚠ 1. If the device does not nebulize normally after taking the above-mentioned procedure, please contact the store where you purchased the device or the nearest dealer.
2. Never operate this device if any of the parts are not working properly or have been damaged.


10. SPECIFICATION

Model	MESH-S600A
Product	Portable Ultrasonic Nebulizer
Method of operation	Active Ultrasonic Vibrating Mesh Technology
Power supply	Internal power supply, Lithium battery only
Power consumption	Approx. 2.0W
Mode of operation	10 mins Auto off (with rechargeable battery)
	Press the button once to turn on the device and enter the normal gear
	Press the button again to enter the low gear
	Press the button again to turn off the device
Vibrating frequency	The nebulization rate of low gear is about 80% lower than the normal gear
	Long press for 3 seconds to cleaning mode, after that it'll turn off automatically
Capacity of medication cup	Approx. 107 kHz \pm 10%
Nebulization rate	Max. 6ml
Particle size	\geq 0.25 ml/min
MMAD	Approx. (74.53 \pm 3.01)% within 0.5 to 5 μ m range
Battery charging voltage and current	2.93 \pm 0.19 μ m
Battery charging time	DC 5.0V, Max. 2A
Battery Non-stop usage time	Approx. 180 mins
Lithium Battery Cycle life	Approx. 120 mins
Durable period	The cycle times is not less than 300 times.
	Durable Period as follows: The product is used to nebulize 0.9% saline 2 times a day, 10 minutes each time at room temperature (23°C). Durable Period may vary depending on environment and solution.
	Main unit: 36 months Medication cup: 12 months
Noise level	\leq 50dBA (acc. to DIN EN13544-1)
Dimension	Approx. L40 \times W75 \times 92 mm
Weight	Approx. 84g
Operating temperature and humidity	10 ~ 40°C, 30 ~ 85 % RH, 800 ~ 1060hPa
Storage and delivery condition	-20 ~ 60°C, 20 ~ 75 % RH, 800 ~ 1060hPa
Safety level	Type BF Class II

11. TECHNICAL DATA

In accordance with EN 13544-1:2007+A1:2009 'Respiratory therapy equipment – Part 1: Nebulizing systems and their components', Annex CC.3, the multistage cascade impactor method was used to measure particle size.

Particle Size	MMAD 2.93 \pm 0.19 μ m (MMAD = Mass Median Aerodynamic Diameter)
Medication Cup Capacity	6 ml max
Noise	Less than 50 dB (1 meter distance)

In order to regulate the requirements for EMC with the aim to prevent unsafe product situations, the EN60601-1-2 standard has been implemented. The nebulizer MESH-S600A conforms to the EN60601-1-2:2007 standard for both immunity and emissions. Nevertheless, do not use the nebulizer close to the strong electrical or electromagnetic fields. This may result in incorrect operation and create a potentially unsafe situation. Guidance and manufacturer's declaration of MESH-S600A

Guidance and manufacturer's declaration-electromagnetic emissions

The **MESH-S600A** use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

The **MESH-S600A** is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies building used for domestic purposes.

Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The MESH-S600A use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The MESH-S600A is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies building used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions IEC 61000-3-3	Compliance	

Guidance and manufacturer's declaration-electromagnetic immunity

The MESH-S600A are intended for use in the electromagnetic environment specified below.

The customer or the user of the MESH-S600A should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge(ESD)IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s)to line(s) ± 2 kV live(s)to earth	± 1 kV differential mode Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips,short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT(>95% dip in UT) for 0.5 cycle 40% UT(60% dip in UT) for 5 cycles 70% UT(30% dip in UT) for 25 cycles <5% UT(>95% dip in UT) for 5 s	<5% UT(>95% dip in UT) for 0.5 cycle 40% UT(60% dip in UT) for 5 cycles 70% UT(30% dip in UT) for 25 cycles <5% UT(>95% dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment.If the user of the MESH-S600A requires continued operation during power mains interruptions, it is recommended that the MESH-S600A be power form an uninterruptible power supply or a battery.
Power frequency(50/60HZ) magnetic field IEC 61000-4-8	3 A/m	3 A/m	The MESH-S600A power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration-electromagnetic immunity

The MESH-S600A are intended for use in the electromagnetic environment specified below.

The customer or the user of the MESH-S600A should assure that is used in such and environment

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150KHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the MESH-S600A including cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d=1.2\sqrt{P}$ $d=1.2\sqrt{P} \text{ 80MHz to 800 MHz}$ $d=2.3\sqrt{P} \text{ 800MHz to 2.5 GHz}$
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5 GHz	3 V/m	<p>Where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in metres(m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and [and] mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MESH-S600A is used exceeds the applicable RF compliance level above, the MESH-S600A should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the MESH-S600A.
- b. Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distance between
Portable and mobile RF communications equipment and the MESH-S600A**

The MESH-S600A are intended for use in an electromagnetic environment in which radiated RF Disturbances are controlled. The customer or the user of the MESH-S600A can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MESH-S600A as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output Power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $D=1,2\sqrt{P}$	$D=1,2\sqrt{P}$	$D=1,2\sqrt{P}$
0,01	0,12	0,12	0,12
0,1	0,38	0,38	0,38
1	1,2	1,2	1,2
10	3,8	3,8	3,8
100	12	12	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
NOTE1 At 80 MHz and 800MHz, the separation distance for the higher frequency range applies.
NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

12. EU DECLARATION OF CONFORMITY

Manufacturer:

Dongguan Aicisy Machinery & Electronic
Equipment Co., Ltd.

Address:

Part B, 3rd F, Block A, Wentang Industrial
Park, Longhua Road, Zhouwu Dongcheng
District Dongguan City, Guangdong Province,
523118 P.R. China

SRN: CN-MF-000026848rial resources.

File No.: AIDISY-DCN-CE-02 Version: V5.0

Whose Single Authorized Representative:

MedPath GmbH

Address:

Mies-van-der-Rohe-Strasse 8, 80807 Munich,
Germany

SRN: DE-AR-000000087

We, the manufacturer, declare at our sole responsibility that following products

Product name	Model	Basic UDI-DI	Intened purpose
Medical Ultrasonic Nebulizer	MESH-S600A, MESH-S600B, MESH-S600C, MESH-S600D	Basic UDI-DI : 69571303MESH001QH	The Medical Ultrasonic Nebulizer is designed to aerosolize liquid medications for inhalation by the patient

meet the provisions of Regulations (EU) 2017/745.

The medical device has been assigned to class IIa according to rule 20 in Annex VIII of MDR 2017/745. It bears the mark

The product concerned has been designed and manufactured under a quality management system according to Annex IX of Regulations (EU) 2017/745.

Compliance of the designated product with the Annex IX of Regulation (EU) 2017/745 has been assessed and certified by the Notified Body

SGS Belgium NV
SGSHouse Noorderlaan
87 2030 Antwerp Belgium

CertificateNo.: CN19/41056
Issuedate: SGS Belgium NV Notified Body 1639
Expirydate: 24 MAY2024

The above mentioned declaration of conformity is exclusively under the responsibility of

Company: Dongguan Aidsy Machinery & Electronic Equipment Co.Ltd
Address: Part B, 3rd F, Block A, Wentang Industrial Park, Longhua Road, Zhouwu
Dongcheng District Dongguan City, Guangdong Province, 523118 P.R. China

2024/05/01

Place, date

A handwritten signature in black ink, appearing to be 'D. Aidsy', written over a horizontal line.

Legally binding signature, Function

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Tvrtka: Dongguan Aidisy Machinery & Electronic Equipment Co.,Ltd