## NebuNova

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Mesh Nebulizer

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User manual

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Firma: Dongguan Aidisy Machinery & Electronic Equipment Co.,Ltd  
Adres: Part B, 3rd F, Block A, Wentang Industrial Park, Longhua Road, Zhouwu  
Dongcheng District Dongguan City, Guangdong Province, 523118 P.R. Chiny

2024/05/01

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## 01. INTRODUCTION

Thank you for choosing Mesh Nebulizer. Our Nebulizer is a handheld device designed to aerosolize medication for respiratory therapy purposes, which considers the difference of your needs between traveling and indoor use.

This device can be operated on an internal lithium battery. Due to the use of a battery and compact size, you are able to continue your treatment while you are traveling. To help save the planet.

As this device is a medical instrument, please read the instruction manual before use, and be sure to follow the instructions of a doctor and use the device correctly.

### Product Features

1. Microporous structure (Mesh) and oscillator are combined in the atomization technology.
2. USB charging for convenient re-charge.
3. Compact body design, handy to carry.
4. Palm-size device, easy to carry.
5. Easy to operate.
6. Easy to clean & maintain.

Medical purpose and principle	The device adopts active ultrasonic vibrating mesh technology, which is designed to nebulize the solution into an aerosol for respiratory therapy purpose.	
Intended user	Legally certified medical experts, such as doctors, nurses, therapists, or patients under the guidance of qualified medical experts. The user should also be able to understand the general operation of this device and the content of instruction manual.	
Intended patients	This device is designed for all ages of patients, except for the patients who are unconscious and not able to breathe spontaneously or having pulmonary edema.	
Recommended Operation Environment	This device is intended for the use in the medical facility, such as hospital, clinic and doctor's office, a general family room, or open-air environment with a roof.	
Durable period	Durable periods are defined as follows. The product is used to nebulize 0.9% saline twice a day, and 10 minutes each time at room temperature (23°C). Durable period may vary depending on environment and solution.	
	Main unit:	36 months
	Medication cup:	12 months
Precautions for use	Warnings and cautions described in the manual should be read carefully.	

## 02. SAFETY PRECAUTIONS

### ⚠ READ ALL WARNINGS AND INSTRUCTIONS BEFORE USING THIS DEVICE.

As with any medical device, this product may become unusable due to a power failure, battery running out, or mechanical impact. Therefore, a backup device is recommended.

When using electrical products, always follow basic safety precautions. Take particular care of children when operating electrical devices.

### ⚠ WARNINGS



- Only use this device for medications prescribed by your doctor.
- The nebulizer is intended for respiratory therapy only, and any other applications of this device are improper and dangerous. The manufacturer holds no liability for any damage caused by improper or incorrect use.
- Do not share your nebulizer with others. It is intended for private use only. In the case that it is shared by multiple people, there is a risk of infectious illness spreading.
- Clean all parts of your nebulizer before and after each use and after extended storage.
- Adult supervision is required when this device is used by children and individuals who require special assistance.
- Never operate this device if any of the parts are not working properly or have been damaged.
- Be sure the device has been properly cleaned before use to avoid possible contamination.
- Do not plug or unplug the adapter into or from the electrical outlet with wet hands.
- Do not drop medication on the main unit or USB port; if it happens by mistake, immediately wipe it off with gauze. If you use the device while it is still wet, it may cause trouble or injury.
- Do not leave the device in the car if it is subjected to significant heat or cold.

- The battery is not replaceable; do not attempt to replace the battery by yourself, as you may damage or overheat it, causing injury. NEVER incinerate the battery, and NEVER immerse the battery in water; otherwise, it may cause the battery to rupture or damage.
- All batteries deteriorate over time if they are not used or charged; do not store the device for a long period of time without charging it periodically. If the device is not used or charged for a long time, the battery may create a hazardous condition.
- Never use the device while it is being charged.
- If you don't use the device for a long time, disconnect the USB charger.
- The device contains sensitive components, including a stainless-steel disk with 3–5 micron holes carefully measured. Do not drop, crush, puncture, bend, heat, incinerate, or put heavy weights on the device or its parts.
- Do not wash the medication cup and mesh with running water or any other liquid.

#### CAUTIONS

- Do not attempt to clean the mesh with any foreign objects; otherwise, it may damage the mesh.
- Do not drop the device; otherwise, it may not function normally.
- Keep the device away from direct sunlight, excessive heat, or cold to avoid damaging the batteries.
- Do not attempt to open, repair, or modify this device.
- Do not drop any liquids on the main unit or the adapter; if it happens by mistake, wipe it off immediately.
- Follow local laws and recycling plans for the disposal or recycling of components, batteries, and packaging.
- If the USB wall charger and USB port do not fit easily, they probably do not match. Make sure that you use the USB wall charger provided and position it correctly.

#### 03. CLASSIFICATION AND EXPLANATION OF SYMBOLS

	Warning / Caution / Note
	Class II equipment per IEC 60601-1
	Type BF equipment per IEC 60601-1
<b>IP22</b>	Protected against foreign objects equal to or greater than 12.5mm in diameter and drops of water falling at up to 15° vertically
	Consult instructions for use
	ON / OFF Button
	Keep dry
	Operating temperature limits: 10°C to 40°C
	Storage and delivery temperature limits: -20°C to 60°C
	Operating humidity limits: 30% to 85% R.H
	Storage and delivery humidity limits: 20% to 75% R.H
	Manufacturer
	Date of manufacture
	Serial number

	Batch code
	Disposal of Electrical & Electronic Equipment (WEEE): Do not treat this product as household waste.
	CE marking, indicating that this device complies with the requirements of the Medical Devices Directive (93/42/EEC)
	Authorized representative in the European Community


#### 04. PACKAGE CONTENTS / SYSTEM OVERVIEW

Check before use



The following items are contained in the package. Please check all parts for visible damage. Replace any damaged parts before you use this device. In the case of missing parts, malfunction or damage, please contact the store where you purchased the unit or the nearest dealer.

1. Main unit


2. Medicine cup module

3. Instruction Manual

4. Mask  
(Adult and Children)

5. Mouth piece

6. USB cable

● USB Port

## 05. HOW TO OPERATE WITH BATTERY

About the battery:

The device has a built-in lithium battery that you cannot replace, also NEVER use the device when it is being charged.

The device is charged through a USB port, just like many cellular phones and portable electronic media devices.

You can only charge the battery by connecting it to a wall outlet using the USB wall charger provided.

NOTE: the USB charger is an optional accessory. Only use external charger that is in accordance with the requirement of IEC60601-1:2012, output: 5V Max 2A.

To charge the battery with a wall outlet:

1. Open the USB cover to expose the USB port.
2. Gently insert the USB cable into the USB port.
3. Gently insert the USB cable into the USB wall charger port.
4. Insert the USB wall charger into a wall outlet.

Understanding battery states:

If the device is not connected to a power source, the ON-OFF indicator light has three possible states:

No light: OFF

Blue light: in operation

Orange light: low power

When the device is connected to a power source, the ON-OFF indicator light has two possible states:

Orange light blinking: is charging

Blue light blinking: fully charged

NOTE: rechargeable batteries do not have unlimited life. Battery life and number of charge cycles vary with use.

## 06. CLEANING AND MESH UNBLOCKING PROCEDURES


**▲ Clean all parts of your portable nebulizer before and after each use and after extended storage.**

1. Set up the medication cup

Step 1: take out the medication cup and main unit.

Step 2: plug in the medication cup to the main unit.

Step 3: make sure the medication cup is contact with main unit property.

2. Fill the medication

Make sure that the indicator light is off before adding the medication.


- Open the medication cup cap.

- Fill the medication cup with medication DO NOT EXCEED THE LIMIT (Max. 10ml).

- Close the medication cup cap properly.

3. Install the mask or mouthpiece into the mask adapter

The device is now ready for use. Refer to the next section on how to inhale.


**⚠ Caution:** For type, dose, and regime of medication, follow the instructions of your doctor or pharmacist.

1. Tilt the unit slightly, as illustrated in the figure.


In this position, the vibrating element is immersed in the medication, and nebulization will start when power is on.

When the vibrating element is immersed in the medication, the unit can be used at any angle.  
Note: In some positions (e.g. upright), nebulization may stop after a short while. In that case, slightly tilt the unit again to re-immers the vibrating element into the medication.

2. Place the mouthpiece into your mouth, or use the mask to cover your mouth and nose.

3. If you push the ON/OFF (  ) button, nebulization will start. It may automatically power off when the medicine or liquid is nearly depleted.

Always be calm and relaxed at the time of inhalation. Breathe in slowly and deeply so that the medication can reach deep into the lungs. Hold your breath briefly, then breathe out slowly while removing the mouthpiece from your mouth. Do not breathe too quickly. Pause if you want to have a rest.

4. After inhalation, always switch off the nebulizer.

5. Clean the unit after each inhalation, yed completely.

1. Tilt the device slightly. When the solution is almost discharged completely, ensure that the residual solution is in contact with the mesh and sprayed completely.

- ⚠** 2. If the solution with high viscosity is used, nebulization may be reduced.

3. If excessive solution accumulates on the mesh, nebulization may stop. In this case, turn off the power and absorb the solution with gauze or a lint-free towel.

4. Do not drop medication onto the main unit or into its USB port. Using the device when it is still wet may cause trouble or injury.

## 07. CLEANING AND DISINFECTION

**⚠** If the device is not cleaned and disinfected correctly and frequently as indicated, microorganisms may remain in the unit and cause risk of infection.

Caution: **NEVER** use the device when it is charging.

**NEVER** pour liquid into it. Liquid may **DAMAGE** the main unit.

If not cleaned immediately after use, the atomized mesh will be blocked, which may lead to the absence of fog in the next use.

### 1. After every use:

Clean the medication cup with mineral water after each inhalation.

1. Pour out the residual solutions in the medication cup.
2. Pour some mineral water into the medication cup.
3. Turn on the device to nebulize the mineral water for 2 to 3 minutes to clean the mesh.
4. During use, keep the main unit in a dry state.
5. Remove the medication cup from the main unit.
6. Wash and rinse the medication cup with mineral water.
7. Shake off excess water and allow parts to fully air dry in a clean place.
8. Use gauze or a clean towel to wipe off stains on the main unit if necessary.
9. Make sure that all cleaned parts are completely dry before storing or using them next time.

### 2. If the atomized mesh is blocked, please refer to the following operations:

Clean the medication cup with Approx. 60% Edible white vinegar After use it for 9 times.

1. Pour out the residual medication in the medication cup.
2. Fill some mineral water into the medication cup.
3. Pour out the mineral water from the medication cup.
4. Fill 3ml of approx. 60% edible white vinegar into the medication cup.
5. Turn on the device to nebulize the 60% white vinegar for 5 to 10 minutes to clean the mesh.
6. Disconnect the adapter and remove the medication cup from the main unit.
7. Wash and rinse the medication cup with mineral water.
8. Shake off excess water and allow parts to fully air dry in a clean place.
9. Use gauze or a clean towel to wipe off the main unit if necessary.
10. Ensure that the medication cup and all other cleaned parts are completely dry before storing or using them again.

1. Keep the battery compartment dry at all times.

- ⚠**
2. Do NOT poke the mesh with fingers, cotton swabs, or other objects.
  3. Do NOT clean parts in a dishwasher.
  4. Do NOT use a microwave to dry any parts.
  5. DO NOT wash the medication cup and mesh with running water or any other liquids.

### Daily disinfections:

It is important to disinfect the medication cup on a daily basis.

#### 1. Disinfection by boiled water:

- Rinse the medication cup with distilled water.
- Boil a saucepan of distilled water.
- **NOTE: Do NOT boil the medication cup directly.**
- Carefully immerse the medication cup in the boiled water for a maximum of 15 minutes.
- Carefully remove the medication cup from the boiled water and shake off excess water.
- Allow parts to be cooled down and fully air-dried on a clean, dry towel, and place them out of the reach of children.
- Make sure that all cleaned parts are completely dried before storing or using them next time.

#### 2. Disinfection by alcohol:

- Rinse the medication cup with distilled water.
- Immerse the medication cup in 75% ethyl alcohol for 1 minute.
- Rinse the medication cup with distilled water again, shake off excess water and allow parts to be fully air-dried on a clean and dry towel.
- Make sure that all cleaned parts are completely dried before storing or using them next time.

#### 3. Clean:

Cleaning & disinfection of nebulizer accessories (including: nebulizer cup, mouthpiece, etc.):

After the use of the nebulizer accessories:

- Clean the nebulizer accessories with distilled water.
- Completely immerse the nebulizer accessories into Approx. 75% medicinal alcohol.
- Soaking for Approx. 30 minutes.
- Wash and rinse the nebulizer accessories with about 50°C warm distilled water.
- Shake out water and dry them with a clean and soft towel, or leave them air-dried in a clean place.
- Make sure that the nebulizer accessories and other cleaned parts are completely dried, then put all these parts in a dry sealed container.

- ⚠** 1. Alcohol is highly flammable. Do not use alcohol in close vicinity of open fire or smoke.

The following item is able to be disinfected by boiled water or alcohol:

Medicine cup module

The parts below are NOT able to be disinfected by boiled water or alcohol:

Main unit

USB cable

⚠ Do NOT rinse or immerse the main unit in any liquid.

#### 08. CARRYING AND STORING:

Close the mouthpiece cover, store the device and medication cup in a dry and clean environment.

- ⚠
1. Do NOT leave or carry the device containing residual liquid in the medication cup.
  2. Do NOT leave the device under direct sunlight, or in highly humid, extreme heat or cold environment.
  3. Keep this device away from fire, high electromagnetic fields and out of the reach of children.

#### 09. TROUBLESHOOTING

If any trouble occurs while you are using the device, please check the following list first.

Problem	Possible Cause	Action
Low atomization	Batteries are low (orange light)	Please charge the battery
	The stains on the electrode cause a fault connection	Use rubbing alcohol to clean electrodes
	The mesh holes have been clogged or stained	Refer to cleaning and disinfection procedure to clean the medication cup page
	The mesh is broken	Replace the medication cup
After turning on the power, power indicator is not on and no mist comes out	Batteries are too low	Please charge the battery
	Fault connection between adapter and main unit	Check and reconnect the USB charger to the main unit
After turning on the power, the power indicator is on but no atomization occurs	The medication cup is not installed properly	Refer to set-up procedure of the medication cup to re-install the medication cup
	The mesh holes have been clogged or stained	Refer to cleaning and disinfection procedure to clean the medication cup page
	The mesh is broken	Replace the medication cup
Power indicator turns into orange light	Low battery power	Please charge the battery

1. If the device does not nebulize normally after taking the above-mentioned procedure, please contact the store where you purchased the device or the nearest dealer.

- ⚠
2. Never operate this device if any of the parts are not working properly or have been damaged.

## 10. TECHNICAL DATA

<b>Model</b>	<b>MESH-S600B</b>
Product	Portable Ultrasonic Nebulizer
Method of operation	Active Ultrasonic Vibrating Mesh Technology
Power supply	Internal power supply, Lithium battery only
Power consumption	Approx. 2.0W
Mode of operation	<b>20 mins Auto off (with rechargeable battery)</b>
	Press the button once to turn on the device and enter the normal gear
	Press the button again to enter the low gear
	Press the button again to turn off the device
Vibrating frequency	Approx. 107 kHz $\pm$ 10%
	Approx. 107 kHz $\pm$ 10%
Capacity of medication cup	Max. 10ml
Nebulization rate	$\geq$ 0.25 ml/min
Particle size	Approx. (74.53 $\pm$ 3.01)% within 0.5 to 5 $\mu$ m range
MMAD	2.93 $\pm$ 0.19 $\mu$ m
Battery charging voltage and current	DC 5.0V, Max. 2A
Battery charging time	Approx. 180 mins
Battery Non-stop usage time	Approx. 120 mins
Lithium Battery Cycle life	The cycle times is not less than 300 times, Durable Period as follows: The product is used to nebulizer 0.9% saline 2 times a day, 10 minutes each time at room temperature (23°C), Durable Period may vary depending on environment and solution.
Durable period	Main unit: 36 months
	Medication cup: 12 months
Noise level	$\leq$ 50dBA (acc. to DIN EN13544-1)
Dimension	Approx. L39 X W59 X H107 mm
Weight	Approx. 90g
Operating temperature and humidity	10 ~ 40°C, 30 ~ 85 % RH, 800 ~ 1060hPa
Storage and delivery condition	-20 ~ 60°C, 20 ~ 75 % RH, 800 ~ 1060hPa
Safety level	Type BF Class II

## 11. EU DECLARATION OF CONFORMITY

File No.: AIDISY-DCN-CE-02 Version: V5.0

**Manufacturer:**  
Dongguan Aicisy Machinery & Electronic  
Equipment Co.,Ltd.  
**Address:**  
Part B, 3rd F, Block A, Wentang Industrial  
Park, Longhua Road, Zhouwu Dongcheng  
District Dongguan City, Guangdong Province,  
523118 P.R. China  
**SRN: CN-MF-000026848trial resources.**

**Whose Single Authorized Representative:**  
MedPath GmbH

**Address:**  
Mies-van-der-Rohe-Strasse 8, 80807 Munich,  
Germany

SRN: DE-AR-000000087

We, the manufacturer, declare at our sole responsibility that following products

Product name	Model	Basic UDI-DI	Intened purpose
Medical Ultrasonic Nebulizer	MESH-S600A, MESH-S600B, MESH-S600C, MESH-S600D	Basic UDI-DI : 69571303MESH001QH	The Medical Ultrasonic Nebulizer is designed to aerosolize liquid medications for inhalation by the patient

meet the provisions of Regulations (EU) 2017/45.

The medical device has been assigned to class IIa according to rule 20 in Annex VIII of MDR 2017/745. It bears the mark

The product concerned has been designed and manufactured under a quality management system according to Annex IX of Regulations (EU) 2017/745.

Compliance of the designated product with the Annex IX of Regulation (EU) 2017/745 has been assessed and certified by the Notified Body

**SGSBelgium NV**  
**SGSHouse Noorderlaan**  
**87 2030 Antwerp Belgium**

**CertificateNo.:** CN19/41056  
**Issuedate:** SGS Belgium NV Notified Body 1639  
**Expirydate:** 24 MAY2024

The above mentioned declaration of conformity is exclusively under the responsibility of

**Company:** Dongguan Aidisy Machinery & Electronic Equipment Co.,Ltd  
**Address:** Part B, 3rd F, Block A, Wentang Industrial Park, Longhua Road, Zhouwu  
Dongcheng District Dongguan City, Guangdong Province, 523118 P.R. China

2024/05/01

Place, date

A handwritten signature in black ink, appearing to be 'F. Function', written over a horizontal line.

Legally binding signature, Function

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Tvrtka: Dongguan Aikisy Machinery & Electronic Equipment Co., Ltd