ÚČINNOST:

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WHITE PEARL**PAP WHITENING GEL
REVOLUTION IN TEETH WHITENING**

PAP is a new, non-invasive teeth whitening technology, which is better for enamel health than most common whitening methods containing hydrogen-peroxide.

PAP also has anti-plaque and long lasting anti-bacterial benefits.

Not only whitens, but also removes deep enamel stains

Suitable for vegans and vegetarians

Neutral pH

Not tested on animals

KIT INCLUDES:

2x Tubes (each 40 ml) of Non-Peroxide PAP whitening gel with applicator
4x Thermoforming Whitening Trays (2 sets)
1x Instruction manual

HIGHLIGHTS:

- Professional at-home tooth whitening
- Fast, effective and gentle tooth whitening
- Provides up to 80 tooth whitening treatments
- Contains ZERO PEROXIDE
- Developed by US dentists
- Strengthens and remineralizes tooth enamel
- Safe on teeth, gums, veneers, caps and crowns
- Manufactured in the USA
- Perfect for red wine, tea and coffee drinkers and smokers
- Smile with confidence, look fresher and younger

**TEETH WHITENING INSTRUCTIONS:**

1.

Apply a thin line of gel half way up the inner front surface of the tray for both the upper and lower whitening trays.

2.

Position the tray inside your mouth and place firmly against the front of your teeth. Remove any excess gel that gets on your gums.

3.

Wear for 30 minutes.

4.

After whitening, remove trays and rinse mouth and trays thoroughly with lukewarm water.

WHITE AND HEALTHY SMILE IN JUST 3 FAST AND EASY STEPS:

1. Mold custom mouth trays to your teeth for a more precise fit, resulting in a better whitening experience.
2. Apply a thin line of gel to front side of each tray. Press tray to both upper and lower teeth for 30 minutes, once or twice a day.
3. The White Pearl Whitening Kit can whiten your teeth up to 5 shades. Use the product until you are satisfied with your results. White teeth could last up to two years depending on personal habits.

WHITENING TIPS:

- ◆ Moisture could reduce the effectiveness of the gel. Ensure both teeth and trays are clean and dry before whitening
- ◆ For 1-2 hours after treatment, avoid eating or drinking and tobacco use
- ◆ Store gels in a cool, dry place and avoid exposure to heat or sunlight

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CUSTOM FORM YOUR TRAYS:

1.

Place the tray over your upper teeth before placing it in boiling water. Practice with your tongue on the roof of your mouth and touching the inside of the tray. Suck the air from the tray. This contours the guard to the inside surface of the teeth. When finished repeat with the lower tray.

2.

Bring water to a boil. The bowl should be at least 4" deep to ensure the trays full immersion. Remove the container from the heat and dip one tray in the hot (but not boiling) water holding the tab, keep the tray in the water for 5 seconds. Make sure the edges do not fold onto each other.

3.

Lift the tray vertically from the water using the tab. Carefully and quickly place the tray around the teeth and then immediately perform the procedure outlined in Step 1. Repeat with the second tray.

4.

Once you are satisfied with the fit, trim the front tab away from the tray with scissors. Make sure it is flush with the front surface of the tray.

5.

Your custom formed trays are now ready for whitening gel.

WHITENING GEL INGREDIENT HIGHLIGHTS:

- ◆ Glycerin delivers teeth whitening agents deep into the enamel, providing fast teeth whitening results
- ◆ PAP, which stands for Phthalimido-Peroxy-Caproic acid, is a new, non-invasive teeth whitening technology, which is better for enamel health than most common whitening methods containing hydrogen-peroxide.
- ◆ Xylitol natural sweetener derived from vegetable sources has low impact on blood sugar, reduces tooth decay and treats dry mouth. It is also effective in remineralizing deeper layers of enamel
- ◆ Mint adds a refreshing taste plus has the medicinal and germicidal benefits of diminishing harmful bacterial growth inside the mouth

WHITENING GEL INGREDIENTS:

Glycerin, Aqua, Hydrated Silica, Xylitol, Phthalimidoperoxycaproic Acid, Sorbitol, Hydroxyapatite, Potassium Citrate, Mentha Piperita (Peppermint) Oil, Ammonium Acryloyldimethyltaurate/VP Copolymer, PVP, Sodium Hexametaphosphate, Piroctone Olamine, Titanium Dioxide, Sodium Lauryl Sarcosinate, Mica, Saccharin Sodium, Arginine, Potassium Hydroxide

CAUTIONS:

Products are not recommended for use by children under 12 years of age. Use for periods of longer than 14 days is to be only under the supervision of a dentist. If irritant ion (such as redness, swelling, soreness) of the gums or the mouth occurs, discontinue use and consult a dentist. Avoid swallowing the cosmetic or part thereof. Avoid direct contact of the active surface of the tooth whitening product with the gums and/or salivary flow. Avoid contact of the product with the eye. Keep gel out of reach of children. Do NOT use if pregnant or during nursing. Do NOT smoke, drink or eat while whitening your teeth.

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