Rénergie Nano-Resurfacer Special Page

ENG

LANCÔME RÉNERGIE NANO-RESURFACER 400 BOOSTER

The Lancôme Rénergie Nano-Resurfacer 400 Booster is a revolutionary skincare device designed to deliver clinical results in the comfort of your home. Enter the new era of skincare.

THE FUTURE OF SKINCARE

This advanced tool utilizes innovative micro-surfacing technology to target key skin concerns such as texture, lines, wrinkles, firmness, and uneven skin tone. It enhances the skin's appearance for a smoother, brighter, and more unified complexion.

PREPARATION PHASE

1. Unpack your Nano-Resurfacer, its charger, 12 extra cartridges and 2 guides.
2. Charge your device (Charging time is three hours).
3. Disconnect the USB charging cable. The device will not work if connected to the electrical grid.
4. To prevent skin reactions, we recommend first testing the device on the inner side of your forearm to ensure your skin does not react irritably.
5. Turn your Nano-Resurfacer on and enjoy the procedure.
6. Dispose of the cartridge after each use.

HOW TO APPLY

1. Thoroughly cleanse and dry your skin.
2. Apply Rénergie H.C.F. Triple Serum to your face and neck. (Press the serum applicator repeatedly until all three components of Rénergie H.C.F. Triple Serum are released).
3. Remove a new, sterile nano cartridge from its protective packaging and insert it into the device. (Each cartridge is intended for single use only).
4. Turn on the NANO RESURFACER | 400 BOOSTER device and set the highest speed - level 3 using the "+/-" button.
5. Place the device with the nano cartridge on your face and move it in circular motions from the center of the face towards the cheekbones to cover the entire skin surface. Then repeat the procedure on the other side of the face.
6. Focus on target wrinkles and treat them with circular motions outwards. Repeat on the other side of the face.

AFTER USE

1. Turn off the device and remove the cartridge. (The device automatically turns off after 5 minutes).
2. Never reuse a cartridge.

WARNING

1. Avoid the lip area, under eyes, eyelids, and direct contact with eyes.
2. Do not go over the same part of the face more than 3 times and do not stay longer than 3 seconds.

THE PERFECT ROUTINE FOR THE PERFECT SKIN

1. Apply Rénergie H.C.F. Triple Serum

An anti-aging serum that uniquely combines Hyaluronic Acid, Vitamin C + Niacinamide, and Ferulic Acid in an innovative triple-dose formula, renewing skin layer by layer to plump, progressively correct wrinkles, and reduce dark spots.

2. Rénergie Nano-Resurfacer procedure

A revolutionary skincare device designed to deliver clinical results in the comfort of your home.

3. Rénergie H.P.N. 300 - Peptide Cream

An anti-aging cream that combines Hyaluronic Acid, over 300 peptides, and Niacinamide to deliver high-performance regeneration, accelerating skin renewal for visible lifting, wrinkle correction, and dark spot reduction.

4. Rénergie H.P.N. UVMUNE Cream SPF50

An anti-aging cream that prevents wrinkles and dark spots with high UV protection, utilizing peptides, Hyaluronic Acid, and Niacinamide to target photo-aging signs caused by the longest UVA rays.

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