



Hair growth brush with light therapy

Model: HC2006



Product description

Healthier-looking hair is within reach. Red LED light therapy. Helps reactivate dormant hair follicles. Scalp massage function.

Recommended use

Frequency: 3–4 times a week Session length: max. 30 minutes per day Compatibility: Can be used alone or with hair serums, minoxidil, or scalp products

Controls

One button (Power on/mode selection) LED indicators on the display for battery/charging status and selected mode Modes: Light only (red + NIR therapy) Bio-Ems (combined therapy) Massage only (scalp relaxation)

Usage description

Press and hold power button for 1,5s to turn on the device. And then short press the power button to switch modes.

Technology

Light Therapy:

Wavelengths: Red laser (660–785 nm) LLLT – low-level laser therapy, Light source: LED emitters, Penetration depth: 0.5–1 mm (scalp tissue), Benefits: Stimulates hair follicles, improves blood circulation, and promotes cell metabolism

Bio-EMS:

Simulates human bioelectricity, stimulates acupuncture points on the scalp, revitalizes hair follicles, regulates sebum production, improves blood circulation, and promotes natural hair growth.

Sonic massage function:

Frequency: 6,000–7,500 vibrations per minute, Benefits: Increases scalp stimulation, promotes nutrient absorption, relieves scalp tension.

How it works

Red LED light stimulates hair follicles and improves microcirculation in the scalp, blue LED light has an antibacterial effect and contributes to scalp health, regular use strengthens hair roots and promotes natural hair growth.

Main advantages

- Non-invasive technology with no side effects.
- Suitable for both men and women.
- Easy and convenient to use at home.
- Elegant and ergonomic design in the spirit of Morgan's quality.

Parameters

- Bristle material: Highly durable ABS with rounded ends
- Weight: approx. 110 g (lightweight ergonomic design)
- Dimensions: approx. (L)215 mm×(W)65 mm×(H)42 mm

Power supply and battery

- Rechargeable battery Type: Lithium-ion, 400 mAh
- Charging time: 3–4 hours
- Usage time: 7–10 sessions (15 minutes each) per charge
- Charging method: USB-C charging cable or magnetic charging base
- Input voltage: 5 V DC, 1 W

Safety precautions

- Clean the scalp before use and keep it dry.
- During use, place the teeth of the comb close to the scalp.
- Can also be used with a hair nourishing solution.
- Do not use for more than 30 minutes per day.
- Do not shine the light into your eyes.
- This is a hair care device, not a medical device.
- Do not disassemble, repair, or modify this device.
- Do not throw it into fire or water.
- After use, wipe the device with a damp paper towel.