

EN

CAUTIONS

- Please read and follow these safety guidelines before using the device.
- Use the device only as described in the instructions. Misuse may cause malfunction or injury.
- For cosmetic use only. Not a medical device.
- Do not use if you are pregnant, breastfeeding, or under 18 years old.
- If you have a skin condition such as lupus, photosensitive eczema, or albinism, consult a doctor before use.
- Do not use if you have photosensitive skin or are taking medication that causes light sensitivity.
- Avoid using the mask on the same day as skincare products that may increase light sensitivity.
- Always ensure the eye protection is correctly positioned and keep your eyes closed during use.
- If you experience pain, irritation, or discomfort, stop using the device immediately.
- The device is not waterproof. Do not immerse it in water or rinse it with water.
- Keep out of reach of children.
- Do not disassemble, modify, or repair the device yourself.
- For best results, use 3–5 times per week for about 10 minutes per session.
- After use, it is recommended to moisturize your skin.

PRODUCT PARAMETERS

Product:	LED 7-Color Face Mask	Power:	4W
Model:	MJ06	Auto Timer:	10 minutes
Battery:	500 mAh Lithium Battery	Charging Time:	2 hours
LEDs:	70 pcs (210 LED chips)	Net Weight:	180 g
Input:	5V = 1A	Material:	ABS + Silicone
Output:	0.3–0.6A		

OTHER LANGUAGES

For manuals in additional languages, please visit:
www.breakoutaid.com/manual
or scan the QR code below.



breakout+aid
love your skin



CLEAR GLOW LED MASK INSTRUCTION MANUAL

Please read the instruction manual carefully before use and keep it properly for future reference.


WHAT'S INCLUDED:

- 1 × LED Light Therapy Facial Mask
- 2 × Eye Comfort Cushions
- 1 × USB Charging Cable
- 2 × Adjustable Head Straps
- 1 × Instruction Manual

PREPARATION BEFORE USE

Please charge the light therapy facial mask fully before first use. During charging, the indicator light on the mask controller will remain solid red, and once charging is complete, the indicator light will turn off. The mask comes with two adjustable nylon straps for comfortable wearing.

PRODUCT EFFICACY



The Light Therapy Facial Mask device is based on the principles of light therapy (LLT), offering seven different colors to address various skin issues, aiming to help users achieve a radiant complexion. It is safe, easy to use, and can be repeatedly used for facial skin maintenance.

In addition to the three primary colors of red, blue, and green light, it also forms four compound lights through different combinations, totaling **seven modes: red, blue, green, yellow (red + green), purple (red + blue), cyan (green + blue), and mixed light (red + green + blue).**

Each light color has different effects and can address various skin concerns. Users can freely choose according to their skin condition.

LIGHTING INTRODUCTION

-  **Red light (wavelength 630nm):** Supports collagen and helps improve the appearance of fine lines. Great for skin that looks tired or needs a rejuvenating boost.
-  **Blue light (wavelength 460nm):** Helps reduce breakout-causing bacteria and calm the look of blemish-prone skin.
-  **Green light (wavelength 520nm):** Helps balance the appearance of uneven skin tone and supports oil control.
-  **Yellow light (wavelength 520+630nm):** Helps brighten the look of dull skin and support a more radiant complexion.
-  **Cyan light (wavelength 460+520nm):** Helps refresh and energize the look of the skin.
-  **Purple light (wavelength 460+630nm):** A combination of red and blue light that helps support clearer-looking skin.
-  **White light (mixed spectrum 460-630nm):** Penetrates deeper to support overall skin rejuvenation and improve the appearance of texture.



HOW TO USE

Step 1 – Cleanse Skin: Wash your face and gently pat dry.

Step 2 – Fit the Mask: Attach the adjustable straps and place the mask comfortably over your face.

Step 3 – Turn On & Select Mode: Press and hold the power button to turn the device on. Short press the button to switch between light modes. Use the mask for 10 minutes. The device will turn off automatically.

Step 4 – After Use: Remove the mask and continue with your regular skincare routine. Please hydrate your skin promptly after use.

Tips: Before using the light therapy facial mask, please shine it on your face for 1 minute to conduct a photosensitivity test. If there are no adverse reactions, you can proceed with the treatment. Depending on your skin condition, choose the appropriate light mode. We recommend using it 3-5 times a week for 10 minutes each session. With regular use over several weeks, this device may help improve the appearance of your skin.



Press and hold the power button to turn the device on/off. Press briefly to switch color modes. The controller emits a single 'beep' when powering on or switching colors, and two 'beeps' when powering off.



For the initial use, please ensure the face mask is fully charged. While charging, the indicator light will be solid red, and it will turn off once charging is complete. Do not use the product while it is charging.

PHOTOSENSITIVITY TEST

Before using this product for the first time, please conduct a photosensitivity test:

1. Place the face mask on your forearm and turn it on, then expose it for 3 minutes.
2. If your skin shows significant redness or itching, it may indicate photosensitive skin. Please stop using it immediately and consult your doctor.
3. After use, slight redness on the skin is a normal reaction and typically subsides within 1-2 hours. If you experience any burning sensation or discomfort during the process, please stop using it immediately.