BEAUTIFLY

Derma Glow

PRO

series

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- www.beautify.eu

- www.youtube.com/beautifyeu

Nominal voltage:	400 mAh
Nominal power:	175 mm x 46 mm x 55.5 mm
Battery capacity:	305 g
Size:	> 100 minutes
Weight:	2 hours
Working time:	4 interchangeable heads, USB cable, case,
Charging time:	interchangeable filters, user manual
The set includes:	
3.7 V	
< 1.6 W	

Safety: Before using the product, read this manual and follow the recommendations contained therein when using the device.

Product description [page 2 - f1]

1. Head / 2. Display screen / 3. Power button / 4. USB charging port / 5. Charging indicator (back of the device)

Interchangeable suction heads [page 2 - f2]

1. Small round suction head - Facilitates the removal of blackheads from hard-to-reach places.
2. Round suction head - Vacuum lifting gives a gentle, stimulating massage effect that improves microcirculation and tightens the skin.
3. Elliptical suction head - can be used on the corners of the eyes, the wings of the nose and the skin around the mouth.
4. Microcrystalline round suction head - Removes dead and calloused skin, acne scars.
5. Microcrystalline suction head.
6. Microcrystalline small round suction head

Charging

1. Charge the device for a maximum of 8 hours.
2. Make sure the device is turned off when charging.
3. You can connect the device to the power supply via a USB cable or directly to a computer.
4. During charging, the reverse charge indicator flashes, and after full charge, the charging indicator lights continuously.

Precautions: Dry your hands before using the device. Do not rinse or immerse the device under water. Perform daily skin care after using this device. Do not use other mechanical facial cleansing devices while using this device. Do not exert too much pressure on the skin during use, as the device itself has a strong suction. The suction cups should be disinfected to avoid the formation and spread of viral, bacterial and fungal skin diseases.

Contraindications: Active purulent acne; rosacea; retinoid therapy; acid treatments, lasers (at least 4 weeks before microdermabrasion); psoriasis; viral, bacterial and fungal skin diseases; vascular skin; tendency to scarring; broken skin continuity (fresh wounds); skin tumors or suspicious moles in the treatment area; facial fillers (at least 4 weeks before microdermabrasion, otherwise the preparation may move); do not wash the device with water, the device is not waterproof.

Warnings:

1. Before using the purchased product, read the product user manual and these safety warnings, the product maintenance and storage instructions, as well as the rules for handling used electrical and electronic equipment. Follow the recommendations contained therein and use the product in accordance with its intended use.
2. The product is not a toy and is not intended for use by children. The product may only be used by adults.
3. Do not use the product if it does not work or works incorrectly, or has been dropped or damaged in any other way. In such a situation, do not open or repair the product yourself. If the product is still subject to liability for defects in the sold item (warranty) or guarantee, contact the seller or guarantor. If the warranty or guarantee period has expired, or the defect is not caused by a cause inherent in the product, return the product for repair to an appropriate service center.
4. Do not use the product near sources of heat, water or moisture, and do not immerse it in water or other liquids.
5. The product is intended for use at room temperature, between 10 and 400 Celsius.
6. Do not expose the product to direct sunlight.
7. After use, store the product in a dry place, at room temperature, between 10 and 400, out of reach of children and animals.
8. If the product becomes dirty, disconnect it from the power source, wipe it with a damp cloth and leave it to dry before storing it. Do not use abrasive materials or corrosive cleaning agents, as this may damage the product.
9. Do not use the product if you notice a damaged power cord or plug, an unpleasant smell (e.g. burning), overheats or smoke is emitted from it. In such a situation, do not open or repair the product yourself. If the product is still subject to liability for defects in the sold item (warranty) or guarantee, contact the seller or guarantor. If the warranty or guarantee period has expired, or the defect is not a result of a cause inherent in the product, return the product to an appropriate service for repair.
10. Do not leave the product unattended while it is switched on.
11. The power source must be compatible with the required voltage of the product.
12. After use you disconnect the product from the power source and allow it to cool down before moving or storing it.
13. Do not twist or coil the power cord, as this may damage it.
14. Accessories for the product that do not come from the manufacturer may cause the product to malfunction or be damaged.

More information about the company and products [page 2 - f3]

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[EN] Rules for handling electrical and electronic equipment: Waste electrical and electronic equipment, as well as batteries and accumulators, must not be disposed of with other waste. This is indicated by the symbol of a crossed-out bin placed on the product packaging. Electrical and electronic equipment, along with batteries and accumulators, may contain hazardous substances, mixtures, and components harmful to the environment and human health. Improper handling or damage to these items may result in harm to health or environmental pollution.

Proper handling of waste electrical and electronic equipment, as well as batteries and accumulators, in households contributes positively to their reuse and recovery, including recycling.

Used electrical and electronic devices, as well as batteries and accumulators, can be returned at local selective waste collection points. Distributors are required to collect household waste equipment free of charge at the point of sale, provided that the used equipment is of the same type and served the same function as the equipment sold. Distributors delivering equipment intended for household use must also collect used household equipment free of charge at the place of delivery, provided that the used equipment is of the same type and served the same function as the delivered equipment. A distributor operating a retail unit with a sales area, as defined in Article 2, point 19 of the Act of March 27, 2003, on spatial planning and development, of at least 400 m² dedicated to the sale of household appliances, is required to accept free of charge, at that unit or in its immediate vicinity, waste equipment from households, provided that none of its external dimensions exceed 25 cm, without the need to purchase new household equipment.

The product was introduced to the market after August 13, 2005.

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[EN] Summary declaration of conformity: The manufacturer, under their sole responsibility, declares that the Derma Glow Pro product complies with the Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, the Directive 2014/35/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, and the Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, along with the Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU regarding the list of restricted substances, and the Directive (EU) 2017/2102 of the European Parliament and of the Council of 15 November 2017 amending Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment. The full declaration of conformity is available at <https://beautify.eu/deklaracje>.

Importer: Beautify Sp. z o.o., Plac Bankowy 2, 00-095 Warszawa, Polska
Model: Derma Glow PRO
Made in China.