BEAUTIFLY

Scrub Mini

Young
series

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- www.beautifly.eu

- www.youtube.com/beautiflyeu

Rated voltage:	AC 100V~240V 50HZ/60HZ
Output voltage:	5V/400mA
Battery:	lithium battery
Weight:	70g
Ultrasound:	25KHZ
Power:	2W
Size:	140*42*20mm
Kit Includes:	user's manual, USB cable, accessories

Safety

Please read this manual before using the product and follow the instructions in it when using the device.

Device description [page 2 - f1]

1. Stainless steel blade
2. Protective cover
- 3 Power button/mode switch
4. USB charging port

Product features

- Provides deep cleansing of the skin.
- Reduces pigment spots and breaks down pigment.
- Smoothes wrinkles and improves skin elasticity.
- Improves blood circulation and lymphatic metabolism.
- Increases the regenerative capacity and permeability of skin cells.
- Improves the condition and appearance of the skin.
- Is easy to carry and navigate.

Charging

Note: Fully charge the battery before using the device for the first time.

- Check the condition of the lithium battery before charging.
- Connect the USB charging cable to the USB interface.
- Connect to the USB interface of your computer, an adapter (sold separately) or the lithium battery of your mobile phone.
- The LED flashes continuously during charging and lights up when charging is complete.

Notes on charging

- Charge the device if it does not display a start-up status when the power button is pressed for a long time.
- During operation, the device will sound a warning tone and then the shutdown state will start. This indicates that the battery level is insufficient. The product must be recharged immediately.
- After 10 minutes of use, the product automatically enters standby mode and, if not used for one minute, switches off automatically.
- It takes approximately 120 minutes to fully charge. (It takes more than 5 hours to charge before first use).
- Once the battery is fully charged, the device can be used for approximately 60 minutes.
- It is recommended to recharge the device periodically.

Operating instructions

Important: The switch lights flash when changing modes.

Cleaning mode

1. Press the power button for 1.5 seconds to switch the unit on.
2. The first mode - cleansing will start automatically when the unit is switched on.
3. Daily care: reduction of blackheads, whiteheads and acne. The ultrasonic vibrations penetrate deep into the muscles, deeply cleansing the skin. Use together with a facial cleanser.

Moisturising mode

1. Press the power button for 1.5 seconds to switch on the device.
2. Press the power button once to start the second mode - Moisturising.
3. Ultrasonic vibrations and negative ion make the active ingredients in cosmetics penetrate deep into the skin.

Lift mode

1. Press the power button for 1.5 seconds to switch on the device.
2. Press the power button twice to activate the third mode - Lifting.
3. Micro-currents generated by ultrasound and EMS stimulate cell activity. This process stimulates collagen and elastin production, resulting in firmer and smoother skin.

Cleaning and maintenance

- Protect the device from moisture, do not place it in the dishwasher or microwave.
- Switch off the device before cleaning.
- Do not use harsh, caustic or volatile thinner cleaning materials.
- Use a soft cloth or cotton swab moistened with cold or warm water to gently wipe the unit.
- Do not rinse the device under running water.

Warnings:

1. Before using the purchased product, read the user manual, these safety warnings, maintenance, and storage instructions, as well as the rules for handling waste electrical and electronic equipment. Follow the provided recommendations and use the product as intended.
2. The product is not a toy and is not intended for use by children. It is designed for use by adults only.
3. Do not use the product if it is malfunctioning, damaged, dropped, or operating improperly. In such cases, do not attempt to open or repair the product yourself. If the product is still under warranty or statutory liability for defects, contact the seller or guarantor. If the warranty period has expired or the defect is not related to the product itself, take it to an authorized service center for repair.
4. Do not use the product near heat sources, water, or moisture, and do not submerge it in water or other liquids.
5. The product is intended for use at room temperature, within a range of 10 to 40 degrees Celsius.
6. Do not expose the product to direct sunlight.
7. After use, store the product in a dry place at room temperature (10 to 40 degrees Celsius), out of reach of children and pets.
8. If the product becomes dirty, disconnect it from the power source, wipe it with a damp cloth, and let it dry before storage. Do not use abrasive materials or corrosive cleaning agents, as they may damage the product.
9. Do not use the product if the power cord or plug is damaged, if it emits an unpleasant smell (e.g., burning), overheats, or emits smoke. In such cases, do not attempt to open or repair the product yourself. If it is still under warranty or statutory liability for defects, contact the seller or guarantor. Otherwise, take it to an authorized service center for repair.
10. Do not leave the product unattended while it is turned on.
11. The power source must comply with the product's voltage requirements.
12. After use, disconnect the product from the power source and allow it to cool before moving or storing it.
13. Do not twist or coil the power cord, as this may damage it.
14. Accessories not supplied by the manufacturer may cause the product to malfunction or become damaged.

More information about the company and products [page 2 - f2]

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[EN] Rules for handling electrical and electronic equipment: Waste electrical and electronic equipment, as well as batteries and accumulators, must not be disposed of with other household waste. This is indicated by the crossed-out wheeled bin symbol on the product packaging. Electrical and electronic equipment, batteries, and accumulators may contain hazardous substances, mixtures, and components that are harmful to the environment and human health. Improper handling or damage may lead to health risks or environmental pollution. Proper household handling of waste electrical and electronic equipment, batteries, and accumulators contributes positively to reuse, recovery, and recycling, including material reclamation.

Waste electrical and electronic equipment, as well as batteries and accumulators, can be returned to local selective collection points. Distributors are required to accept waste equipment from households free of charge at the point of sale if the waste equipment is of the same type and served the same function as the equipment being sold. Distributors delivering household equipment to consumers must also accept waste equipment from households free of charge at the delivery location, provided the waste equipment is of the same type and served the same function as the delivered equipment. Distributors operating retail units with a sales area of at least 400 m², as defined in Article 2, point 19 of the Spatial Planning and Development Act of 27 March 2003, dedicated to household equipment sales, are required to accept small waste equipment from households (with no external dimension exceeding 25 cm) free of charge at the retail unit or in its immediate vicinity, without the requirement to purchase new household equipment. This product was placed on the market after 13 August 2005.

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[EN] Summary declaration of conformity: The manufacturer, under its sole responsibility, declares that the Scrub Mini Young device complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, as well as Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU regarding the list of restricted substances. The full text of the declaration of conformity is available at <https://beautify.eu/deklaracje>.

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