BEAUTIFLY

CollagenBoost

PRO

series

DermaPen
DermaPen

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beautifly.eu

@beautifly.eu

beautifly.eu

/beautiflyeu

@beautifly.eu

f1

f2

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- www.beautifly.eu
- www.youtube.com/beautiflyeu

Input: 100 - 240V 50/60 Hz
 Output: 5V=1A
 Number of modes: 5
 Included: Dermapen device, two massage rollers, charger.

Security:

Before using the product, read this user manual and follow the recommendations contained therein when using the device.

Device description [page 2 - f1]:

1. Charging input.
2. Speed display.
3. LED lighting button.
4. EMS intensity display.
5. EMS intensity button.
6. Power/speed button.
7. Scale point.
8. Scale adjustment ring.
9. Needle insert.
10. Induction base.

Functions:

The device combines microneedling technology with microcurrents (EMS) and light therapy. There are 2 different types of needle cartridges designed for the device.

Microneedling technology improves the appearance of the skin and reduces the signs of aging. They stimulate the skin using microneedles to obtain collagen and produce elastin. Micro channels created in the surface of the skin help in the absorption of skin serum.

Microcurrents (EMS) can increase cell activity, stimulate muscle movement, restore skin elasticity and accelerate blood circulation in capillaries. The device has 3 different levels of EMS intensity. Use the appropriate intensity depending on the treatment and skin type. By pressing EMS, the intensity changes with each press.

Level 1: 1 light on.

Level 2: 2 lights on.

Level 3: 3 lights on.

The EMS function closes when all lights are turned off.

Press the LED light button to turn the light on or off. The color of the light changes each time you press the button.

Red light: reduces inflammation and stimulates collagen production.

Blue light: good for sterilizing and improving acne.

How to use the device:

1. Place the needle cartridge in place and lock it in place device by turning the cartridge clockwise.
2. Turn the needle cartridge counterclockwise to remove it.
3. Operate the device only when the needle cartridge is properly assembled. Press and hold the power button for 3 seconds to turn on/off the device.
4. When you turn on the device, the speed changes once each time you press the power button.
5. Use the appropriate needle size according to different treatment areas and skin problems.
6. Adjust the needle depth by turning the scale adjustment ring. After starting the treatment, adjust the needle depth. The scale number corresponding to the scale point

shows the current needle depth.

Landing:

1. This device is powered by a built-in battery.
2. Charging the device takes approximately 3 hours.
3. When the device battery drops to 10%, the speed number flashes slowly.
4. When the battery is low, the device will turn off automatically and the speed on the display will flash quickly for one minute.
5. The device shows "L" symbol when restarting after automatic shutdown.
6. The device shows "0" symbol while charging.
7. The device shows "F" symbol when charging is complete. The device will automatically resume operation after being removed from the base.

Recommendations:

1. Used needle cartridges should be placed in the garbage container.
2. If the needle cartridge is accidentally contaminated during installation (dropped on the floor, package open/damaged), discard it and purchase a new needle cartridge.

Contraindications:

- active inflammation and purulent processes during the procedure
- pregnancy and breastfeeding
- unregulated diabetes
- unregulated thyroid diseases
- tendency to develop keloid scars
- allergy to any of the ingredients of the substances used for the treatment
- in people with vascular skin
- use of drugs that reduce blood clotting
- active cancer diseases
- use of anticoagulant drugs (e.g. aspirin)
- fever
- if Botox has been administered to the treated area in the last three months

Warnings and recommendations:

1. Please check the device and charger before use.
2. Make sure the device is turned off when not in use.
3. Keep the device away from direct sunlight and high temperature or moisture.
4. Keep the device out of the reach of children.
5. Turn off the device after use or before cleaning and maintenance.
6. Do not use the device on the eyelids, thyroid area, lips, ears, breasts or genital area.
7. Do not repair, disassemble or modify the controller and body. This may cause a malfunction.
8. When using a device that has not been used for a long time, check whether it is working properly and safely.

More information about the company and products [page 2 - f2]:

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[EN] Rules for handling electrical and electronic equipment: Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled. The device was placed on the market after August 13, 2005

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[EN] Summary declaration of conformity: The manufacturer hereby declares under its sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances.

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Importer: Beautify Sp. z o.o.

Plac Bankowy 2, 00-095 Warszawa, Polska

Model: CoLagenBoost