medivon®

User manual
EYE MASSAGER

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1. O produkcie

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1. About the product

The modern, intelligent Horizon Pro eye massager is designed to massage the temples and other acupuncture points around the eyes. The innovative device offers pneumatic, vibrating and hot compress massage. These treatments relieve tired eyes and relax the user both physically and mentally. The product has relaxing music pre-loaded, which increases the comfort of the massage. The device looks like glasses, has a fashionable and practical design, and is easy to carry. The product is equipped with a built-in rechargeable battery and an intuitive backlit LCD screen. The integrated and compact device will serve everyone, it can be used anytime and anywhere, whenever you want. Allow yourself and your eyes to relax.

2. Purpose

The product is aimed at students, people who strain their eyesight at work and work long at the computer, writers, accountants, people sitting long in front of a computer or TV screen, the elderly experiencing eye fatigue, people who often stay awake late or otherwise abuse their eyesight.

3. Functionality

Pneumatic massage: pressure massage of acupuncture points around the eyes, consisting in alternately filling the chambers with air (pressure) and emptying the chambers (reducing pressure).


Vibration massage: Simulating impacts, shiatsu and other traditional massage techniques - provides muscle relaxation.

Hot compress: increasing the temperature of the device accelerates blood circulation around the eyes, bringing relief and relaxation.


Integrated and compact design; easy adjustment of the band will provide comfort to every user.

The material, which is in direct contact with the user's skin, prevents dirt and sweat stains, is removable and washable, making it easier to use and ensuring comfort.


4. The appearance and elements of the product

1. Clasp adjusting the pressure
2. Rubber band
3. Key on / off / mode change
4. Music key
5. USB charging interface
6. LCD screen

Velcro Tape

Removable and Washable Cloth

5. Preparation before use

Built-in rechargeable battery

To start the product, briefly press the "ⓘ" button. Check the battery charge indicator. If the battery indicator is blinking, it means that the battery needs charging - use the included USB cable for this purpose. It is recommended to use the device when the battery is fully charged.

Product fit and use

Attach the removable material on the inside of the device as shown in the diagram, use three Velcro straps to do this. Then put on the massager, adjust its position to the face and eyes, adjust the rubber headband to the head so that you feel comfortable during the massage.

6. Buttons and symbols on the screen

1. Buttons:

- "ⓘ": Start / switch mode; long press to turn off;
- "🎵": Music function key

2. Symbols (on the screen):

- "⚙️": Mode functions;
- "🕒": Timing;
- "🎵": Playing music;
- "🔋": Battery level;
- "Air": Indicates that the pneumatic massage function is working;
- "Heat": Indicates that the heating function is working;
- "Vib.": Indicates that the vibration function is activated.

7. Function buttons

1. On / Off button / mode change "ⓘ"

To start the product shortly press the "ⓘ" button, the LCD screen will turn on, the product will start running in the default Mode 1 with the following options:

Turn on the pneumatic massage

Turn on the hot compress function

- Turn on the music
- The remaining runtime is 15 minutes.
- The LCD screen shows "air" and "heat".

· Each time the button "ⓘ" is pressed, the operating mode can be changed (there are 5 of them). The device has the following modes:

Mode 1: pneumatic + heating

Mode 2: all functions

Mode 3: pneumatic + vibrating

Mode 4: pneumatic

Mode 5: heating

- Hold the "ⓘ" button for 2 seconds to turn off the device.

2. Music button "🎵"

- When the music is off, press the button once to start playing;
- Pressing this button while playing music will go to the next track;
- Press the button twice in succession while playing music to adjust the volume; By default, the volume is set to a high level. The volume level is switched from high to medium to low in sequence.
- When playing music, long press the music key to stop playing the music;

8. Uploading your own songs

Connect the product to the computer with the USB cable (as shown below) then copy your favorite music to the device. The product only supports MP3 format music.

NOTE: The product turns off when connected to another device with the USB cable.

9. Care and cleaning

To clean the product, wipe it with a slightly damp soft cloth or towel. Do not use volatile preparations. Store the product in a dry and ventilated place and do not expose it to direct sunlight, high temperature or moisture.

10. Maintenance of the removable material:

After use, wipe sweat and other removable cloth stains with a slightly damp cloth. The material can be washed according to the following instructions:

- The material can only be washed by hand, do not machine wash. · ·
- The water should be warm when cleaning. do not scrub or rub vigorously.
- Do not wring.

- Do not use a decolorizing agent;
- Do not dry mechanically;
- Do not iron;
- Hang freely to dry.

11. Battery charging instructions

Low battery level is signaled during massager operation by a vibration signal repeated three times, the battery charge level icon then displays "□". Use the USB cable and connect the device to another charger to charge the battery; while charging, the battery icon will blink. When the battery is fully charged, the battery level icon will display "▢". Charging the battery takes approx. 2 hours.

- A fully charged battery can be used for 4-5 days (15 minutes each day).
- Normally, the service life of a rechargeable battery is ≥ 500 full charges and discharges, i.e. about 3 years.
- The service life will be shortened if the battery is not used for a long time. It is recommended to charge the battery at least once a month.

12. Precautions

Please read this instruction manual completely before operating the product.

- Clean the cover of the massager before use.
- The product is not intended for children and people without the ability to express opinions.
- It is forbidden to use the massager by people with hypersensitive or heat-insensitive skin.
- It is forbidden to use the product by pregnant women and patients with heart disease, dysarteriotonia, malignant tumors, vascular diseases of the skull and other acute diseases and people under

the supervision of a doctor without prior consultation.

- Stop using the massager immediately if you feel uncomfortable or notice changes in your skin during the massage. Consult a doctor as soon as possible.
- Do not use or only use the product under medical supervision if you have vision problems. It is forbidden to use the device by people after eye surgery and patients suffering from glaucoma, cataracts, retinal detachment and other eye diseases.
- Contact lens wearers must remove them before wearing the massager.
- Do not use the Product in places with high temperature, flammable or humid, as well as in places exposed to electromagnetic radiation.
- Do not disassemble the device without authorization.

Before use, make sure the product is not damaged or wet to reduce the risk of failure.

13. Scientific Knowledge: Avoid "Computer Vision Syndrome"

I. After a long time in front of the computer, we feel discomfort.

Ophthalmologists and optometrists call this "computer vision syndrome (CVS)". The details of CVS are as follows:

① Eyestrain will increase if the person looking at the screen for a long time does not blink. When viewing content and images on the screen, the computer user blinks as little as possible to avoid missing information.

Blinking is critical to protecting your eyes. It allows you to change the position of the eyeball, reduce the tension of the lens and the iris muscle, moisturize the eyes and prevent the feeling of burning, itching and similar eye ailments.

② The "McCullough Effect" occurs when you see green words on a black screen, and this is a normal reaction. Although harmless to the body, this symptom may last for several days. The same effect will be given by a piece of paper stuck to a white wall, if, when looking at the paper, we see a pink edge around it.

③ Information for bifocal glasses wearers:

The optimal position for the display should be at or below the chest height with a distance of 12-18 inches. If the screen height is inappropriate or the distance is too short, the user will bend his neck incorrectly or hold his head at a rather unnatural angle, which may cause dizziness, neck and shoulder pain, or a general feeling of being unwell. Operating the computer strains the eyes and reduces the production of lubricant and enzyme in the eyes. By limiting the number of blinks, it exposes you to dry and tired eyes, visual disturbances, headaches and other symptoms.

II. After popularizing the computer, operating it is almost an indispensable part of modern life. Many people have to watch a computer screen all day as part of their work. Practicing this leads to eye fatigue and even diseases such as glaucoma. Therefore, people who spend a lot of time in front of screens should definitely pay attention to eye care. Main assumptions:

① Avoid using the computer continuously for long periods. Allow 5-10 minutes of rest for each hour of computer work. The recommended exercise is to focus your eyesight on a distant object, e.g. a building / tree outside a window.

② Maintain the correct posture while working. The computer screen should be straight ahead or slightly below the height of your eyes. The optimal reference point is the center of the screen, which is 20 ° below the horizontal line of sight.

③ Maintain a sufficient working distance. The eyes should be at least 60 cm from the computer screen.

- ④ Create and maintain perfect working conditions. The light in the workplace should be soft; the computer screen should have adequate brightness and display the data correctly. Do not work in front of the screen in a dark environment as the contrast of black and light on the screen could damage your eyes. The height of the desk and chair should suit the user and the height of the monitor.
- ⑤ Consult your ophthalmologist if you experience dry or red eyes, burning, blurred vision, heavy eyelid or foreign body invasion, eye swelling / pain, or headache despite following eye care recommendations.

Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out basket on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may lead to health damage or environmental pollution during further disposal or recycling of the equipment.

Armare S.A. Hereby declares that this massage device complies with the European Directive RE 2014/35/EU The full text of the Declaration of Conformity is available at the following internet address:
www.medivon.pl/declarations

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