medivon.

User manual

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Device Manual

Fully charge the power unit before using it for the first time and after a long period of non-use. After charging the battery, attach the power unit to the blindfold at the connector location.

After holding the power button for 2 seconds, the massager will start and allow you to start the massage. By repeatedly pressing the POWER button, you change the massager's operating modes:

- 1st vibration + heating
- 2nd vibration
- 3. heating

The set with the device also includes a cooling band. Just put it in the freezer a few hours before use. After cooling down, you can put it in the eye massager for relief and relaxation.

Recharge the device when the battery level is low. Do not fully discharge the battery. The supply voltage is 5V/1A. Use the factory charger for charging. Do not use the device while it is charging.

After the massage is completed, the device can be turned off by pressing the POWER button again.

The massager is ideal for people who need deep relaxation or improve blood circulation. With systematic use of the massage, you will reduce the symptoms of fatigue, muscle swelling or chronic pain of unknown origin. The construction of the massager is very light and suitable for everyone.

The time of the massage should be adjusted to your own preferences. The device must not be used by persons with injuries or who have undergone surgery in the last 3 months. Do not use the device in high temperatures and damp rooms. Absolutely keep away from water. Also, people with hypertension, hypotension or severe strokes cannot use the device. People with skin problems should consult a doctor before use.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment.

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The manufacturer hereby declares that the Medivon Hypnos device complies with Directive 2014/30 / EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65 / EU of the European Parliament and of the Council of 8 June 2011 on restriction of the use of certain hazardous substances in electrical and electronic equipment. The full text of the declaration of conformity is available at the internet address:
www.medivon.pl/deklaracje

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska
Model: Medivon Hypnos