medivon.

User manual

PL | EN | DE | FR | ES | PT | IT | CZ | SK
SL | EE | LT | LV | HU | RO | RU | NB | FI
EL | HR | SV | DA | NL | TR | BG | AR

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User Manual

The device is ready for immediate use after purchase. Make sure the device is charged, then start the massager and select the appropriate operating mode according to your preferences. When the massager is turned off, short press the POWER button to turn it on. The first default operating mode will be launched. To turn off the device, hold down the POWER button again. When the device is turned on, short cyclical presses of the power button change the predefined operating modes with various functions of the device.

If the device is discharged, the battery indicator will turn red. Connect the device to a power source to charge the battery. A full battery will be indicated by a green LED. The supply voltage is 5V / 2A.

The massage icon pressed cyclically changes the intensity. Three massage strengths are available to suit the user's preferences.

Pressing the heating button activates the heating mode. There are 3 heating levels. The pressing cycles alternate in sequence on the display. The fourth press turns off the heating function.

Using the vibration button, as in the case of heating, we change its intensity cyclically. A fourth press turns off the vibration function.

When the heating or compression massage function is turned on, EMS mode can be activated. The EMS+ and EMS- buttons are used to control the intensity of electromassage. The vibration and EMS frequencies have been automatically adjusted by specialists for the best possible use.

The strength of the massage should be adjusted to your own preferences. The device must not be used by persons with injuries or who have undergone surgery in the last 3 months. Do not use the device in high temperatures and damp rooms. Absolutely keep away from water. The device cannot be used by people with hypertension, hypotension, during pregnancy, with a pacemaker or after severe strokes. People with skin problems should consult a doctor before use.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after 13.08.2005.

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The manufacturer hereby declares that the Medivon Knee Pro device complies with Directive 2014/30 / EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65 / EU of the European Parliament and of the Council of 8 June 2011 on restriction of the use of certain hazardous substances in electrical and electronic equipment. The full text of the declaration of conformity is available at the internet address: www.medivon.pl/deklaracje

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska
Model: Medivon Knee Pro