medivon.

## User manual

Anti-cellulite massager

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## Device manual

1. Connect your device.
2. Install the appropriate head for the desired treatment.
3. Set the intensity using the dial depending on the area you will be working on and its sensitivity to massage.
4. After using the massager, turn the switch to off.
5. Disconnect the device.

- Finger head - stimulates a professional massage performed with the fingers, suitable for use on cellulite.
- Roller head - gives immediate results by eliminating cellulite, works against sagging, stimulates blood flow, tones muscles, performs lymphatic drainage, relieves back pain, sciatica and lumbar pain, warms muscles, relaxes, relieves stress, etc.
- Ball head - suitable for delicate areas. Activates and firms muscles, improves circulation in the abdominal area, eliminating excess toxins.
- Wavy head - firms and improves blood flow, suitable for massage of thighs and buttocks.
- Flat head - ideal for a relaxing massage of the neck, shoulders and back.

For greater comfort, wear a fabric cover that cannot be used with creams. Intensity adjustment allows you to adjust the intensity of the massage depending on your needs. Use the low setting to massage delicate areas (breasts, legs with varicose veins). Use the highest one to massage less sensitive areas of the body.

To clean the massager, use only a soft cloth and neutral detergent.

## Warnings and Safety Instructions

Store the device at room temperature, in a dry and shaded place, out of reach of children. Do not use the device if you are pregnant or suffer from heart problems or have any physical injuries. If you are unsure, consult your doctor before use. Before using, wait one hour after your last meal. Keep the massager away from children - it is not a toy! Do not use the device on the same area of the body for more than 15 minutes during one massage session. Do not use the device in the bathroom or in humid places. If the device does not work properly, immediately disconnect the power supply and contact the service center. Do not use the device if you feel unwell. Before replacing the head, make sure the device has stopped working. If the power cord is damaged, do not attempt to start the device. The cable must be replaced by the manufacturer, its service agent or a similarly qualified person in order to avoid danger to life or health.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after 13.08.2005.

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The manufacturer hereby declares that the Medivon Cello device complies with Directive 2014/35/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65/EU of the European Parliament and of the Council of June 8, 2011 on restrictions on the use of certain hazardous substances in electrical and electronic equipment. The full content of the declaration of conformity is available at the following internet address: [www.medivon.pl/deklaracje](http://www.medivon.pl/deklaracje)

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska  
Model: Medivon Cello