User manual

PL | EN | DE | FR | ES | PT | IT | CZ | NL
SV | SK | SL | EE | LT | LV | HU | RO | RU
NR | FI | EL | HR | DA | BG | TR | AR

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User Manual

The device is suitable for immediate use after purchase. Make sure the device is charged, then put it on your head and adjust the position of the massager using the rear dial to your liking. When the massager is turned off, press and briefly hold the power button to turn it on. The first - default massage mode will be launched. Press the button to cycle through the operating modes according to your preferences.

The following is a diagram of the cyclical sequence of modes and functions:

- 1) Air massage + Warming massage
- 2) Air massage + Vibration
- 3) Air massage
- 4) Air massage + infrared light

If the device is discharged, connect the device to a power source to charge the battery. Full battery will be indicated by a blue icon. The supply voltage is 5V / 1A (USB-C).

To play music, connect to your mobile device via Bluetooth. To do this, open the Bluetooth settings and search for our head massager in the list of devices. Once connected, we can control playback using the Bluetooth button. A short press plays/pauses the music, a double press plays the next song, a longer press on the Bluetooth button starts Bluetooth pairing with the devices. You can also answer phone calls using the Bluetooth button.

The massager is ideal for people who require deep relaxation, want to improve the condition of their hair or improve scalp circulation. With systematic use of massage, you will reduce the symptoms of fatigue, insomnia and chronic headaches. The device must not be used by people with injuries or who have undergone surgery in the last 3 months. Do not use the device in high temperatures or humid rooms. Always keep away from water and children. The device cannot be used by children, people with hypertension, hypotension or after severe strokes. People with skin problems should consult a doctor before use. Store in a dry and shaded place, at room temperature, out of reach of children.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after August 13, 2005.

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The manufacturer hereby declares that the Medivon Hebe Pro device complies with Directive 2014/35 / EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, Directive 2014/30 / EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65 / EU of the European Parliament and of the Council of 8 June 2011 on restriction of the use of certain hazardous substances in electrical and electronic equipment. The full text of the declaration of conformity is available at the internet address: www.medivon.pl/deklaracje

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska
Model: Medivon Hebe Pro