

<div data-bbox="188 171 265 209"> </div> <div data-bbox="188 238 561 288"> <p>Usmile Sonic Electric Toothbrush U3 User's Manual</p> </div> <div data-bbox="351 317 398 717"> </div>	<div data-bbox="625 174 728 192"> <p>Quick Guide:</p> </div> <div data-bbox="625 198 976 276"> <p>Install the brush head The usmile trademark is on the front. Please avoid installing backwards. Otherwise, it will cause the product to be damaged or not work properly!</p> </div> <div data-bbox="702 288 942 470"> </div> <div data-bbox="625 499 719 516"> <p>Clean teeth:</p> </div> <div data-bbox="625 522 951 540"> <p>1. Press the mode button to select the brushing mode.</p> </div> <div data-bbox="694 550 921 717"> </div>	<div data-bbox="1062 174 1430 212"> <p>2. Wet the brush head, squeeze the toothpaste, put it in the mouth, press the power button, and start cleaning the teeth.</p> </div> <div data-bbox="1062 218 1238 302"> </div> <div data-bbox="1062 311 1430 349"> <p>3. In the process of brushing, press the function switch to adjust the vibration intensity.</p> </div> <div data-bbox="1062 355 1165 448"> </div> <div data-bbox="1062 457 1430 541"> <p>4. Scientific brushing time is 2 minutes, which will vibrate once every 30s to remind you to switch to another zone. After brushing for 2min, it will vibrate three times to remind you and turn off automatically.</p> </div> <div data-bbox="1062 547 1148 631"> </div> <div data-bbox="1062 637 1430 719"> <p>Reminder: For those who use for the first time or who have sensitive gums, it is recommended to use the care mode first. And gradually increase the intensity to obtain a better adaptation process.</p> </div>	<div data-bbox="1499 174 1867 212"> <p>Thank you for choosing usmile sonic electric toothbrush.</p> </div> <div data-bbox="1499 218 1867 294"> <p>Be sure to read this manual carefully before using this product in order to maximize its excellent performance in tooth cleaning and care. Please keep this manual for future reference.</p> </div> <div data-bbox="1499 355 1567 375"> <p>Contents</p> </div> <div data-bbox="1499 410 1705 448"> <p>Safety Instructions Using the Product Properly</p> </div> <div data-bbox="1499 457 1713 602"> <p>Battery Status and Charging Install usmile Special Brush Head Adjusting the Brushing Mode Brushing the Teeth Properly Brushing your teeth zone Cleaning and Placement Travel with usmile</p> </div> <div data-bbox="1499 614 1696 672"> <p>Environmental Protection Product Information Service and Warranty</p> </div>	<div data-bbox="1936 174 2304 253"> <p>Safety Precautions To avoid risks of breakdown, electric shock, injuries, fires, and equipment or property damage, be sure to observe the following safety precautions.</p> </div> <div data-bbox="1936 259 2304 719"> <ol style="list-style-type: none"> 1.This product has a built-in rechargeable battery, and should not be thrown into fire or exposed to the sun, otherwise the electrolyte may leak out, and the battery may be overheated or explode. 2.Do not use this product while it is being charged. 3.Do not use this product while the charging port is not properly covered. 4.Do not refit or repair this product yourself. 5.Do not put the charging cable attached to this product in a humid environment or water. 6.Keep this product out of reach of children and infants. 7.Do not use this product for any purpose other than oral care. 8.Prevent this product from falling and impacts. 9.Do not use this product while taking a bath, shower or hot spring. 10.If you are suffering from a severe periodontal disease and being treated, or you worry about any oral discomfort, please consult your dentist before use this product. 11.Do not clean this product with any paint thinner, benzene, alcohol or any other chemical. 12.Slight bleeding may occur when you have just begun to use this product, because your gums are not used to such irritation. </div>	<div data-bbox="2373 174 2741 404"> <p>Bleeding would usually stop one or two weeks later. If bleeding lasts for more than two weeks, your gums may be abnormal, and you should suspend using this product and consult your dentist. 13.In the first few days of using this product, your oral cavity may itch, which is a normal reaction of the surface of your oral cavity and tongue because it is not used to the high-frequency vibration of the brush head. This feeling will disappear after several days of continues use. 14.If this product is damaged, rusty, cracked, etc., please stop using it immediately.</p> </div> <div data-bbox="2373 439 2501 457"> <p>Product Structure</p> </div> <div data-bbox="2416 468 2724 725"> </div>	<div data-bbox="2810 174 3161 717"> </div>	<div data-bbox="3247 174 3444 212"> <p>Using the Product Properly Battery Status and Charging</p> </div> <div data-bbox="3247 218 3615 297"> <p>This product has been charged before factory release, and is ready for use. When the battery indicator is flashing red, means battery is low, and should be charged as soon as possible.</p> </div> <div data-bbox="3247 317 3598 532"> </div> <div data-bbox="3247 582 3615 684"> <p>This product is charged with a Type-C standard charging interface that features high compatibility and fast charging. You can use a USB charging adapter or even power bank that meets the standard to charge this product, but using an inferior charging adapter may damage this product.</p> </div>	<div data-bbox="3684 174 4052 273"> <p>When charging, open the charging port cover on the bottom of the handle, and connect the charging cable attached to this product to a USB charging adapter meeting the standard for charging. And put this product in a dry environment for charging.</p> </div> <div data-bbox="3718 294 4001 387"> </div> <div data-bbox="3684 407 4052 617"> <p>During charging, the battery indicator will flash in white. When the battery is fully charged, the indicator will remain illuminated in white. When the battery is about to run out, the indicator will remain illuminated in red, and this product cannot operate, and should be charged immediately. When the battery completely has run out, the indicator will not flash in white until several minutes after the charging. After fully charging, be sure to tighten the charging port cover before using this product.</p> </div> <div data-bbox="3718 637 4001 719"> </div>
<div data-bbox="188 783 479 800"> <p>Install the usmile Special Brush Head</p> </div> <div data-bbox="188 806 561 969"> <p>When this product is used, please use a dedicated brush head for usmile. When the brush head has been used for three consecutive months, or is linting or forking, please replace it to avoid injuring your gums and teeth. When install this product, insert the brush head into the driving shaft slightly hard. And there should be a gap of 0.5mm-1mm between the brush head and the handle to ensure that the brush head can vibrate normally.</p> </div> <div data-bbox="205 998 561 1278"> </div>	<div data-bbox="625 783 848 800"> <p>Adjusting the Brushing Mode</p> </div> <div data-bbox="625 806 993 841"> <p>You can press the function switch keys continually to adjust this product to the brushing mode that you need.</p> </div> <div data-bbox="651 861 968 1324"> </div>	<div data-bbox="1062 783 1190 800"> <p>Brushing Modes:</p> </div> <div data-bbox="1062 806 1430 1179"> <p>Before the electric toothbrush is started, press the function switch key to select mode. Cleaning mode: Brush the teeth with a normal force for 2min, which will vibrate once every 30s to remind you to switch to another zone. Whitening mode: Brighten the tooth surface with polishing and powerful cleaning alternating frequency conversion work, every 30s reminder for zone change, and it will automatically turn off after 2min. Soft mode: Brush the teeth with a weaker force for 2min, which will vibrate once every 30s to remind you to switch to another zone. After brushing for 2min, it will vibrate three times to remind you finishing your brush, and turn off automatically. Massage mode: Brush your teeth with a comfortable massage frequency conversion for 2min, every 30s reminder for zone change, after brushing for 2min, it will vibrate three times to remind you finishing your brush, and turn off automatically.</p> </div> <div data-bbox="1062 1202 1268 1220"> <p>Brushing Gear Description:</p> </div> <div data-bbox="1062 1225 1430 1304"> <p>Three levels of intensity, weak/medium/strong, are set to satisfy soft, moderate, strong, etc. different user needs. And in the process of brushing, press the function switch key to adjust the vibration intensity.</p> </div>	<div data-bbox="1499 783 1713 800"> <p>Brushing the Teeth Properly</p> </div> <div data-bbox="1499 806 1867 946"> <p>When using the usmile sonic electric toothbrush to brush your teeth, do not press the head hard to keep the bristles in contact with the tooth surface, and move slowly laterally for cleaning. It is recommended that you keep the bristles and tooth surfaces as shown in the figure when brushing the outer and inner sides of your teeth angle of inclination to obtain better</p> </div> <div data-bbox="1542 966 1842 1237"> </div>	<div data-bbox="1936 783 2047 800"> <p>Start Brushing</p> </div> <div data-bbox="1936 806 2304 905"> <p>It is recommended that you divide your teeth into four zones, and brush each zone for 30s. After this product is started, it will give a frequency change prompt every 30s. When your teeth have been brushed for 2min, three consecutive frequency changes will be made as the prompt.</p> </div> <div data-bbox="1936 966 2304 1278"> </div>	<div data-bbox="2373 783 2553 800"> <p>Cleaning and Placement</p> </div> <div data-bbox="2373 806 2741 905"> <p>After brushing, please press the switch button on the handle to stop the product. Clean the brush head and handle with water. After drying with a towel, stand it upright on a dry surface, or put it in a mouthwash cup.</p> </div> <div data-bbox="2373 934 2741 1301"> </div>	<div data-bbox="2810 783 2947 800"> <p>Travel with usmile</p> </div> <div data-bbox="2810 806 3178 905"> <p>When you're going to travel with usmile sonic electric toothbrush, please remove the brush head and place it in the dust-proof tube of the brush head provided with the product. Please put the handle in your luggage properly to prevent any nearby object from pressing the button on it.</p> </div> <div data-bbox="2810 926 3178 1310"> </div>	<div data-bbox="3247 783 3444 800"> <p>Environmental Protection</p> </div> <div data-bbox="3247 806 3615 928"> <p>usmile sonic electric toothbrush has a built-in rechargeable lithium battery, please do not throw it away mixed with general household garbage. Send it to your local battery recycling facility or send it back to the usmile After-Sales Service Department. In this way, you're contributing to environmental protection.</p> </div> <div data-bbox="3247 952 3401 969"> <p>Product Information</p> </div> <div data-bbox="3247 975 3615 1118"> <p>Product: usmile Sonic Electric Toothbrush U3 Model: U3 Rated voltage: DC≐ 3.6V; Rated power: 1.2W Manufacturer: Guangzhou Haoxingwan Technology Co., Ltd. Address: Room718, No.9 Linhe West Road, Tianhe District, Guangzhou Origin: Made in China</p> </div> <div data-bbox="3247 1255 3615 1319"> <p>usmile is the registered trademark for oral care products exported by Guangzhou Haoxingwan Technology Co., Ltd. All rights reserved.</p> </div>	<div data-bbox="3684 783 3838 800"> <p>Service and Warranty</p> </div> <div data-bbox="3684 806 4052 1010"> <p>For any product fault or damage as conformed by our After-sales Service Dept. have occurred during normal operation within one year after purchase, you will be entitled to free warranty. Free warranty services exclude the brush head, damages arising from the use of unauthorized spare parts, damages arising from misuse, abuse, omission, modification or unauthorized repair, normal wear and tear, including fissures, scratches, wears and discolorations, and damages arising from force majeure.</p> </div> <div data-bbox="3684 1033 3915 1074"> <p>When you need to contact us, Please email us via info@mysmile.com</p> </div>