# Ritual peeler

## Reclar Plus

•

•

•

•

•

RECLAR

**Thank you for trusting the Ritual Reclar Plus peeler.  
Beauty is now in your hands.**

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
Ultrasound water peeling effectively removes sebum, old skin cells and is great for deep cleaning of pores. A frequency of 28–30kHz ensures the maximum cleansing effect while being gentle towards the skin.

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# Cleansing ritual

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**Result:** clean skin free of old skin cells and impurities, ready to absorb active substances from your favourite cosmetic product.

**You will need:** a quiet place with a mirror, removed make-up, Ritual Reclar Plus peeler, 5 minutes of your time.

**Method:** Turn on the Peeling feature and move the peeler from your temples towards your nose and up to the centre of your forehead, then from your mouth towards your ears and from the tip of your nose upwards. Continuously moisturise your skin, for example using a wet cotton pad. Booster Peeling function can be used for extra cleansing of the T-zone or other areas with severely clogged pores. To activate this function, press the button when the device is off and hold it for 5 seconds, until you hear a “beep” and the lights are cross-lit. The booster peeling mode lasts for 1 minute.

We recommend performing the cleansing ritual 1 – 2 times per week, depending on your skin type.

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
# Healing ritual

## Léčivý rituál

**Result:** healed and calm skin, removed scars and spots, and improved effectiveness of your cosmetics as a bonus.

**You will need:** your favourite cream, serum or mask, Reclar Ritual Peeler Plus, a couple of minutes.

**Method:** Apply your favourite cream, mask or serum and turn on the G-Blue feature. Turn the front side of the peeler towards yourself and slide the tip from the centre of your face to the sides. Keep touching the sides to keep the blue light on.

We recommend performing the healing ritual regularly, a couple of minutes a day are sufficient.

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
The blue light heals acne and sterilises tissues to prevent more inflammation. It also aids the healing of minor wounds and scars. Galvanisation improves the effects of your cosmetic products and helps keeping the skin flexible.

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# Rejuvenating ritual

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**Result:** smooth and bright skin, tightened pores and activated collagen production.

**You will need:** your favourite cream, serum or mask, Reclar Ritual Peeler Plus, a couple of minutes.

**Method:** Apply your favourite cream, mask or serum and turn on the G-Red feature. Turn the front side of the peeler towards yourself and slide the tip from the centre of your face to the sides. Keep touching the sides of the peeler to keep the rejuvenating red light on.

We recommend performing the rejuvenating ritual regularly, a couple of minutes a day are sufficient.

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
The red light activates collagen, reduces wrinkles and pigment stains and closes pores. Galvanisation improves the effects of your cosmetic products and helps keeping the skin flexible.

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# Quick start

## Rychlý start

### 1) Charge the device fully using the supplied cable.


Do not apply excessive force when plugging the cable into the device. Charging is confirmed by a red light. Blue light means the device is fully charged.

### 2) Take off the blade cover.

### 3) Turn the device on by holding the control button.

### 4) Pick your desired feature:

 Peeling – press the button once, for Booster peeling mode press the button when the device is off and hold it for 5 seconds, until you hear a “beep” and the lights are cross-lit.

 G-Blue – press the button twice

 G-Red – press the button three times

Note: To activate the light therapy you need to keep touching the electrodes – side metal bars – throughout the procedure. If you are applying G-Blue or G-Red onto a different person, physical contact with that person is necessary to close the circuit.

### 5) Turn the device off by holding the control button.

Clean the blade gently with an attached cloth after every use or use a disinfectant.

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
PEELING FEATURE/FUNKCE PEELING

# Ultrasound skin cleansing

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Apply once or twice per week, depending on the skin type. The time for this cleansing is set to a maximum of 5 minutes. The booster peeling mode lasts for 1 minute. Afterwards, the feature automatically ends. We do not recommend exceeding 5 minutes/ 1 minute per one application.

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
**Step 1:** Remove make-up and wash your skin. Your skin should remain moist throughout the procedure, so we recommend having a cotton pad and some water ready so you can moisturise throughout.

**Step 2:** Turn the device on by holding the button. To activate the Peeling feature, press the button once to turn on the indicator light of the desired feature. If your device is in standby mode, you can just press the button once. To activate Booster peeling mode, press the button when the device is off and hold it for 5 seconds, until you hear a “beep” and the lights are cross-lit.

**Step 3:** Keep the peeler’s back side turned towards you and gently place the blade onto moist skin. Move the peeler from your temples towards your nose and up to the centre of your forehead, then from your ears towards your mouth and from the tip of your nose upwards. Keep your face moisturised.

Peeling can dehydrate your skin so we recommend applying your favourite hydrating skincare product afterwards.

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
G-BLUE FEATURE/FUNKCE G-BLUE


# Galvanisation with blue light therapy

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Suitable for daily use. The time for one therapy is set to a maximum of 5 minutes. Afterwards, the feature automatically ends. We do not recommend exceeding this time.

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G-Blue Feature  
Funkce G-Blue

**Step 1:** Apply your favourite skincare product (cream, mask, serum) on a clean skin. Your skin must be sufficiently moist.


**Step 2:** Turn the device on by holding the button. To activate the G-Blue feature, press the button twice to turn on the indicator light of the desired feature. If your device is in standby mode, you can just press the button twice.

**Step 3:** Turn the front side of the peeler to yourself and slide the tip from the centre of your face to the sides. Keep touching the sides of the peeler to keep the healing blue light on.

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G-RED FEATURE/FUNKCE G-RED


# Galvanisation with red light therapy

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G-Red Feature  
Funkce G-Red

**Step 1:** Apply your favourite skincare product (cream, mask, serum) on a clean skin. Your skin must be sufficiently moist.



**Step 2:** Turn the device on by holding the button. To activate the G-Red feature, press the button three times to turn on the indicator light of the desired feature. If your device is in standby mode, you can just press the button three times.

**Step 3:** Turn the front side of the peeler to yourself and slide the tip from the centre of your face to the sides. Keep touching the sides of the peeler to keep the rejuvenating red light on.

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# FAQ/Časté dotazy

## Can I use the peeler while showering?

The device is not waterproof. The device can be damaged if it gets wet.

## How to properly clean the blade?

You can clean the blade with our Reclar cloth or using a soft cotton fabric.

## How to disinfect the device?

Use any alcohol-based disinfectant.

## Is it normal if the blade moves?

The blade is deliberately not fixed in place so that the device can function properly while vibrating. Do not pull at the blade with force, just wipe it gently.

## Where should I put the peeler after using it?

Keep the peeler in a dry place. We do not recommend keeping it in a bathroom or other damp room.

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### **Can I use the peeler while pregnant?**

You can use the Peeling feature. Both G-Blue and G-Red features are not recommended during pregnancy.

### **Is the light harmful to the eyes?**

Keep to the instructions to prevent direct shining of the LED light to your eyes.

### **What do I do when the lights are off?**

To close the galvanising circuit and for the light to come on you need to hold the side electrodes (metal bars) while using the G-Blue and G-Red feature.

### **Can I use the peeler on someone else?**

If you are applying G-Blue or G-Red onto another person, physical contact with that person is necessary to close the galvanising circuit. The Peeling feature does not require it.

### **Is it necessary to use special cosmetics products?**

You can use the Reclar Peeler with any skincare products you want.

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### **Can I use the individual features repeatedly?**

We recommend keeping to the times given for each feature and not exceeding it.

### **How do I charge the peeler?**

Use the supplied cable along with any USB charger. Plug the charging cable in gently. We recommend charging the device fully every time.

### **How do I know the battery is low?**

The device provides sound signaling during use when the battery is low, along with the lights going on.

### **When is the device fully charged?**

Charging is indicated by a red light. When the device is fully charged a blue light will appear.

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### **Can I use the peeler while it is charging?**

The device is not usable while charging.

### **The device is not charging, what do I do?**

If your device is not charging do not try to push the cable into the device. The charging cable may be damaged. In that case, use a different compatible cable.

### **For more info on how to use the device, check out the video manuals at:**

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[reclar.com/instructions](https://reclar.com/instructions)

reclar.com