



## Quick Start Guide

# pmd | kiss

### ANTI-AGING LIP ENHANCING TREATMENT

CAUTION: Please read this manual carefully before using

---

#### GETTING STARTED:

Mode 1: Standard Speed for a quick treatment

Mode 2: Low Speed for a more intense targeted treatment

#### CUSTOMIZE YOUR TREATMENT:

**kissfoliator™ Kisser:** The kissfoliator removes dead skin cells using anti-bacterial aluminum oxide crystals for smooth kissable lips.

**Select Your Kisser:** Two sizes of Kisser tips are available. Test both tips and decide which one is the best fit for your lips.

**hydrakiss™ Lip Mask:** For added lip hydration, and to enhance the anti-aging results of the treatment

#### HOW-TO USE

Instructions:

1. Insert Kissfoliator tip into PMD Kiss
2. Power on PMD Kiss. Start with clean dry lips, gently glide PMD Kiss across the lips while device pulses. For a more intense treatment, try applying pressure or scrubbing. Use as desired on dry, chapped lips.
3. Apply Smart Lip Plumping Serum to clean, dry lips.
4. Swap PMD Kiss Tips. While Smart Lip Plumping Serum takes effect, remove Kissfoliator tip and insert Kisser tip.
5. Power on PMD Kiss. Place Kisser tip firmly and directly onto lip creating a tight seal for suction. Allow PMD Kiss to pulse and glide across your lips until desired result is achieved.
6. Apply hydrakiss lip sheet mask to lips for a minimum of 5 minutes to enhance anti-aging results.

#### Charging:

Charge device for 2-3 hours before first use.

#### INCLUDES:

PMD Kiss  
Charging Cord  
Kiss Serum  
Kissfoliator™  
2 Kisser Tips  
2 hydrakiss™ Lip Sheets



The PMD Kiss comes with a One-Year Limited Manufacturer Warranty.

## FAQ's

### What are the long term results?

When you use the PMD Kiss and Smart Lip Plumping Serum daily for 4-6 weeks, your lips will stay full from increased collagen production. The appearance of lip lines will begin to fade and lips will appear healthy and youthful. Results vary with every user.

### What are the short-term results?

You will see immediate plumping, and natural lip color will be revitalized in minutes. Lips will appear fuller and more rosy for 2-4 hours. If you continue to use the PMD Kiss daily, results will increase and the appearance of lips will remain full! The PMD Kiss is safe to use as often as desired.

### How long does the plumping last?

You will see an instant plump immediately after using the PMD Kiss and Smart Lip Plumping Serum. Your lips will stay plump for several hours. With continued daily use, lips will stay full and results will be significantly more visible.

### Can I use PMD Kiss if I have lip injections?

We suggest consulting your doctor if you have injections, or waiting for them to dissolve before using the PMD Kiss.

#### Caution and Warnings:

Use only as directed. Not following instructions may lead to rash, redness, itching or even injury. If this occurs, discontinue use and consult a medical professional. Product is intended for external use only. Avoid sensitive areas, including any area with broken, irritated, or itching skin. Keep product out of reach of children. PMD is not responsible for malfunction or injury if the product is not used as intended. Not intended for clinical or commercial use. The effects of the PMD Kiss are not always immediately visible and results may vary. If you have any questions regarding application, please reach out to our experts at 888-445-4532.

---

## #BRILLIANTCONFIDENCE

To learn more visit: [eu.pmdbeauty.com/kiss](http://eu.pmdbeauty.com/kiss)

### FOLLOW US

Youtube  
@pmd

Instagram  
@pmdbeauty

Pinterest  
@pmdbeauty