

Electric toothbrush

Contents of the package:

1. Toothbrush - electronic device
2. spare toothbrush head - replacement
3. USB to USB-C charging cable

First of all, thank you for choosing our product. Our company provides a reliable electric toothbrush to help you brush your teeth and keep your oral health.

The electric toothbrush can slow down the accumulation of plaque, effectively remove stains from teeth, tea and coffee, and play a supportive role in protecting your mouth and improving the condition of your gums.

Important safety information:

1. At the beginning of usage of our product, if you feel any irritation or numbness, please do not worry, it is normal due to the high-speed movements of the toothbrush bristles. This sensation will gradually reduce and disappear within 2 days.
2. Please do not disassemble, repair or modify the product.
3. Please do not hit or drop the product, if the product is damaged, do not use it.
4. Before connecting to the power supply, please check whether the voltage is the same with the voltage of the toothbrush.
5. This electric toothbrush is only compatible with our specified sonic toothbrush head. Do not attempt to connect other second heads, the safety and stability of the head attachment is not guaranteed.
6. If you have had oral surgery in the last 2 months or your gums bleed after using the electric toothbrush in the period for more than 1 week, please consult using it with your professional dentist.
7. If you have a heart rhythm controller or other electronic devices installed in your chest, please consult usage with your dentist before use.
8. This product is for cleaning teeth, gums and tongue only, do not use it for other purposes.
9. Place the product out of reach of children, the toothbrush cannot be used as a toy, there is a risk of permanently damaging the device.
10. This product is for indoor use only. Do not use it in direct sunlight, do not leave it under direct sunlight for long periods of time. It should be placed and stored in a dry and ventilated place at room temperature.

For users with a sensitive oral cavity!

1. Persons who have had recently oral surgery related to periodontology.
2. Persons suffering from other non-standard oral conditions such as insertion dentures, or other non-fixed dentures
3. Persons with sensitive and fragile oral cavity with frequent bleeding gums.

If you are one of these 3 points, please do not use this product or consult toothbrush usage with your dental professional before using it. This product is not suitable for persons with disabilities, persons with mental disorders and children under 6 years of age unless using this product under the supervision and guidance of an adult.

Care of electronics in the equipment

The device has sensitive electronics in it. Disconnect the power brush when the product is not charging or when you are cleaning it. Make sure to avoid dust and dirt on the USB cable end, otherwise it will affect charging. This means that charging may not take place or will not be continuous and will drop out as dust will be introduced when touching the USB connectors.

First usage of the toothbrush

Switching on: a short press of the button switches the brush on and the first time you switch it on, the gentle mode is the default.

Switch mode: Within 5 seconds of switching on, press the power button briefly and you can switch working modes in this way.

Custom mode: In custom mode, long press the button without release to adjust the vibration intensity of the toothbrush. Select the appropriate intensity for you and immediately release the button. Then the current intensity will be adjusted.

Smart mode memory: After using any mode for more than 5 seconds, the toothbrush will intelligently record that mode and store it in memory. The next time you switch it on, it will automatically set that mode.

Power Off: Press the button briefly after use for 5 seconds, or the brush will turn off automatically after 2 minutes of use.

Charging

1. If the toothbrush is discharged after unboxing, please charge it. After it is fully charged, you can charge it at any time. You need to wait for about 3-4 hours to fully charge.
2. When the battery has less than 10% capacity, the charging indicator will turn red to alert you to the need to charge the toothbrush immediately.

3. This product is charged using a USB cable.
4. Charging warning: green light is on and pulsing while charging. When fully charged, the green light is on continuously. You need to wait for approximately 3-4 hours for a full charge.
5. During charging, DO NOT press the power button. After disconnecting the USB cable, you can use it.
6. When the toothbrush is in use and you insert the USB cable, the toothbrush will stop operating and start charging.
7. When fully charged, disconnect the USB cable and reinstall the waterproof cover to prevent the toothbrush from coming into contact with water.

Notice

1. If the toothbrush does not turn on or does not work properly, return it to the manufacturer for repair or replacement and do not attempt to disassemble it yourself. This will prevent accidents.
2. If the toothbrush will not be used for a long time, please charge it at least once every three months. This will prevent the battery capacity from being destroyed.
3. Do not charge the toothbrush with a USB cable different from the standard USB cable of this product to prevent the battery from overheating.
4. It is normal for the battery to become warm during charging. If the temperature does not exceed 40 degrees Celsius.
5. When the low battery indicator comes on, please charge the brush as soon as possible to avoid damaging the battery and affecting its vibrancy.
6. While charging, please place the brush out of the reach of children to prevent accidents.

Dentist's recommendation for the correct way to brush teeth

Cleaning time

Oral hygiene experts recommend: use an electric toothbrush to brush your teeth for two minutes, not less.

During brushing

A short pulse will alert you to change the brushing site of your mouth. The pulse will sound every 30 seconds to clean the four parts of your mouth for optimal results.

Brush head replacement

Dentists recommend replacing the brush head every 2-3 months for optimal oral hygiene results. It should only be used with a compatible (original) hello coco toothbrush head.

Cleaning

Wipe the brush only with a soft damp cloth, do not use rough and sharp objects, brushes or sponges. Then allow the brush to dry naturally. Please do not immerse the brush in water for a long time.

Brush head

After use, clean the brush head, rinse it properly and let it dry naturally. Please do not tear or pull the bristles of the brush head, as they can easily become deformed. It is recommended to replace the brush head every 2-3 months, after a period of use the brush head will become rougher and this will reduce the cleaning effect.