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Symbol  znamená UPOZORNENIE

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Massager eye device**eta 3353****INSTRUCTIONS FOR USE**

Dear customer, thank you for purchasing our product. Please read the operating instructions carefully before putting the appliance into operation and keep these instructions including the warranty, the receipt and, if possible, the box with the internal packing.

I. SAFETY INSTRUCTIONS

- Consider the instructions for use as a part of the appliance and pass them on to any other user of the appliance.
- The appliance is only designated to be used in household and for similar purposes. It is not designed to be used in beauty parlours or for any other commercial utilisation!
- Children and persons with limited physical or mental capabilities or lack of experience may only use this appliance if they are supervised or have been instructed of the use of the appliance in a safe manner and understand possible risks. Unattended children must not perform cleaning and maintenance intended to be performed by users. Children must not play with the appliance.
- Never immerse the appliance in water!
- Children must be kept away from the appliance and its adapter.
- The appliance has no charging source of its own, it is only designed for charging via USB (an adapter with USB output, PC with USB output with output voltage of 5.0 V, and output current of at least 500 mA). Before making any connection, please pay attention to instructions of the manufacturer of the respective power source.
- For additional protection we recommend installing a protective switch (RCD) in the electrical circuitry of the bathroom, with the nominal actuating current not exceeding 30 mA. Ask an inspecting technician or an electrician for advice.
- The appliance and the connected adapter must not be used at places where it could fall into bath tub, wash basin or swimming pool. If the appliance fell into water anyway, do not take it out! First of all, unplug the adapter from the socket and only then take

the appliance out. Then bring the appliance to a specialised service centre to check its safety and proper functionality.

- If the appliance is used in a bathroom, it must be disconnected from the mains after use by unplugging the power cord from the wall socket, because the appliance poses a risk if near to water, even when switched off.
- Always unplug the appliance from charger if you leave it unattended and cleaning.
- The mains adapter supplied with the appliance is only designated for this appliance. Do not use it for other purposes. Only use the attached mains adapter for this appliance. Do not use any other adapter for charging (for example an adapter from a different appliance).
- Never plug the adapter in the electrical socket or unplug it with wet hands.
- Never immerse the adapter into water or other liquids (not even partly) and do not wash it in running water!
- The appliance gets hot when charging, which is absolutely normal.
- Charge the battery at a temperature close to the room temperature.
- Do not dispose the battery into fire. There is a risk of explosion!
- Never charge the appliance in direct sunlight or close to heat sources, and if the temperature is below 10 °C or above 40 °C, adversely affects the accumulator service life.
- Prevent hanging the power adapter cord loosely over desktop edge where children may reach it.
- The power cord of the adapter must not be damaged with sharp or hot objects or open fire, and must not be immersed into water or bent over sharp edges.
- Before use make sure to remove glasses, contact lenses, and/or jewellery.
- When handling the appliance, please proceed so that no personal injury or damage to the appliance will occur.
- Never use the appliance if it does not correctly operate, fell on the ground or suffered any damage.
- Protect the appliance from dust, chemicals, excessive humidity, toiletries, aggressive liquid cosmetic preparations and intensive temperature changes, and do not expose it to direct sunlight.
- Do not insert the appliance in any bodily orifices.
- Never use the appliance for other than the intended purposes as described in these instructions for use!
- **WARNING:** There is a risk of injury in case of incorrect use (not in accordance with the manual) of the appliance.
- The producer does not undertake any responsibility for damage caused by improper use of the appliance (e.g. injury etc.) and its guarantee does not cover the appliance in the case of failure to comply with the safety warning above.

SPECIAL SAFETY INSTRUCTIONS

- This appliance has a heated surface. Persons insensitive to heat must be careful while using this appliance.
- The following persons should not use this device: pregnant women or while nursing or during menstruation; persons allergies; persons with facial nerve disorders; persons with severe skin diseases on the face, person with eyelash implants.
- When using the red-light therapy, close the eyes and avoid looking directly at the light source.
- Make sure to use the appliance only on face, not on the eyelids to avoid injury to the eye.
- Do not use the appliance:
 - if your while wearing contact lenses.
 - if your skin is oversensitive, damaged or exematic

(e.g. in the case of fresh scars, small injuries, burns, erythema, swollen cheeks, skin diseases...). To avoid dissemination of micro inflammations, do not use the massager roller including cases when your skin suffers acne. We do not recommend using the massager roller on injection-treated skin provided that in the given area you have undergone surgical intervention or plastic facial rejuvenation (such as e.g. Botox, facial fillers, thread lifts etc.).

- On areas which may present a health hazard (varicose veins, moles, birthmarks).
- If you have any tattoo on your face or permanent makeup.
- If you feel pain or experience any unpleasant feeling during the using, or if any other health problem appears, immediately switch off the appliance, and consult its further use with your practitioner or dermatologist.
- The appliance is not intended to cure any skin disease symptoms. Please consult use of the appliance with your doctor, where necessary.
- Do not use the appliance on children, incapacitated persons or persons unable to express their feelings.

Note

ETA a.s. is not a licensed provider of health and cosmetic care and does not have expertise in the field of treatment, diagnostics of the health condition and related effects or suitability for the use of appliances that may affect the health condition in any way. During any exercise, massage therapy or cosmetic intervention, there is a theoretical possibility of physical injury or, for example, worsening of the symptoms of existing diseases. By using this appliance, you agree that you use it at your own risk and assume all risks of injury or deterioration. You also agree to waive any claims (whether known or unknown) against ETA a.s. arising from improper use of this appliance or through negligence. Before using any appliances that may affect your health, first consult your doctor about the suitability of their use.

II. DESCRIPTION OF FEATURES (PIC. 1)

A – Massager body

- A1 – ON/OFF switch + mode switch 
- A2 – indicator light 1 (LIGHT)
- A3 – indicator light 2 (SONIC)
- A4 – indicator light 3 (HOT)
- A5 – massage head
- A6 – light
- A7 – power USB slot
- A8 – protective cover

B – USB cable: (Use AC power adapter with USB or mobile power)

Mode on/switch button control A1

- function 1: indicator 3 is on = the head is heating up,
- function 2: indicator 2 is on = head is vibrating,
- function 3: indicator 1 is on = red light is on,
- function 4: indicators 2 and 3 are on = the head is vibrating and heating up,
- function 5: indicators 1 and 3 are on = red light is on and the head is heating up,
- function 6: indicators 1 and 2 are on = red light is on and the head is vibrating, head,
- function 7: indicators 1 and 2 and 3 are on = red light is on, the head is vibrating and warming up.

III. INSTRUCTIONS FOR USE

Remove all packaging materials and take out the appliance.

CHARGING THE APPLIANCE (Fig. 2).

- 1) Make sure the appliance is off. While charging, the appliance cannot be in use. If the appliance is turned on and you start charging, the appliance shall switch off.
- 2) Put the provided charge USB cable **B** into the charge slot **A7** of brush body, connect the other charger compatible AC wall adapter or mobile power.
- 3) The operation of the battery charging is signaled by the gradual switching on and off of the indicators 1, 2, 3.
- 4) When the battery is fully charged, the indicator stays on when charged. Unplug the adaptor from the wall socket first and only then unplug the connector from the case.
- 5) For a full charge, the battery needs to be charged for approximately 2 hours. Depending on the using mode used, a fully charged battery will provide approximately 30 minutes of operation.

Caution

- Flashing indicators 1, 2, 3 during appliance operation indicate that the battery is flat.
- Do not exceed the recommended time of charging during regular charging!

Apply a layer of eye cream to the skin around the eyes (do not apply any substances to the skin if you use the light therapy mode). Press and hold (3 sec.) the button **A1** to turn on the appliance (a beep will sound). Choose the mode you want to use depending on the needs of your skin. Gently press the massage head **A5** against the skin, slowly move it or apply slight pressure from the inner to the outer corner of the eye. The recommended time of use is 1.5 minutes for each eye. Use daily for better results. After use, switch off the appliance and continue normal skin care.

Notes

- The appliance always switches on in function 1 setting.
- The appliance switches off automatically after 5 minutes of operation. If you change the function during the operating mode, the 5 minute interval will start from the beginning.

IV. MAINTENANCE

For maximum hygiene, clean the massage head **A5** after each use (wipe the massage head surface with a wet soft cloth). Do not use any abrasive or aggressive cleaning agents! Do not immerse the massager body **A** in water. Store the appliance at a safe, dry and dust-free place, out of reach of children and incapacitated persons.

Guarantee of 24 months does not apply to decrease in the accumulator's effective capacity caused by its use or age. Effective capacity gradually decreases depending on the method of use of the accumulator.

**Expert service must do the large maintenance or to intervene into inside parts of the appliance!
Non-observance of the instructions of the manufacturer cancels the right for guarantee repair!**

V. ENVIRONMENT

If the dimensions allow, there are labels of materials used for production of packing, components and accessories as well as their recycling on all parts. The symbols specified on the product or in the accompanying documentation mean that the used

electrical or electronic products must not be disposed of together with municipal waste. For proper disposal, hand them over at special collection places where they will be accepted free of charge. Suitable disposal of the product can help to maintain valuable natural resources and to prevent possible negative impacts on the environment and human health, which could be the possible consequences of improper waste disposal. Ask for more details at the local authorities or at a collection site. Fines can be imposed for improper disposal of this type of waste in agreement with the national regulations. Dispose of the battery safely by means of collection networks intended for this purpose. Never dispose of the batteries by burning!

Removing the battery

Remove the battery from the appliance only when completely discharged (i.e. leave the appliance running until it stops) and after the end of the appliance's service life! Open the pedicure set. Gradually disconnect and remove the batterie. The battery should not be removed and is not accessible during normal operation!

VI. TECHNICAL DATE

Protection class of the appliance III.

The manufacturer reserves the right to change the technical specifications for the respective models. Expert service must do the large maintenance or to intervene into inside parts of the appliance! Non-observance of the instructions of the manufacturer cancels the right for guarantee repair!

WARNINGS AND SYMBOLS USED ON THE APPLIANCE, PACKAGING OR IN THE INSTRUCTIONS MANUAL: *HOUSEHOLD USE ONLY.*

TO AVOID DANGER OF SUFFOCATION, KEEP THIS PLASTIC BAG AWAY FROM BABIES AND CHILDREN. DO NOT USE THIS BAG IN

CRIBS, BEDS, CARRIAGES OR PLAYPENS. THIS BAG IS NOT A TOY.

The  symbol indicates a WARNING.

- Please read the instructions manual

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© DATE 26/5/2022