# UBrush!

Smile, because UBrush!

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**UBrush!** Enterprises B.V.

## **EN Instruction for use**

1. Install the UBrush! and place the battery.
2. Choose the right brush size which you want to use between your teeth.
3. Press the LED light button (bottom button in gray section).
4. Look in the supplied mirror when using the UBrush!, so you can see what you are doing.
5. Place the small brush into the space between the teeth and press the upper vibration button (the electric interdental brush makes about 167 vibrations per second) and make five back and forth movements. Because of the vibrations the gums are massaged and more plaque is absorbed. You can see this on the small brush. In the meantime you can clean the brush with water or a wipe if there is plaque or blood on the brush. Your gums may bleed because it is inflamed. It is important to use the brush 1 to 2 times a day to complete your oral hygiene.
6. After you finish with the UBrush! it is wise to clean the brush and dry it with a wipe. When you don't do this bacteria will multiply. After cleaning place the cap back on the brush for a second use.
7. Replace the brush when the brush is getting thinner.

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## UBrush! Enterprises B.V.