

USER MANUAL

BeautyRelax BR-2400 Neck massager

SHIATSUSOFT



CONTENT:

NOTICE	4
IMPORTANT NOTICE	4
PACKAGING CONTENT	5
CONTRAINDICATIONS	5
INTRODUCTION OF THE DEVICE	6
MAIN ADVANTAGES	6
TECHNICAL PARAMETERS	6
HOW TO USE	7
BATTERY CHARGING	8
CLEANING AND MAINTENANCE	8
DISPOSAL OF USED ELECTRONIC EQUIPMENT	8
CERTIFICATION	8
IMPORTER	8
ES DECLARATION OF CONFORMITY	9

Dear customer, thank you for your trust in purchasing our product. Please read all the detailed information in the instruction manual carefully before using it for the first time. Please pay particular attention to the safety warnings. Please use the device exactly as instructed. All of these steps will help to maximize the use of all functions of the instrument. In this manual we will try to give you the best possible idea of the technology used and how to use the instrument to best advantage.

If you encounter any uncertainties while reading the instructions, please do not hesitate to contact us at the following email address: mh@santadria.cz

Keep the instructions for use carefully for later reference.

NOTICE

- Do not attempt to disassemble or repair the unit yourself in any way
- store and use the device in a dry, closed environment
- (it is forbidden to use the device in the bathroom)
- it is forbidden to use the device if you have had surgery on the cervical spine or any other part of the body for which the device is designed
- it is forbidden to use the device if you have serious chronic health problems with the cervical spine and if your neck area is very sensitive
- do not use the device for too long; 5 to 10 minutes a day is enough for pleasant relaxation and muscle relaxation; use for longer than the specified time could cause headaches (in case of application to the cervical spine)
- if you are unsure, consult your GP as a last resort; the device does not require any major mechanical maintenance
- disconnect the unit from the power adapter before cleaning the controller; clean the surface of the controller only with a clean soft cloth
- Do not use aggressive cleaning agents; the surface of the controller could be damaged; do not leave the packaging material lying around; it could become a dangerous toy for children
- Handle the product with care; do not subject it to shocks, knocks or impacts; permissible ambient temperature: 5°C to +50°C
- damage caused by neglect of these safety warnings and instructions for use will void any claim to the warranty of the appliance; we shall not be liable for consequential damage caused by breach of these conditions
- we do not accept any liability for personal injury caused by improper or incorrect handling of the product; in such cases, the product warranty is void



IMPORTANT NOTICE

It is strictly FORBIDDEN to use the massager while driving in the car. In the event of an accident, you run the risk of injuring your spine or any other area that would be massaged by the massager.

PACKAGING CONTENT

- BeautyRelax BR-2400 SHIATSUSOFT massage device
- USB charging cable
- User manual

CONTRAINDICATIONS

Under no circumstances should the device be used by persons affected by any of these conditions:

- Epilepsy, viral disease
- Surgery of the body part for which the device is intended
- In a state of high fever
- Severe skin disease at the site of use of the device
- Heart disease, hypertension
- Patients who have undergone surgery with transplantation
- Persons who require constant supervision may use the device only after consulting a physician



INTRODUCTION OF THE DEVICE

BeautyRelax ShiatsuSoft Massager for instant muscle relaxation in the cervical spine. The unique design of the massage mechanism faithfully follows the movements of the masseur's fingers and is complemented by the infrared light heating function. Two areas of moving massage mechanism with a total of 6 heads. Beneficial heating function. Two directions of massage movement. Timer. Wireless use. Built-in Li-ion battery. The device is very quiet (0-20 db) and very light (weighs 53 g, equivalent to the weight of an iPhone). Fully recharges in 1 hour. The cushion is inflatable and therefore practically storable.



MAIN ADVANTAGES

- Immediate relaxation of stiff muscles
- Ideal for use in the neck area
- Highly effective Shiatsu massage simulates the hand movements of a masseur
- 6 massage heads
- Two directions of movement of the massage heads
- Easy one-button operation
- 2 heating levels (low, high)
- 3 speed levels
- Timer: automatic shutdown after 15 minutes
- Wireless use thanks to built-in rechargeable battery
- Highly breathable and pleasant to the touch material (nylon)
- Charging: built-in Li-ion 1400mAh battery
- Charging voltage: 5V
- Power: 5W
- Dimensions: 20x18x6cm
- Weight: 53g

TECHNICAL PARAMETERS

- Charging: built-in Li-ion 1400mAh battery
- Power: USB cable (included)
- Dimensions: 20x18x6cm
- Weight: 53g

HOW TO USE

- 1.** Always ensure that the battery is sufficiently charged. Only use the cable included in the package for charging. Using any other type of adapter may cause fatal damage to the battery and void the warranty.
- 2.** To switch on the device, long press the **ON/OFF** button. The device has a compressor in it that automatically inflates the collar. Once the collar is inflated, the sensor stops pumping air and starts the movement of the massage heads.
- 3.** With a short press of the **ON/OFF** button you can adjust the speed of the massage. A total of 3 massage speeds are available.
- 4.** With the help of the handles you can regulate the force that the massage head exerts on your neck.
- 5.** Use the HEAT button to turn the collar heating on or off. Short press - low level of heat exposure. Long press - higher level of heat.
- 6.** The device has a timer. It switches off after a few minutes. The recommended massage duration is 10 to 15 minutes. Do not overdo it. Also, do not put too much pressure on the neck with the handles.
- 7.** You can switch off the device at any time by long pressing the **ON/OFF** button.
- 8.** When you are finished using the device, you can simply blow out the collar. Just open the valve on the side of the device. It must be carefully closed before further use to prevent air leakage.

BATTERY CHARGING

The device has a built-in rechargeable battery. A sufficiently charged battery enables the continuous operation of all functions of the device. The battery status is indicated by the battery symbol on the display.

Charging procedure:

1. When the battery is low, the battery symbol on the display will flash.
2. Use the included USB charging cable for charging.
3. Connect the other end to the socket on the side of the collar.
4. The progress of charging is indicated by a flashing red light.
5. Once the charging process is complete, the flashing of the indicator light will stop.
6. The length of charge depends on the state of discharge and the overall condition of the battery.

Please note: Only use the included adapter to charge. Using a different type of adapter may cause fatal damage to the battery and void the warranty.

CLEANING AND MAINTENANCE

Do not use any aggressive cleaning agents. Ideally, use a slightly damp cloth for cleaning. Moisture must NEVER penetrate directly into the device. A cloth or dry wipe is ideal for cleaning. Under no circumstances use any sharp objects that could damage the surface of the device during maintenance. The device is not waterproof.

DISPOSAL OF USED ELECTRONIC EQUIPMENT



This symbol on the products or in accompanying documents means that used electrical and electronic products must not be added to normal municipal waste. For proper disposal, refurbishment and recycling, take these products to the designated collection points. Alternatively, in some European Union or other European countries, you can return your products to your local retailer when purchasing an equivalent new product. By properly disposing of this product, you will help conserve valuable natural resources and help prevent potential negative environmental and human health impacts that could result from improper waste disposal. For further details, contact your local authority or nearest collection point. This symbol is valid in the European Union. If you wish to dispose of this product, please obtain the necessary information on the correct disposal method from your local authorities or your dealer.

CERTIFICATION



The product is approved for use in European Union countries and bears the CE mark. The declaration of conformity and all documentation is available from the importer.

IMPORTER

SANTADRIA s.r.o. , Mojmírovo nám. 75/2 Brno, 612 00
IČO: 29195454 | DIC: CZ29195454

If you have any questions about the product, please do not hesitate to contact us:
Email: mh@santadria.cz Tel: 737 196 160

Printing errors, features and changes in technical parameters reserved.

ES DECLARATION OF CONFORMITY

The undersigned, representing the company

SANTADRIA s.r.o.

Mojmírovo nám. 75/2, Brno 612 00, ID: 29195454

hereby declares, as an importer into the EU, that the following product with the trade name

SHIATSUSOFT massage machine

marketed under the model name: **BR-2400**

EAN: **8594166119239**

belonging to the **BeautyRelax** brand group of products

is in compliance with the following European standards:

According to EMC Directive 2014/30/EU

EN 55014-1:2017

EN 55014-2:2015

Based on compliance with these directives, the product has been CE marked.

The product also complies with RoHS Directive 2002/95/EC (RoHS). All materials are health-safe. They do not contain any heavy metals such as cadmium and other toxic substances.

For safe use, it is essential to strictly follow the instructions in the instruction manual included in the package.




SANTADRIA s.r.o.
Mojmírovo nám. 75/2, 612 00 Brno
DIČ: CZ29195454 Tel.: 737 196 160

Martin Hrbek, Managing Director

Place of publication: Brno

Date: 9.6.2022