BRAUN

Skin i·expert

2131137/V-25

www.braun.com

PL7

Braun Infolines

English	4	UK	0800 783 7010
Polski	12	IE	1 800 509 448
Český	22	PL	801 127 286
Slovenský	31	CZ	800 113 322
Magyar	40	SK	800 333 233
Hrvatski	49	HU	(06-1) 451-1256
Slovenščina	58	HR	01/ 6690 330
Türkçe	67	SI	802 822
Română	76	TR	0850 220 0911
Ελληνικά	85	RO	0800 890 194
Български	95	GR	800 801 3458
Українська	105	BG	008001184973
		UA	0 800 505 000


Type 6036

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CyDen Ltd., Block A,
Bay Studios Business Park,
Fabian Way, Swansea SA18QB, UK

Skin tones and body hair types

✓ = suitable / odpowiedni / vhodné / vhodné / megfelelő / prikladno /
 = Kullanim için uygun / primerno / adecvat / κατάλληλο / подходящий /
 підходить


✗ = not suitable / nieodpowiedni / použití není vhodné / nevhodné / nem megfelelő /
 = nije prikladno / kullanim için uygun değil / ni primerno / neadecvat / ακατάλληλο / не e
 подходящий / не підходить

English

Our products are designed to meet the highest standards of quality, functionality and design. Before using the device, please read the use instructions carefully, including all contraindications, warnings and safety information, and keep them for future reference.

Using professional light technology, the Braun Skin i-expert Intense Pulsed Light (IPL) device is intended for the permanent reduction in unwanted hair.

Ideal body areas include:

Women	Men
<p>Legs, arms, underarms, intimate area*, face (below cheek bones), fingers, toes</p> <p>CAUTION DO NOT USE... around the eyes, along the forehead, on scalp, nipples</p>	<p>Shoulders and below (chest, back, underarms, arms, stomach, intimate area*, legs, fingers, toes)</p> <p>CAUTION DO NOT USE... on face, neck, scalp, nipples</p>
<p>*Intimate area: VIO</p>  <p>Bikini area between legs buttocks and skin around anus (peri-anal)</p> <p>DO NOT USE... in mucous membranes such as labia minora, vagina or on the anus</p>	<p>*Intimate area:</p>  <p>Area above penis buttocks and skin around anus (peri-anal)</p> <p>DO NOT USE... on penile shaft, scrotum or anus</p> 

IPL is not suitable for everyone. It is most effective on light to medium skin tones and dark hair colors. Please, refer to the skin tone and hair color chart on page 3 to see if the device is suitable for you.


Description

- 1 Treatment window
- 2 SensoAdapt™ skin tone sensors
- 3 LED Power bars
- 4 Mode button with indicator lights
- 5 Treatment button
- 6 Handle
- 7 Base station
- 8 Mains cord
- 9 Glass filter (integrated UV filter)
- 10 Treatment head
 - a Smart Flex head*
 - b Standard head *
 - c Wide head*
 - d Intimate head*
 - e Mini precision head*

* not with all models / color can differ

How does it work?

Braun Skin i-expert works beneath the skin surface by targeting the melanin in the hair follicle, helping to break the cycle of hair re-growth. With continued, periodic use, the light energy helps prevent unwanted hair from reaching the skin surface.

The intelligent SensoAdapt™ skin sensors (2) continuously read your skin tone before every flash, automatically adapting the light intensity for best efficacy and safety. This allows treatment of appropriate skin tones, but prevents treatment if your skin tone is too dark.

Contraindications

Do not use the device ...

- if your skin is darker than skin tone 5 shown on the skin tone chart on page 3. Since dark skin absorbs more light energy, treating very dark skin may cause discomfort/pain and adverse effects (e.g., burns, blisters, discoloration or scarring) and injure your skin.
- if you are pregnant, lactating, or under the age of 18 as this device has not been tested with these individuals.
- if you have a history of skin cancer or precancerous lesions (e.g., nevi or a large number of moles).
- If your skin has been artificially or naturally tanned recently. Your skin may be sensitive following sun exposure and particularly susceptible to side effects of IPL treatment (e.g. burn, blister, discoloration or scar to your skin). Also, avoid unprotected exposure to direct sunlight after your IPL treatment.

Certain conditions can also prohibit the use of this device.

Do not use the device if any of the following applies:

- You have a chronic medical condition.
- You have a chronic skin disease (e.g. psoriasis or vitiligo) or damaged skin (e.g. sunburn, cuts, open wounds or active infections) in the areas you wish to treat.
- You have varicose veins in the area you wish to treat.
- You have a known sensitivity to sunlight (photosensitivity) or are taking medication that may make the skin more sensitive (e.g. Retin A, Accutane and/or other topical retinoids).
- You have had a professional skin peel treatment on the area you wish to treat.
- You have any of the above, this device may injure your skin or make existing conditions worse. You may experience side effects such as burns, blisters and skin color changes (hyper- or hypo-pigmentation) or scarring.

If you are unsure whether this device is safe for you to use, please consult your physician or dermatologist. If you take any kind of medication regularly and over a long period of time, please consult your physician for any potential impact on skin sensitivity to light.

Warnings

Do not use on ...

- scalp, eyebrows or anywhere near or around the eyes. Use in these areas may permanently damage your eyes.
- a tattoo or permanent makeup in the treatment area.
- nipples, genitals or around the anus. These areas may have a darker skin color and/or greater hair density and using the device in these areas may cause discomfort/pain or injure your skin.
- dark brown or black spots such as freckles, birthmarks, moles or warts in the treatment area.
- beard / facial hair for men. This may lead to permanent or uneven results which may produce unwanted changes in the facial characteristics.

It is not recommended to use the device on areas of skin where dermal fillers are present.

Do not use more than once a week on the same area. Damage of skin may occur after prolonged or repeated treatment on one site.

Eye safety

Do not look directly into the glass filter (9) or attempt to activate the device toward the eyes. To protect the safety of your eyes and those of individuals nearby, the device will only activate when both skin tone sensors are placed against the skin. If at any time the device activates when both sensors are NOT in contact with the skin, STOP use immediately and contact your retailer.

Possible eye injury (potentially leading to loss of vision) or skin injury if instructions are not followed. Protect the eye from exposure. Do not use the appliance over the eye lids or close to the eye. We recommend you look away from the device when pushing the treatment button.

Precautions

- The device is intended for a single user.
- It is not intended for use by people under 18 years old.
- People with reduced physical, sensory or mental capabilities or lack of experience and knowledge, may use the device if they have been given supervision or instruction concerning the safe use of the device and understand the hazards involved.
- Children should be supervised to ensure that they do not play with the device. Always keep the device out of children's reach.
- There are a number of potential risks - exposure to the light output causing burns or eye damage, electrocution, strangulation from the cord.

Electrical and fire safety

Regularly check the device and cords for visible signs of damage. In case of damage or cracks, stop using the device. Electric shock may result from using a damaged device. In case of defects, have the device checked and repaired by an authorized service center. Unqualified repair work can lead to extreme hazards for the user.

CAUTION

Always unplug the appliance after use and in case of a power failure to avoid risk or damage to the equipment.

Do not subject the device to heavy shocks to avoid risk of damage to the lamp.

Do not position the device so that it is difficult to disconnect from the mains supply. Do not leave the device unattended when it is switched on to avoid the risk of fire or burns.

WARNING

If the wall socket used to power the device has poor connections, the mains plug becomes hot. Make sure you plug the device into a properly installed wall socket to avoid the risk of fire and burns.

No modification of the equipment is allowed. Do not use any other component nor manipulate the device. This may damage the unit and cause it stop working.

DANGER

Never attempt to repair and open the device, as this may expose you to dangerous electrical components and to pulsed light energy, either of which may cause serious bodily damage and/or permanent eye injury.

WARNING: Water and electricity are a dangerous combination!

To avoid risk of electrocution, do not use this device in wet surroundings (e.g. in a bathroom, near a shower or swimming pool); do not let water run into the device. To avoid electric shock, never use a wet device.

Do not use if the device is damaged (e.g. cracked housing, cord damage etc). Internal parts contain electric charges that are hazardous.

Regularly inspect the treatment window (1) for damage:

Do not use if the glass filter (9) is cracked or missing. Do not use if the built-in reflector is missing or not in correct place.

Do not touch the glass filter during or directly after use since it will get hot.

Do not use petroleum-based or flammable cleaning agents (alcohol, acetones) on the device as these pose a fire risk.

Do not use abrasive cleaning agents. These may damage the device or scratch the glass filter.

Do not attempt to use this device other than for removing unwanted hair.

Do not use the device on a surface other than skin. Do not press the treatment button when the treatment window (1) is not placed in full contact with the skin.

Do not block the cooling vents. The device may overheat and pose an electric shock risk or fire hazard.

Under continued use the temperature of the device may reach a maximum of 60°C. After treatment, remove the device from the skin.

How to use

Preparing the treatment area

Remove any visible hair on the skin before using the device. You may use your preferred hair removal method. To keep the treatment window clean and intact, ensure there is no hair remaining on the skin surface. Clean the treatment area and pat it dry.

To power on the device:

Use the mains cord (8) to connect the base station to an electrical outlet. Press any button (4 or 5) to activate the device. The LED power bars (3) will light up once and a fan starts running. This means the device is ready for use. During operation, the device will heat up slightly.

Automatic light adaptation

When the skin tone sensors (2) have skin contact, the device reads your skin tone and selects the appropriate light intensity which is indicated by white LEDs on the Power bars (3). These confirm that the device is ready to flash.

	light skin tone → high light intensity
	dark skin tone → low light intensity with fast flash rate
 no LEDs	no skin contact → device will not flash
 red LED	invalid skin tone → device will not flash

Energy settings

Push the mode button (4) to reduce the light intensity. You can toggle through the following 3 settings:

Power mode (3 lights)

Balanced mode (2 lights)

Sensitive mode (1 light: lowest light intensity with fast flash rate)

Balanced and sensitive modes are recommended for sensitive areas and first-time users. With reduced light intensity, more treatments may be required to see the desired effects.

The power bars flash to confirm a change in the light intensity setting.

When holding the mode button (4) longer than 3 seconds, the device is deactivated.

After 10 minutes without use, the device is automatically deactivated.


Press any button to re-activate it.

Skin Test

Before your first treatment, we recommend you test your skin for reaction on a small area of skin where you plan to treat. Perform two flashes next to each other as described below in chapter «Treatment methods», section «STAMPING». Wait 15 minutes to examine the test area. If the skin appears normal, begin your first treatment. If a skin reaction develops, switch to sensitive mode and repeat the test. Wait for another 15 minutes. If the skin reaction persists, wait for another 24 hours. If skin reaction is present after this period, discontinue using the device and consult your physician.

Treatment methods

The area to be treated must be hair-free, clean and dry. Before pushing the treatment button, look away from the device to avoid the bright light. In rare cases, a user may find the flash too bright. However, there is no safety risk associated with this. For increased comfort, we recommend treating in a well-lit room to help reduce the perceived brightness.

<p>GLIDING for fast and efficient treatment of large areas (e.g. legs).</p> 	<p>STAMPING Spot-by-spot treatment for targeted flash release (small, hard-to-reach areas)</p> 
<p>1. Place the treatment window firmly on your skin, ensuring the skin sensors have skin contact. The power bars (3) show white LEDs</p>	
<p>2. Press and hold the treatment button (5). The device will begin to flash.</p>	<p>2. Press and instantly release the treatment button (5). The device will emit a flash.</p>

3. Move the device <u>slowly and continuously</u> from the bottom to the top of the treatment area. Make sure to always keep full skin contact so that the device flashes continuously.	3. Lift off the device and place it onto the next adjacent spot to be treated. Wait for the power bars to show white LEDs. Then press and release the treatment button.
Once the gliding stroke is completed, lift the device off the treatment area. Repeat steps 1-3 to treat the next adjacent area in a parallel stroke.	Continue treatment spot by spot until you have completed the entire area you wish to treat.

Make sure you avoid missed areas by placing the treatment window directly next to the previously treated spot without leaving any space between treated areas. Avoid overlapping treatment or treating the same area twice as this may increase the chance of producing adverse reactions.

If you feel intense pain, stop treatment immediately and consult your physician before using the device again.

Where to use the treatment heads

- a SmartFlex: all (including curved) areas
- b Standard: all areas
- c Wide: large body areas
- d Intimate: small, curvy, hard-to-reach areas
- e Mini precision: smallest, curvy, hard-to-reach areas

How to treat sensitive areas

Before treating, remember to perform a skin test. Remove any visible hair and cosmetics, lotions or creams. You might want to confirm full hair removal using a mirror.

Face (females only)

When treating the upper lip area, avoid contact with the lips and the nostrils. The central area below the nostrils might be sensitive.

Do not use on permanent / semi-permanent lip liner on the lip margin.

TIP: You can press your lips together to ensure that the light does not activate on the upper lip itself.

Hair growth under the chin can be denser and coarser. This can sometimes lead to slight skin reddening following a treatment. This redness is temporary and should completely disappear after 24 hours.

Intimate area (refer to page 4 for suitable treatment areas)

Stand, sit or crouch down to reach all areas to be treated.

For your initial treatment, we recommend using the intimate head with the stamping method and the sensitive mode as described above. Please select your preferred energy setting that will provide the best comfort for continuous treatment.

NOTE:

- If redness or a burning sensation occurs and lasts for several minutes, apply a cool cloth/pad to soothe the area.
- If any pigment changes in the skin, e.g., lightening or darkening, are observed, stop using the device. Seek medical advice to determine the cause which may be due to an underlying medical condition.
- Since the intimate area tends to have darker skin color, the device may not flash due to safety reasons even if the device has good contact with the skin.

Treatment plan

Start-up phase (weekly or bi-weekly)	1 treatment every 1-2 weeks (4- 12 weeks)
Maintenance phase (monthly or bi-monthly)	continue treating to your needs, e.g. every 1-2 months.

The typical hair growth cycle varies between individuals and could take up to 16 months depending on body area. Since only hairs in their growing phase are susceptible to treatment with light, it is important to do multiple treatments to achieve continuously smooth skin. For this reason, we recommend a start-up phase with one treatment every 1-2 weeks over a period of 4 - 12 weeks.

For the facial area or those with darker skin tones, you may need to extend the start-up phase to achieve your best results.

In case you no longer see hair growing, you can switch to the maintenance treatments.

After treatment

Once you have completed your treatment, unplug the mains cord from the electrical outlet.

After treatment, we recommend using sunscreen SPF > 15 whenever treated areas may be exposed to the sun. We recommend avoiding exposure to artificial light sources, such as tanning booths or solariums.

Avoid any activity or skin products that may irritate your skin for 24 hours after treatment. This includes using hot tubs or saunas, wearing tight-fitting clothing, using bleaching creams, perfumed products or exfoliation products. The above may cause irritation, such as itching and redness, in the treated areas.

Cleaning and storage

After each treatment, inspect the treatment window (1), the skin tone sensors (2), the glass filter (9) and the built-in reflectors for damage. Do not touch or clean the glass filter immediately after use to allow it to cool down. Wipe the treatment window with a dry, lint-free cloth if any marks or black spots are visible. Do not rinse or immerse the device in water or any other liquid.

To ensure maximum performance, keep the heads (10) with their built-in reflectors free from dirt, hair or other debris. Do not clean in water, but use a dry cloth for cleaning them carefully.

Store the device in a cool, dry place. Make sure the treatment window and skin sensors are protected from damage.

Possible side effects

Some people will experience slight discomfort (e.g. heat or redness of the skin) when using the device – this is normal. You may feel heat and tingling when activating the flash, followed by redness in the skin. This type of reaction usually disappears within 24 hours. Adverse reactions are usually immediate or appear within 24 hours. In rare cases, they may take up to 72 hours to appear. If you experience abnormal skin reactions please contact your physician.

Side effect	How to assess and react
Warm feeling or tingling sensations during treatment which typically disappear after a few seconds to a minute and decrease with continued use.	This is expected and normal for all IPL treatments. You can keep on using the device as instructed.
Discomfort or pain is intense during treatment or persists after a treatment.	Stop using the device and consult your physician before using it again.
Skin redness during or after treatment which disappears within minutes to several hours.	This is expected and normal for all IPL treatments. You can keep on using the device as instructed once skin redness has disappeared.
Skin redness which does not go away within 24 - 48 hours after treatment.	Stop using the device and consult your physician before using it again.
In very rare cases:	
Swelling and redness around the hair follicle that disappears within two or three days.	You can keep on using the device as instructed once swelling and redness around the hair follicle has disappeared.
Swelling and redness around the hair follicle that does not disappear within two to three days.	Stop using the device and consult your physician before using it again.
Temporary changes in skin color (lightening or darkening).	If your skin color changes, stop using the device and consult your physician.
Increase in hair growth in areas of treatment (more likely to be seen with subjects of Mediterranean descent).	If this increase persists over several treatments, stop using the device and consult your physician.

In case of adverse events please contact your national competent authority and the device manufacturer.

Trouble shooting

The following troubleshooting guide may help you solving initial problems you may experience with this device.

Problem	Solution
Power bars do not light up when treatment window is placed on the skin.	Make sure the device is connected to an electrical outlet. Press any button (4 or 5) to activate it. Ensure the skin tone sensors are in full contact with the skin. If needed, move the device in a way that the skin tone sensors get skin contact.
The device doesn't work on small areas.	Use a smaller head if available with your device. Make sure both skin tone sensors are in contact with the skin.
Flashing sequence slows down.	Device is in heat-protection mode. To get back to normal speed, allow device to cool down.
On the power bar, only the first LED illuminates red.	Invalid skin tone reading. Your skin is either too light or too dark. The device will not flash.
Centre mode button light is flashing.	Device is deactivated. Press any button to re-activate it.
First LED on power bars flashes red.	Head (10) has been taken off. Re-attach it.
The power bars show an unfamiliar pattern of one red LED shining along with individual white LEDs.	Device error. Device will not operate. Unplug and reconnect. If error persists, contact a Braun Service Centre.
All white LEDs flash.	End of product lifetime. Maximum number of flashes (400,000) has been reached.
The device is broken, cracked or appears defective.	Do not use if the handle, the base station or the mains cord is damaged. If you are in any doubt about the safety or suspect it is damaged in any way, you must not use it.

Disposal

The product contains recyclable electronic waste. To protect the environment, do not dispose of it in the household waste, but take it to appropriate local collection points.

Labels and symbols

The base station is marked with the following symbols:

	Type BF Applied Part		Serial number
	Class II equipment		WEEE waste disposal mark of the European Union
	Optical radiation		Conformity mark of the European Union
	Follow instructions for use		Manufacturer
	Not suitable for use in a bathtub or shower		Keep Dry

Technical specifications

Repetition rate:	every 0.5 – 0.9 seconds (depending on setting), continuous operation
Max optical output:	6 J/cm ²
Pulse length:	2.5 +/- 0.5 msec (full discharge FWHM)
Pulse profile:	partial free discharge with total pulse duration = 2.0 ms to 8.0 ms dependant on skin tone
Wave length:	510 - 2000 nm
Power line input:	100 - 240V~, 50-60Hz, 2.0 - 1.0 A
Operating temperature:	5°C - 40°C
Optimum operating temperature:	15°C - 30°C
Operating humidity:	up to 93% R.H. non-condensing
Operating pressures:	700 hPa to 1060 hPa
Treatment area (spot size):	
SmartFlex/Standard head:	3 cm ²
Wide head:	4 cm ²
Intimate head:	1.5 cm ²
Mini precision:	1 cm ²

Transport Conditions:

Temperature range:	- 25°C to 70°C
Humidity:	up to 93% R.H. non-condensing
Pressures:	500 hPa to 1060 hPa

The device includes 400,000 flashes and is intended for use up to 6 years by a single user.

This device complies with the EU Radio Equipment Directive 2014/53/EU and the UK Radio Equipment Regulations 2017.

The device contains a radio module which operates in 2.401 to 2.473 GHz band with a maximum power of 20mW (13 dBm). Hereby, Braun GmbH declares that the radio equipment for type 6036 is in compliance with Directive 2014/53/EU and UK Statutory Instruments.

The device is intended to connect to a network with the following characteristics.

- Open or Private network with internet access
- WLAN 802.11 b/g/n
- Frequency band: 2.4 GHz
- Channel 1-11
- WPA/WPA2/WPA3

The correct operation of the device does not rely on the connection to the network.

Further information can be found in the terms and conditions of the mobile app.

The full text of the EU and UKCA Declarations of Conformity can be found in the following URL:

<https://uk.braun.com/en-gb/ce>

This device complies with the following EMC and Radio standards and latest amendments applicable to Household and Medical Devices.

EN 55014-1:2021, EN 55014-2:2021, EN 61000-3-2:2014, EN 61000-3-3:2013, IEC 60601-1-2:2014 (Ed 4.0), EN 62311:2008, EN 62233:2008, ETSI EN 301 489-17 (2020-09) V3.2.4, ETSI EN 300 328 V 2.2.2 (2019-07), FCC CRF 47 Part 15B, 15.107, 15.109 and Part 15C, ICES-003 Issue 7.0.

This device also complies with the following Safety Standards and latest Amendments applicable to Household and Medical Devices.

IEC 60335-1:2010/AMD2:2016, IEC 60335-2-113: 2016, IEC 60601-1:2005/AMD1:2012, IEC 62304:2006/A1:2015, IEC 60601-1-11:2015, IEC 60601-1-6:2010/AMD1:2013, IEC-62366-1:2007/AMD1:2014, IEC 60601-2-57:2011, IEC 60601-2-83:2019, IEC 62471:2006.

This device is classified as Risk Group Exempt for all hazards

For Spectral Irradiance, please refer to Annex 1 on page 117.

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