

Electric toothbrush for children

Instructions for use

#### IMPORTANT SAFETY PRECAUTIONS

Please read the safety instructions before using this product. These safety instructions will help you to use this product use this product safely and correctly, but also prevent you from injuring yourself or someone else.

When using electrical equipment, especially in the presence of children, the basic safety precautions should always be observed, including the following:

#### READ ALL INSTRUCTIONS BEFORE USING THE WARNING:

To reduce the risk of burns, fire or personal injury:

This product is not intended for people with physical, sensory or lack of experience and knowledge, unless unless they have been given adequate supervision or instruction in the use of the product, by an adult who is responsible for their safety.

Children should use under supervision so that they do not play with the product it is not a toy.

Use this product only for its intended use and as directed.

written in this manual. This product is designed for cleaning only teeth and gums. Do not use any accessories or the brush head for anything other than that recommended by the manufacturer.

This toothbrush is a personal care device and is not intended to use by multiple patients in a dental office

If your toothpaste contains peroxide, baking soda or other bicarbonate (common in whitening toothpastes), then thoroughly clean the brush head and handle to avoid cracking the plastic parts of the product.

Keep the product away from hot surfaces.

If discomfort or pain in the mouth occurs, discontinue use the toothbrush and contact your dentist.

Do not clean this product in a dishwasher.

Do not use the toothbrush outdoors or in an aerosol environment.

(spray) or where oxygen is used

Never drop or insert any objects into any opening.

For indoor home use only.

If you do not use the toothbrush for a long time, please change the batteries, so as not to damage it in case of "leaking" batteries.

#### HEALTH WARNING

If you have had surgery in the mouth or gum area in the last two months, please consult your dentist before using the toothbrush.

Consult your dentist if you experience excessive bleeding during use brush or if slight bleeding persists after one week of use.

Consult your physician before using this toothbrush or implanted device manufacturer if you have questions about joint use of this product and a pacemaker or other implanted medical device.

Consult your physician if you have any medical problems.

#### KEEP THESE INSTRUCTIONS

The children's electric toothbrush is designed for children up to 5 years of age. It is not intended for adults.

For children under 8 years of age, the dentist recommends that brushing be under supervision by a parent or adult caregiver.

#### HOW TO USE

1. Put the toothpaste on the head of the toothbrush.
2. Insert the toothbrush head into the mouth.
4. Press the ON button on the toothbrush lightly.  
and start brushing your teeth.

Recommendation:

Brush your teeth evenly

- . Top on the outside and Bottom on the outside
- . Top on the inside and Bottom on the inside
- . Guide the brush head slowly from one tooth to the other.

Brush one spot for a few seconds before moving to the next tooth.

Don't push the brush too hard or scrub too hard, simply allow the brush to work itself.

While brushing, lightly press the OFF button to stop brushing.

#### NOTES:

2 minutes is long enough,

for your teeth to be cleaned properly. If you want a more thorough cleaning, you can brush your teeth for longer than 2 minutes.

If you use an electric toothbrush for the first time, it may bleed slightly.

The bleeding should stop after a few days, before you get your gums on the electric to get used to the electric toothbrush. If the bleeding persists, consult your dentist.

#### MAINTENANCE

Switch on the toothbrush and rinse the brush head under running water.

Clean the brush with a slightly dampened cloth or paper towel.