BEAUTIFLY

TouchSense

Body

series


8-IN-1 ANTI-CELLULITE MASSAGER

USER MANUAL


 /beautifuly.global  /beautifuly_global

 /@beautifuly.global  /@beautifuly.global


f1

f2

f3

f4

[youtube.com/@beautifly.global](https://www.youtube.com/@beautifly.global)

[facebook.com/beautifly.global](https://www.facebook.com/beautifly.global)

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Rated voltage: DC 5V
Rated current: 1A
Power: 16W
Timer: 15 minutes
Wavelength range: 620 - 640nm

Safety:

Read this user manual before using this product and follow the instructions contained therein when using the device.

Device description [page 2 - f1]:

1. Power/battery level indicator.
2. Power switch/power control.
3. Charging connector.
4. Device body.
5. Massage head.

Spare parts description [page 2 - f2]:

1. Flower head.
2. Wave head.
3. Roller head.
4. Flat head.
5. Ring head.
6. Brush head.
7. Linear head.
8. Abrasive head.

Operation [page 2 - f3]:

1. Hold the device in your hand. Grasp the cover and turn it counterclockwise to remove it.
2. Locate the two small gaps between the head and the body of the device. Grasp these areas and pull to remove the device.
3. Replace the massage head with the desired one, align it with the desired location, and then press down.
4. Replace the mesh cover and turn it clockwise, aligning it with the threads, until it is securely tightened.

Maintenance and Storage

1. Use a soft cloth or mild detergent to clean the device. The massage heads can be washed directly under running water.
2. Keep water, liquids, and moisture away from the device, as they may damage the internal electronic components.

Warnings

1. Do not use the device until one hour after a meal.
2. Keep the device out of reach of children.
3. Do not use the device in the bathroom or in areas with high humidity. Do not immerse the device in water. Keep away from liquids and heat sources. Keep the ventilation openings free from dust and debris.
4. Store in a cool and dry place. Avoid direct sunlight.
5. Ensure the device is turned off before replacing the massage heads.
6. Do not loosen the screws or attempt to disassemble the device.
7. Do not use any massage head without the mesh cover near hair.
8. The device should not be used by: people with skin conditions, such as active inflammation (e.g., acne pustulosa, eczema), bacterial, viral, or fungal infections (e.g., herpes), open wounds, burns; people with autoimmune diseases; People with cancer, especially in the area of the lesion or undergoing treatment (e.g., skin cancer); pregnant women; people with circulatory system diseases, such as hypertension, varicose veins, or thrombosis; people with sensory disorders, such as diabetic neuropathy; people with medical implants, such as a pacemaker; people taking photosensitizing medications, such as certain antibiotics or retinoids; and children. Particular caution should be exercised by the elderly, those who have undergone surgery, and those with chronic pain of unknown cause.
9. The vibrating massager with a red light function is designed to improve the user's daily comfort by combining gentle vibrations and warm light. Regular use of the device can contribute to: a feeling of relaxation in tense areas

of the body, improved well-being and relaxation after physical exertion or a day at work, promoting muscle relaxation, improving the perceived elasticity of the skin, and giving the skin a more rested and radiant appearance. The emitted red light creates a pleasant sensation of warmth, which can promote relaxation and comfort during use. The device's vibrations provide a gentle surface massage that can be perceived as soothing and relaxing.

10. Although the device has been designed with safety in mind and is equipped with mechanisms to limit potential hazards, its use may involve some residual risk. The following may occur when using the device: temporary redness or irritation of the skin, discomfort with prolonged or excessive use, excessive skin heating at the application site, aggravation of existing conditions (e.g., in the event of improper use), and photosensitivity reactions in susceptible individuals.

11. To ensure safe use of the device: use the device according to the instructions (duration, frequency, body areas), do not use on damaged or diseased skin, discontinue use if pain, burning, or discomfort occurs, consult a specialist if you have a chronic condition or have any doubts, and do not use the device in a manner other than recommended by the manufacturer.

12. If you have any doubts, we recommend consulting a doctor or other qualified professional.

13. The red light vibration massager is a reusable product, intended for long-term use. The expected safe use period is 2 to 5 years, depending on the frequency of use and proper maintenance. After this time, the effectiveness and safety of the device may gradually decrease.

14. To maintain the safety and proper operation of the device, the user should: regularly clean the device according to the instructions, inspect its technical condition (e.g., casing, cables, light-emitting components), do not use the device if it is visibly damaged, store the product in a dry and clean place, and use only original or recommended accessories.

15. Recommended frequency: once daily or 3–5 times a week, maximum frequency: twice daily (with an interval of at least 6–8 hours).

16. If the recommended parameters (time, frequency, or intensity) are exceeded, the following may occur: skin irritation (redness, burning, dryness), overheating of superficial tissues, increased skin sensitivity to touch or light, discomfort or pain at the treatment site, worsening of existing skin or muscle complaints, and temporary deterioration of well-being due to overstimulation.

17. The skin surface to which the treatment is applied should be relatively smooth and free of major irregularities (e.g., scars, nodules, skin lesions).

18. Each treatment session should cover a separate area of skin, the size of which does not exceed the maximum recommended treatment area. Overlap between treatment areas should be minimal and limited to small skin sections to avoid excessive accumulation.

and energy. Treatment areas should be spaced evenly between uses.

19. Regularly inspect the device's condition (housing, cables, light sources). Discontinue use if there is visible damage or malfunction. Store in a dry and safe place, protecting it from moisture and mechanical damage.

20. Cleanse the skin before treatment and leave it dry; use gentle, neutral cosmetics after treatment. Some products may affect the effectiveness of the treatment or increase the risk of adverse reactions. Products containing retinoids, AHABHA acids, or high concentrations of vitamin C may increase the skin's sensitivity to red radiation and vibrations, which may lead to irritation or redness. Products with a high alcohol content may dry out the skin and increase the risk of discomfort during the treatment. Warming or cooling products may intensify the sensation of heat or alter the sensation of vibrations.

21. Before using the device, consult your doctor if you are taking medications that may affect your skin's response to light. Photosensitizing medications (e.g., some antibiotics, antihistamines, and cancer medications) may increase your skin's sensitivity to light, leading to a risk of burns or irritation. Anti-inflammatory medications and topical steroids may reduce your skin's resistance to irritation and affect your skin's response to red light. Anticoagulants (oral or topical) may increase the risk of bruising or minor bleeding during intense vibration massage.

22. When using the device, consider the effects of other radiation sources on the skin to avoid overexposure and an increased risk of irritation or discomfort. Tanning beds and UV lamps. Using UV radiation concurrently with treatment may lead to a cumulative light effect and an increased risk of skin damage. Other devices emitting red or infrared light. Using them simultaneously or in close succession may increase the risk of overheating or discomfort.

23. This device is intended for home use in a dry environment with an air temperature between 10 and 30°C. Do not use the device in places where direct sunlight is directed at the skin.

24. Wear appropriate eye protection if you are exposed to radiation due to reflection, misuse, or improper handling of the product.

25. Do not use the brush head and the matte head with the mesh cover simultaneously.

Use and Operation:

1. Fully charge the device when the battery is low.

2. Press the On/Off/Power Adjustment button to turn the device on or off.
3. Briefly press the On/Off/Power Adjustment button to change the speed as needed. There are 9 levels to choose from.
4. For first-time use, it is recommended to start with a lower massage speed and gradually increase it as you become accustomed.
5. To keep the device in optimal condition, it is recommended to use it for no longer than 15 minutes at a time; the device will automatically turn off after 15 minutes of continuous use.
6. Turn off the device after use.

Charging Instructions:

1. Before first use or after a long period of inactivity, fully charge the device.
2. Before charging, turn off the device and use the included USB cable.
3. The battery indicator flashes during charging and remains solid when fully charged.
4. Once fully charged, disconnect the USB cable from the charging port.
5. If the device is not used for an extended period, remember to fully charge it at least once every 3 months.

More information about the company and products [page 2 - f4]:

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