medivon®

User manual

Gun massager

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1. Device operation

The device is very easy to use - a simple push of a button is all it takes to start a relaxing full-body massage, allowing you to unwind after a hard day at work or effectively massage tense muscles after a workout.

Once activated using the power button on the back of the unit, you can cycle through the massage modes, gradually increasing the gears.

To switch the device off, press the POWER ON / OFF button again.

Before use, it is recommended to fully charge the device battery using the charger. Only use the charger supplied for this purpose. A fully charged battery will be indicated by a backlight. It is recommended not to fully discharge the device to extend the battery life. It is prohibited to use the massager while the battery is charging. There is a risk of damage or fire.

Note: only replace the massage heads when the device is switched off. Simply pull them off the device by pulling (they do not have a screw thread).

2. Massage and settings

Start by choosing the right massage head for you and your needs. Make sure it is firmly attached and not damaged. Then start the device using the power ON/OFF button.

The massage force is adjustable by the user. Depending on the force with which you press the working massager against your body, you will get different final results. It is recommended to use the massager gently during the first massage.

To change the gear of the device, touch the switch and select accordingly: mode 1 for warming up the muscles, mode 2 for relaxing the fascia, mode 3 for removing lactic acid, mode 4 for a deep massage, mode 5 for a professional massage. To switch off the device and end the massage, press the Opower OFF button.

3. Dangers

1. Do not use spare parts that are not approved by the manufacturer.
2. Non-professional service personnel should not dismantle or replace parts of the appliance without authorisation.
3. If the power cord is damaged, it must be replaced by appropriate professional personnel to avoid danger.
4. Do not use or store the product in a place with high humidity or high temperature.
5. Unplug the power cord when you have finished massaging.
6. Do not use the product for longer than 15 minutes at a time.
7. If any side effects occur, stop using the massager immediately.
8. If the temperature is too high or the product starts to smell bad or give off a burning odour, stop using it immediately and disconnect the power supply.
9. Do not carry the product by pulling on the power cord to prevent the power cord from coming loose.
10. This product is not intended for use by persons who lack physical, sensory or mental ability or experience and knowledge (including children).
11. People who are not sensitive to heat should use the heating function with caution - preferably under the supervision of a carer.

4. Warnings

If you have any of the following conditions, consult a doctor before using:

- pregnant women, children, or the elderly,
- individuals suffering from skin infections or having skin redness and burns,
- patients using a pacemaker or similar implanted medical devices,
- patients with heart diseases or malignant tumors, etc.
- people who are under treatment or experiencing discomfort,
- individuals suffering from osteoporosis,
- people with high or irregular blood pressure.

5. Cleaning and maintenance

Turn off the product and unplug the power cord before cleaning. Wipe the surface of the product with a dry, soft cloth. Do not clean the product and accessories with harsh detergents such as benzene, thinners, and store in a cool, dry place when not in use. Also avoid exposure to direct sunlight or humid environment, otherwise it may cause color fading or product damage.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after 13.08.2005.

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The manufacturer declares under his sole responsibility that the Medivon Relix device complies with Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, as well as Directive 2011 of the European Parliament and of the Council /65/EU of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to the Directive of the European Parliament and 2011/65/EU as regards the list of restricted substances. The full text of the declaration of conformity is available at the internet address www.medivon.pl/deklaracje

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska
Model: Medivon Relix