
medivon.

Medivon Body Pulse Go

Model: AST-303J

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
User manual EN

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B

A


C


D


E

F

A

H

G

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This manual applies to the Medivon Body Pulse Go product. Before using the device, read this manual and keep it in a safe place. Technical data: Power supply: 5 V / 1 A | Rated power: 0.5 W Battery capacity: 90 mAh.

1. Intended Use of the Product

The AST-303J device is a medical neuromuscular stimulator (EMS/TENS). It is intended for home and clinical use for:

- Relieving muscle pain and neuralgia.
- Stimulating weakened muscles in the rehabilitation process.
- Improving local blood circulation.

The device is intended for adults. It does not require specialized medical training for home use.

2. Safety Precautions

To ensure proper and safe use of the product, read the user manual carefully before use. The safety precautions described in this manual are intended to ensure safety, proper use of the product, and to prevent injury to the user and others, as well as damage to the property.

Safety precautions regarding use: If you feel anything abnormal on your skin or in your body, stop using it immediately and consult a doctor. Failure to do so may result in an accident or deterioration of health. Using the product at a very high level or for a long time can cause muscle fatigue, tingling, or skin redness. While using the product, remove all jewelry or accessories such as necklaces, earrings, or belts. Failure to do so may result in an accident, injury, or damage. During use, keep electronic devices such as wristwatches, pedometers, mobile phones, etc., away from the body. Failure to do so may cause malfunction of the electronic device. Do not heat or cool the product before use, as it may deform it. Do not use in a way that exposes the product to sudden temperature changes, e.g., pouring hot water on it, using it in a sauna, etc. Clean the area where the product will be used. **Pay attention to the following points when handling gel sheets:** Keep away from direct sunlight, heat, and moisture. Store in a cool, dark place. Do not leave gel sheets outside and do not store them with dirty objects. If gel sheets get dirty, they may stop adhering properly, electrical stimulation may weaken, you may feel pain or suffer burns. Although gel sheets may start to peel off more easily after 25-30 uses (depending on frequency of cleaning, skin condition, or storage conditions), cleaning them helps extend their life. If a gel sheet accidentally sticks to clothing, do not pull it off by force. Moisten the fabric on the other side where the sheet is stuck, and it will come off easily. Applying gel sheets: Do not use damaged, worn, or dirty gel sheets.

This could result in burns. Important: Gel sheets are applied in the same way to both the large gel pad and the small gel pad. If the gel sheets are attached incorrectly, they may lose their adhesion to the skin, stop working normally, or wear out faster than usual. Gel sheets are consumables and should be purchased separately.

3. Danger and Warning

Never use simultaneously with the following types of medical electrical devices:

- Implanted medical devices such as pacemakers.
- Medical electrical equipment used for life support, such as a heart-lung machine.
- Wearable medical devices such as an ECG.

The following people should not use this product. Failure to follow this may result in an accident or problems with the body or skin:

People with heart disease and women who are pregnant or immediately after childbirth.

Everyone with the following conditions should consult a doctor before using this product. Failure to do so may result in an accident or problems with the body or skin:

- People prone to allergies; People with fever (38°C or higher); People with perception disorders or skin problems (infection, inflammation, or redness); People requiring acute (painful) care, such as spinal curvature, sprains, muscle strains, etc.; People with blood flow disorders; People taking medications; People with bone diseases such as osteoporosis; People experiencing any abnormal cranial nerve sensations (e.g., epilepsy, etc.); People with abnormal blood pressure; People sensing abnormalities in heart function or cranial nerves; People whose skin is irritated by gel sheets; People experiencing non-specific physical anomalies; People with perception disorders due to significant peripheral circulation disorders, e.g., in the course of diabetes; People undergoing treatment in a medical facility.

Even if none of the above conditions apply to you, consult a doctor if you have any concerns about using the product.

Please do not use this product on the following body parts; it may result in an accident or problems with the body or skin:

round the heart, such as the chest, shoulder blade, or parts of the neck upwards (head, face, etc.). On wounds or mucous membranes, such as eyes, mouth, or genitals. On body parts where depilation was performed on the same day. On body parts where plastic surgery was performed. On body parts with metal, plastic, or silicone implants. Spine. Tattooed areas. Abdomen during menstruation.

Please do not place the device in any other place not specified in this manual.

- Infants and small children who are unable to understand the instructions should not use this product. This can result in an accident or problems with the skin and body.
- Keep the button battery and battery cover in a place inaccessible to children to prevent them from being swallowed. If any of these elements are swallowed, contact a doctor immediately.
- Do not disassemble, modify, or attempt repair. This may result in an accident or injury.
- Do not expose the product to shocks, such as dropping or impact. This can cause its damage. In such a case, further use may result in an accident or injury.
- Do not touch metal electrodes. This may cause electric shock.
- Do not use in combination with any other electrical or cosmetic device. This may result in malaise or nausea.

Caution

While using the product, the following safety precautions must be followed. Failure to follow them may result in an accident, electric shock, injury, device damage, or contamination.

- Do not heat the product or gel sheet in a microwave oven. Do not use for purposes other than intended and inconsistent with the instructions in the manual. Do not use the product on small children. Do not allow children or animals to play with the product. Do not spill any liquids on the product. Do not use while exercising or driving a car. Do not use in wet places such as the bathroom.
- If pain, discomfort, roughness, skin dryness, rash, redness, itching, burning, inflammation, or another unusual situation occurs after use, stop use immediately and consult a doctor. Further use may lead to worsening symptoms.
- Please do not share the device with other people. If you notice a strange smell or sound, or if a foreign body such as water gets inside, turn off the device immediately. Failure to do so may result in a short circuit, fire, or device damage.

When handling the battery, the following instructions must be followed. Failure to follow them may result in an accident, injury, or product failure.

- In the event of leakage from the button battery, immediately stop use and contact the distributor of the purchased product. If fluid from the button battery gets into the eyes, it may cause damage or blindness. Do not rub them, rinse immediately with plenty of clean water for at least 1 minute, and then contact a doctor. If fluid from the battery gets on the skin or clothing, it may damage the skin. Immediately rinse the skin with plenty of clean water. When installing the button battery, make sure it is inserted the right way. When the button battery is exhausted, it must be removed immediately.

When handling gel sheets, the following instructions must be followed. Failure to follow them may result in burns.

Please do not use damaged, worn, and/or dirty gel sheets. Do not bend or fold the gel sheet; otherwise, it may be damaged. Stick the gel sheet to the skin; if a gap occurs at the application site, the current will not flow uniformly across the entire surface of the sheet, which may cause a feeling of a strong stimulus. Turn off the product before removing the gel sheet. When using the product, always stick the gel sheet to the electrodes. Do not touch the electrodes directly when the power is on.

4. Identification of parts

Set AST-303J includes:

- Ⓐ 1. Main unit (controller).
- Ⓑ 2. Electrode pads (hydrogel).
- Ⓒ 3. USB charging cable.
- 4. User manual

5. Charging

The device has a built-in rechargeable battery.

1. Connect the USB cable to the main unit.
2. Connect the other end of the cable to a certified medical power supply (5V).
3. The LED indicator will inform you about the charging process. After a full charge, disconnect the device.

WARNING: Do not use the device while it is connected for charging.

6. User manual and method of use

Step 1: Clean and dry the skin at the application site. Ensure that there is no sweat or lotions on the skin. Absolutely do not turn on the device if the gel does not adhere to the skin – there is a risk of burning.

- Ⓓ Step 2: Connect the main unit to the electrode pad (via snaps).
- Ⓔ Step 3: Remove the protective film from the electrodes (keep it for later).
- Ⓕ Step 4: Stick the device directly to the skin in the desired location.
- Ⓖ Step 5: Press and hold the power button to turn on the device.
- Ⓗ Step 6: Select mode and intensity.

Step 7: After the session ends (default 15 minutes), turn off the device, peel it off the skin, and apply the protective film to the electrodes.

Areas of use: Abdomen and waist: If you do not feel stimulation immediately, move the product 1 to 2 cm up or down. Turn off the power while repositioning the product.; Arms: attach the product so that it wraps around the inside of the arm, and the main unit is facing down; Legs: attach the product so that it wraps around the thigh (front/inside), and the main unit is facing down.

7. Operating modes and intensity adjustment

The device offers various stimulation modes: Mode 1-2: Surface stimulation (relaxation). Mode 3-4: Deep tissue massage (pain relief). Mode 5-6: Intensive muscle work (rehabilitation). Always start the session at the lowest intensity level, gradually increasing it until a comfortable tingling sensation is felt.

Intensity adjustment: There are a total of 25 intensity levels available (from 0 to 25). During the first use, apply a low level; once you get used to it, adjust it as desired.

MODE	TYPE	FUNCTION
M0	Automatic	The following ten modes loop. By muscle stimulation, the device helps to strengthen and firm the abdomen, arms, buttocks and thighs.
M1	Flapping	Relax muscles and relieve fatigue
M2	Acupressure	Dredge meridian, harmonizing qi and blood
M3	Shiatsu	Dredge meridian, harmonizing qi and blood, stopping pain, relieving rheumatism chills
M4	Kneading	Relax muscles and relieve fatigue
M5	Tuina	Restores the physical strength, Relieve fatigue, Improving the body's blood circulation, Strengthens the human body disease-resistant ability
M6	Tapping	Dredge meridian, harmonizing qi and blood, stopping pain, relieving rheumatism chills
M7	Massage	Promoting blood circulation, pain relief
M8	Scrapping	Dredge meridian, harmonizing qi and blood stasis
M9	Accupuncture	Dredge meridian, harmonizing qi and blood, stopping pain, relieving rheumatism chills
M10	Squeezing	Dredge meridian, harmonizing qi and blood

8. Storage and Maintenance

- **Attaching to the storage sheet:** Do not allow the gel sheets to come into contact with each other, e.g., by folding the product. Do not press the gel pads too hard onto the protective sheet, as this may damage the gel surface. Please do not bend the storage sheet; this causes deformation and makes normal use impossible.
- **Placement in the storage bag:** Keep the storage sheet and storage bag clean. When storing the product, avoid hot or humid places (such as the interior of a car) and direct sunlight. When the device is not in use, store it in a safe place to avoid tripping and injury. Keep out of reach of children and away from pets. Please do not place the device near fire. If the device is not used for a long time (a month or more), it should be charged at least once to avoid shortening the battery life.
- **Disposal method:** Dispose of the device in accordance with local government regulations. Dispose of the product if cuts or cracks appear on it after long-term use.
- **Maintenance:** Before cleaning the product, make sure the power is turned off; cleaning while the device is on (ON) may cause electrical stimulation. Do not rinse the product with water. Do not wipe the product with paint thinner, benzene, or alcohol. Make sure your hands are clean when handling the gel sheets. **How to care for this product:** Gently wipe sweat or moisture off the product with a soft cloth, avoiding scratches. If the device is dirty, moisten a cloth with lukewarm water with the addition of 3–5% detergent, wring it out well, and wipe the dirt off the surface. Rinse the cloth in water, wring it out well, and wipe off any remaining detergent. Allow to air dry in a well-ventilated, shaded place. **How to care for gel sheets:** If the adhesive surface of the gel sheet becomes dirty and stops adhering properly, it should be cared for according to the instructions: Apply a few drops of water to your fingertips and wipe the adhesive surface (gel side) to remove dirt. Shake off excess water and leave the sheet with the adhesive (gel) side up to dry in a well-ventilated, shaded place. Allow to air dry until there is no moisture on the surface, then attach the product to the storage sheet. Do not clean the sheets for too long or too often. This may damage the adhesive surface (gel side) and shorten its life.

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PRODUCENT: Shenzhen OSTO Medical Technology Company Limited

Adres siedziby: 2104, Unit4, Building2, Zhongxin Longsheng Square, NO. 219 Shenshan Road (Longgang Section), Longxin Community, Baolong Street, Longgang District, 518116 Shenzhen City, Guangdong Province, Chiny

Dodatkowy adres: Building 3, No. 3, Huitai South Road, Sandong Central Park, High-tech Industrial Park, Huicheng District, 516025 Huizhou City, Guangdong Province, Chiny

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Shenzhen OSTO Medical Technology Co., Ltd.	CE Technical File	File No	A014	Version	V1.4
	Declaration of Conformity	Effective Date	2025/02/24	Page	1 / 1

DECLARATION OF CONFORMITY

According to Art. 19 of Regulation (EU) 2017/745 on Medical Devices

MANUFACTURER: Shenzhen OSTO Medical Technology Company Limited

Main Site: 2104, Unit4, Building2, Zhongxin Longsheng Square, NO. 219 Shenshan Road (Longgang Section), Longxin Community, Baolong Street, Longgang District, 518116 Shenzhen City, Guangdong Province, China

Additional Site: Building 3, No. 3, Huitai South Road, Sandong Central Park, High-tech Industrial Park, Huicheng District, 516025 Huizhou City, Guangdong Province, China

Medical Device: Transcutaneous Electronic Nerve and Muscle Stimulator

Models: AST-303J

Classification - Annex IX: class IIa, rule 9

Conformity assessment Route: Annex IX of the Regulation (EU) 2017/745

Basic UDI-DI:695094080001GT

EMDN Code:Z12062801

SRN:CN-MF-000022960

We, Shenzhen OSTO Medical Technology Co., Ltd., here with declare that the stated medical devices meet the transposition into national law, the provisions of the Regulation (EU) 2017/745 on Medical Devices (MDR).And the DOC is issued under the sole responsibilities of the manufacturer.
All supporting documentation is retained at the premises of the manufacturer.

Intended use: This Device is intended to be used to pain relief

Standards applied: EN ISO14971, EN1041, EN 60601-1, IEC 60601-1-6, EN 60601-1-2, IEC 62304, EN 62366-1, IEC 60601-1-11, MEDDEV 2.7.1, EN ISO 10993-1.

Notified Body: Intertek Medical Notified Body AB
Torshamnsgatan 43, Box 1103 SE -164 22 Kista Sweden

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Shortened warranty terms

This Guarantee is provided by Armare S.A., a company governed by Polish law, with its registered office in Warsaw, ul. Grzybowska 87, 00-844, Poland. You can find us in the official Polish register of companies under the KRS number: 0000730590, NIP: 5213833614 and REGON: 380721860. We guarantee that the equipment you have purchased from us (hereinafter referred to as the "Product") will work properly. If something goes wrong for a reason inherent in the Product, this warranty protects you. The warranty starts from the day you receive the Product and lasts:

24 months, if you are a consumer or a natural person running a business and the purchase of the Product is not of a professional nature for you. Otherwise, the warranty period is 12 months.

If the law in your area requires a longer warranty, the warranty period will be extended to meet local law requirements.

In addition, if the Product does not comply with the contract, you have the right to legal protection from and at the expense of the seller. This warranty does not affect these rights - it does not limit, suspend or exclude them.

Detailed warranty conditions, how to submit complaints, what happens when you find a defect, what are exceptions to the warranty, when the warranty expires or is extended, and what are the rules regarding the processing of personal data, you can find on our website: <http://medivon.pl/warranty>

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EN | Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution. Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment. Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m² devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market after 13 August 2005.

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Shenzhen OSTO Medical Technology Co.,Ltd.
2104, Unit4, Building2,Zhongxin Longsheng Square,
NO.219 Shenshan Road(Longgang Section),Longxin
Community, Baolong Street, Longgang District, 518116
Shenzhen City, Guangdong Province, China
E-mail: info@sztosto.com

Phoenix Medtech GmbH
Koenigsberger Strasse,
11 64839 Muenster-Hessen,
Germany
Contact : Kevin Han
Tel : +49 6071 4977513
E-mail: contact@phoenix-medtech.de

Model: AST-303J
Basic UDI-DE: 695094080001GT

Date of preparation of the user manual: April 2026

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Importer: Armare S.A.