
medivon.

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



User manual **EN**

LED therapy lamp

Light therapy at home

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COMPREHENSIVE USER MANUAL (PANEL A300)

Control and Display

1. Light mode: click cyclically: RI (both), R (red), I (infrared), OFF (off).

2. Time (Timer): click cyclically: 5, 10, 15, 20, 25, 30 min or ON (continuous operation).

3. Digital display: shows currently selected mode and remaining time.

4. Brightness: click cyclically, select power: 25, 50, 75, 100%.

5. Pulsation: click cyclically, select speed: 10, 20, 40 Hz or ON (steady light).

6/7. LED Indicators: control lights inform which type of light is currently on.

Accessories: safety goggles (8), door hook (9), power cord (10), blocking system - height adjustment (11), connecting screws (12), mounting cable (13), connection cable (14).

Mounting Instructions:

A. Screw the threaded posts into the top of the lamp, tighten the cable bracket to the posts; if connecting multiple lamps, check the instructions part.

B. Hang the door hook on any standard size door.

C. Attach the cable bracket to the carabiner, hang the other end of the set on the door hook, use the regulator to set the height of the lamp.

D. To raise the lamp, pull the cable supporting the bottom of the device; to lower the lamp, press the lock lever, adjust the height to your needs.

How it works and how to use the therapy?

A300 series lamp supports muscle regeneration, reduces joint pain and improves skin appearance.

Red light (630 nm): acts on the skin surface and improves cell vitality.

Red light (660 nm): acts on the skin and collagen production.

Near infrared (810 nm): penetrates deep, soothing muscle fatigue and pain.

Near infrared (830 nm): supports tissue regeneration and reduces inflammation.

Near infrared (850 nm): reaches deeper - to muscles

and joints. Deep near infrared (1064 nm): supports deep tissue healing and chronic pain relief.

Note: infrared diodes are invisible to the eye, which is a natural phenomenon. Preparation: expose the area of the body, as clothing blocks the light. Optimal results are achieved at a distance of 15 cm to 20 cm from the panel.

Standard exposure time is 10 to 20 minutes per area (do not exceed 20 min in one session). Effects: consistency is key (1-2 times a day). Pain relief can occur quickly, but improvement in skin condition usually requires 8 to 12 weeks of regular treatments.

Target users: people aged 12 and over can use adult mode, starting with low intensity for adaptation, children under 12 must use under full adult supervision with session time halved; contraindications: infants 0-12 months cannot use, patients with epilepsy should avoid use as light can trigger seizures; users requiring caution should consult a doctor, this applies to patients with photosensitive diseases such as lupus or porphyria, pregnant women who should not illuminate the abdomen and lumbar region, cancer patients and people with severe diabetes or diabetic foot.

Residual risk description and control methods

The device is safe when used correctly according to the manual; possible risks are skin overheating, redness or dryness caused by too small distance, too long use or sensitive skin; maintain proper distance and stop use if irritation occurs; do not look directly at the diodes and use safety goggles.

The device's lifespan is about 5 years with normal use; it should be cleaned regularly. Before cleaning, disconnect the power and wait 30 minutes, clean with 75% alcohol or a soft cloth, do not allow liquid to enter the device, leave to dry, do not submerge or wash under pressure, store in a dry and cool place.

Operating parameters, intensity, time and frequency

Minimum effective intensity is at least 30 mW/cm^2 , maximum effective intensity up to 200 mW/cm^2 , minimum therapy time 10 minutes daily, maximum 30 minutes daily, minimum frequency 3 times a week, maximum 7 times a week and no more than once a day on the same area.

Unsafe radiation levels

Intensity above 200 mW/cm^2 , time over 30 minutes, more than one session daily on the same area and distance less than 3 inches can lead to skin burns.

Wavelength range is 630 nm to 1060 nm

Maximum treatment area A300 up to 900 cm^2

The LED panel should be held parallel to the body surface to ensure uniform illumination.

The panel should be held parallel at a distance of 8 cm to 30 cm.

Safety features. The A300 model has a cooling system and zero heat emission.

Stability and durability: LED diodes have a lifespan of about 100,000 hours, the device works stably at temperatures from 0 to 40 degrees in humidity from 10 to 75 percent, altitudes below 2000 meters, do not use in damp rooms.

Interactions with cosmetics and drugs: Caution should be exercised when using preparations containing capsaicin, salicylic acid or heating ingredients, as well as with drugs such as tretinoin, tetracycline, sulfonamides and other photosensitizing substances.

Avoid sun and UV lamps one hour before and after the treatment, after sunbathing wait 48 hours.

Required for people with photosensitive diseases, pregnant women, oncological patients and people with severe diabetes.

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EN | Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution. Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment. Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m² devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market after 13 August 2005.

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Model: A300
Basic UDI-DI: 697600641LEDLightK5

Date of preparation of the user manual: March 2026

Declaration of Conformity

Manufacturer:	Shenzhen Sungrow LED Technology Co., LTD 3F, 2th Factory Bldg(Phase II) Jiulonghai Sea & Tech Ind Park, Guangming Hightech Ind Park, Shenzhen, Guangdong Prov, China 518060
SIRN #:	CE-MP-000047979
European Representative:	Bernovis Sociedad Limitada Calle de Almansa 35, 1D, Madrid 28039 Spain Tel.: +34 658 396 730 SIRN # ES-AR-000001202
Trade Mark:	S&GOW
Product name:	LED Therapy Light
Model:	ESPL13460/ESPLUS760/ESPLUS1500/ESPRO340/ESPRO730/ESPRO1300 /ESPRO31000/ESPRO4800/ESRO/EP90/E1500/E1800/ E2400/ E4500/A300/ A900/ A1 200/P40MAX, MH160PRO/SL280 /B6-G422/BE105/ BE105P/ BE120/ BE120P/ 101/ 102 / BE400/ BE700 /D05
Basic UDI-DE:	697606641EDLgHMS
Intended Use:	LED Therapy Light is intended to provide Red and Infrared light to relieve pain in the body of adults. Red light is used to relieve skin surface pain. Infrared light provides topical heating for the purpose of elevating tissue temperature, for the temporary relief of minor muscle and joint pain, arthritis and muscle spasms, relieving stiffness, promoting the relaxation of muscle tissue and to temporarily increase local blood circulation. Red + infrared light is used to relieve both skin surface pain and muscle pain.

MDR Classification: Class I, According to Rule 1, Annex VIII of Regulation.

Conformity Assessment Route: EU DECLARATION OF CONFORMITY following the Annex II + Annex III +

Article 19 of MDR (EU) 2017/745.

We herewith under our sole responsibility declare that the above mentioned products meet the transposition into national law, the provisions of the following EU Regulation and Standards. All supporting documentations are retained under the premises of the manufacturer.

The manufacturer is exclusively responsible for the declaration of conformity. General applicable regulations, directives (Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.

Applied standards, common specification, guidance:

EN ISO 15223-1:2021,
EN ISO 20417:2021,
EN ISO 14971:2019,
IEC 62366-1:2015+A1:2020,
EN ISO 10993-1:2010,
EN ISO 10993-5:2009,
ISO 10993-10:2010,
EN ISO 13485:2016,

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Shortened warranty terms

This Guarantee is provided by Armare S.A., a company governed by Polish law, with its registered office in Warsaw, ul. Grzybowska 87, 00-844, Poland. You can find us in the official Polish register of companies under the KRS number: 0000730590, NIP: 5213833614 and REGON: 380721860. We guarantee that the equipment you have purchased from us (hereinafter referred to as the "Product") will work properly. If something goes wrong for a reason inherent in the Product, this warranty protects you. The warranty starts from the day you receive the Product and lasts: 24 months, if you are a consumer or a natural person running a business and the purchase of the Product is not of a professional nature for you. Otherwise, the warranty period is 12 months. If the law in your area requires a longer warranty, the warranty period will be extended to meet local law requirements.

In addition, if the Product does not comply with the contract, you have the right to legal protection from and at the expense of the seller. This warranty does not affect these rights - it does not limit, suspend or exclude them.

Detailed warranty conditions, how to submit complaints, what happens when you find a defect, what are exceptions to the warranty, when the warranty expires or is extended, and what are the rules regarding the processing of personal data, you can find on our website: <http://medivon.pl/warranty>

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Übersetzung erstellt durch den Importeur Armare S.A., ul. Grzybowska 87, 00-844 Warszawa, Polen
ARMARE S.A.