
medivon.

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
User manual EN

LED Therapy Light

Deep Tissue Recovery Series

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Product description and operation

LED therapy lamp – a medical device for phototherapy using 18 LED diodes with a total power of 54 W. The device generates therapeutic light (red and near infrared) that penetrates tissues and joints, supporting regenerative processes. The product is intended for phototherapy and is not a hyperthermic bulb. Thanks to its mobile design, it can be used in various home conditions. The device emits light that activates metabolic pathways, leading to increased ATP (cellular energy) production. Red light (660 nm) accelerates cellular metabolic reactions, supports regenerative processes and skin condition, while near infrared (850 nm) penetrates deeper into tissues, affecting muscles, joints and other structures. The effects include temporary relief of muscle and joint pain, reduction of post-exercise soreness, support of tissue regeneration, improvement of blood circulation and overall skin condition.

Remote control operation

- a – digital display – shows the current operating time settings and the selected brightness level of the device.
- b – timer (clock icon button) – switches the exposure duration; there are 6 levels to choose from: 5/10/15/20/25/30 minutes. The device will switch off automatically after the set time has elapsed.
- c – mode button (lightning icon) – used to switch between two operating modes: constant light (Constant Light Mode) and pulse mode at 10Hz (Pulse Flicker Mode).
- d – brightness adjustment (+/- button) – allows you to choose one of 5 brightness levels, from L0 (lowest) to L4 (highest).
- e – power button (Power icon) – a short press switches the device on or off completely.

Technical data and device parameters

Device lifetime: 30,000 hours, wavelength range: 660–850 nm, maximum treatment area: 200 cm², no additional therapeutic benefit above 200 mW/cm², stable operation at 0–40°C in dry conditions, operating temperature: 0–40°C, humidity: 10–75% RH. Recommended usage parameters: distance 7–15 cm, treatment time 5–25 minutes per area, frequency 1–2 times daily or 4–6 times weekly, treatment duration 3–6 weeks. Increased risk conditions: usage time over 30 minutes, more than 2 treatments per day on the same area, distance less than 7 cm – may cause skin burns.

Instructions for use and safety

Keep the device as parallel as possible to the skin surface to ensure even irradiation, maintain a distance of 7–15 cm and avoid overlapping treatment areas. Do not use in humid environments such as bathrooms or saunas, do not cover the light source during operation, do not use during sleep, as this may cause overheating or skin damage. Excessive use may cause redness, burning, dryness, allergic reactions and photosensitivity, therefore follow the recommended usage parameters. Wear protective glasses during use and do not look directly at the light source. Avoid sun exposure or UV radiation 1 hour before and after treatment, and if skin irritation occurs after sun exposure, wait 48 hours before using the device.

Assembly and installation

Installation using a mobile stand: 1. Unfold the mobile stand, 2. Set the desired height, 3. Mount the socket on the stand, 4. Screw in the lamp and adjust the tilt angle (up to 180°).

Cleaning, disinfection and maintenance

Disconnect power and wait 15 minutes for the device to cool down, wipe the housing with 75% alcohol or a soft cloth, allow to dry for 10 minutes before reuse, do not immerse, spray or pressure-wash the device. Store in a cool and dry place away from high temperature and humidity.

Medical notes and contraindications

The device is intended for home use and is safe when used as instructed, however skin overheating, slight redness or dryness may occur – especially at short distances, prolonged use or sensitive skin, as well as temporary eye discomfort when looking directly at the light source. Do not use on known or suspected tumors, as red light therapy may stimulate cell proliferation and potentially accelerate tumor growth or metastasis. Do not use on active infections. Use on open wounds only after initial healing or under medical supervision. Use caution with cosmetics and medications such as capsaicin, salicylic acid, retinoids, tetracyclines, sulfonamides and photosensitizing substances.

User groups and medical warnings

Children under 12 must use under adult supervision with reduced treatment time, users over 12 may use normally, starting with low intensity is recommended. Infants (0–12 months) must not use the device. Medical consultation is required for photosensitive diseases, pregnancy (no irradiation of abdomen and lower back), cancer conditions and severe diabetes.

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EN | Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution. Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment. Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m² devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market after 13 August 2005.

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Model: D05 | Marka: SGROW
Basic UDI-DI: 697600641LEDLightK5

Date of preparation of the user manual: March 2026

Declaration of Conformity

Manufacturer:	Shenzhen Sungrow LED Technology Co., LTD 3F, 2# Factory Bldg(Phase II) Bangkai Sci & Tech Ind Park, Guangming Hightech Ind Park, Shenzhen, Guangdong Prov. China. 518000
SRN #:	CN-MF-000047979
European Representative:	Riomaviv Sociedad Limitada Calle de Almenara 35, 1D, Madrid 28039 Spain TEL: +34 638 398 230 SRN # ES-AK-000001202
Trade Mark:	SGROW
Product name:	LED Therapy Light
Model:	ESPLUS3607/ESPLUS750/ESPLUS15007/ESPLUS3000/ESPRO3607/ESPRO750/ESPRO15007/ESPRO3000/ESPRO4800/ES1007/ES15007/ES18007/ES24007/ES42007/AS3007/AS6007/AL24007/AL42007, ESOPRO/SLE807/BE-G127/BE1037/BE10397/BE1207/BE12097/LO37/LO87/BE1039/BE1209/DO3
Basic UDI-DE:	697600641EED1GHTK5
Intended Use:	LED Therapy Light is intended to provide Red and Infrared light to relieve pain in the body of adults. Red Light is used to relieve skin surface pain, Infrared light provides topical heating for the purpose of elevating tissue temperature, for the temporary relief of minor muscle and joint pain, arthritis and muscle spasm, relieving stiffness, promoting the relaxation of muscle tissue, and to temporarily increase local blood circulation. Red + Infrared light is used to relieve both skin surface pain and muscle pain.

MDR Classification: Class I, According to Rule 1, Annex VIII of Regulation.

Conformity Assessment Route: EU DECLARATION OF CONFORMITY following the Annex III +

Article 19 of MDR (EU) 2017/745.

We herewith under our sole responsibility declare that the above mentioned products meet the transposition into national law, the provisions of the following EU Regulation and Standards. All supporting documentations are retained under the premises of the manufacturer.

The manufacturer is exclusively responsible for the declaration of conformity general applicable regulations,

directives: Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.

Applied standards, common specification, guidance:

EN ISO 15223-1:2021,
EN ISO 20417:2021,
EN ISO 14971:2019,
IEC 62366-3:2015+A1:2020,
EN ISO 10993-1:2010,
EN ISO 10993-5:2009,
ISO 10993-10:2010,
EN ISO 13485:2016,

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Shortened warranty terms

This Guarantee is provided by Armare S.A., a company governed by Polish law, with its registered office in Warsaw, ul. Grzybowska 87, 00-844, Poland. You can find us in the official Polish register of companies under the KRS number: 0000730590, NIP: 5213833614 and REGON: 380721860. We guarantee that the equipment you have purchased from us (hereinafter referred to as the "Product") will work properly. If something goes wrong for a reason inherent in the Product, this warranty protects you. The warranty starts from the day you receive the Product and lasts: 24 months, if you are a consumer or a natural person running a business and the purchase of the Product is not of a professional nature for you. Otherwise, the warranty period is 12 months. If the law in your area requires a longer warranty, the warranty period will be extended to meet local law requirements.

In addition, if the Product does not comply with the contract, you have the right to legal protection from and at the expense of the seller. This warranty does not affect these rights - it does not limit, suspend or exclude them.

Detailed warranty conditions, how to submit complaints, what happens when you find a defect, what are exceptions to the warranty, when the warranty expires or is extended, and what are the rules regarding the processing of personal data, you can find on our website: <http://medivon.pl/warranty>

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ARMARE S.A.