
Dr. Mayer

Rechargeable

Sonic Power Toothbrush

Operating Instructions

Model: GTS2065

READ ALL INSTRUCTIONS BEFORE USING

IMPORTANT SAFEGUARDS

The Power toothbrush is safe and effective when used according to directions. Like any electrical appliance, however, basic safety precautions should be followed, especially when children are present.

DANGERS-WATER

To reduce the risk of burns, electrocution, fire, or injury to persons:

- Do not place or store the charger base where it can fall or be pulled into a bathtub, sink or toilet.
- Do not place or drop the charger into water or other liquid.
- Do not reach for a product that has fallen into water or other liquid. Unplug immediately.
- Do not use while bathing.
- Unplug charger base before cleaning and assure that it is dry before plugging it into an electrical outlet.

DANGERS-VOLTAGE

To reduce the risk of burns, electrocution, fire, or injury to persons:

- Do not plug the charger base into an outlet with a voltage other than specified on the bottom of the charger.
- This product is designed to operate within a range of 100-240 volts.
- Voltage converters do not guarantee compatibility.
- Do not force the plug into an outlet.
- If the plug does not fit easily into the outlet, discontinue use.

DANGER-WARNINGS

To reduce the risk of burns, electrocution, fire, or injury to persons:

- This product is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities or lack of experience and knowledge, unless they have been given adequate supervision or instruction concerning use of the product by an adult responsible for their safety.
- Children should be supervised to ensure they do not play with the product.
- Use this product only for its intended use as described in this manual. This product is designed to clean your teeth, gums and tongue only. Do not use any attachments or brush heads other than those recommended by the manufacturer.
- The toothbrush is a personal care appliance and is not intended for use on multiple patients in a dental office or institution.
- If your toothpaste includes peroxide, baking soda or other bicarbonate (common in whitening toothpastes), thoroughly clean your brush head and handle to avoid cracking the plastic.
- Never operate this product if it has a damaged cord, plug or charger base.
- Do not use this product if it is not working properly, if it has been dropped or damaged or dropped into water. Immediately discontinue use if the product appears damaged in any way.

- Keep the toothbrush, charger base and power cord away from heated surfaces.
- If you experience discomfort or pain in your mouth, discontinue use of this product and contact a dentist.
- Do not attempt to remove the battery.
- Do not clean the toothbrush in the dishwasher.
- This appliance has a polarized plug (one blade is wider than the other). As a safety feature, this plug will fit in a polarized outlet only one way. If the plug does not fit fully in the outlet, reverse the plug. If it still does not fit, contact a qualified electrician. Do not attempt to defeat this safety feature.
- Do not use outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered.
- Toothbrushes can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Cleaning and user maintenance shall not be made by children without supervision.
- Children shall not play with the appliance.
- The battery must be removed from the appliance before it is scrapped.
- The appliance must be disconnected from the supply mains when removing the battery.
- The battery is to be disposed of safely.
- If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.
- Regarding how to clean the tooth brush surfaces in contact with hand and skin, please refer to section <MAINTENANCE>.

MEDICAL-WARNINGS



- If you have had oral surgery or gum surgery within the past two months, consult your dental provider before using the product.
- Consult your dental provider if excessive bleeding occurs while using this product or minor bleeding continues to occur after one week of use.
- Consult your physician or implant device manufacturer prior to using this product if you have questions concerning the combined use of this product with a pacemaker or other implanted medical device(s).
- Consult your physician prior to using the toothbrush if you have any medical concerns.
- The toothbrush has been tested and is certified as per FCC VOC and CE regulations.
- Always unplug this product immediately after using.
- Never use while sleeping or drowsy.

SAVE THESE INSTRUCTIONS

GETTING STARTED

Charging the toothbrush

Be sure to turn off the toothbrush before charging. plug the power plug into the outlet. Place the handle in the charger. The charger indicator light will be lit, indicating the handle is charging. After the battery is fully charged, charge indicator light will turn solid green, indicating the toothbrush is fully charged.

Notes:

- 1 It takes at least 24 hours to charge the battery fully.
- 2 The toothbrush must be charged for 24 hours, if:
 - a. this is the first time using the toothbrush;
 - b. the toothbrush has not been used for more than three months;
- 3 To keep the battery fully charged, we recommend that you keep the toothbrush in the charger when not in use.
- 4 During charging, the toothbrush will get warm, but this is normal.
- 5 The operating duration will become shorter as the battery gets older. If the toothbrush operates only for 1 short time after being fully charged, this means that the battery is nearing the end of its life cycle. Contact service center for support.

HOW TO USE

- 1 Slide the brush head onto the handle. Push brush head down until it locks into place on handle.

Note: There is small gap between color ring and the handle.(Fig 1)


Fig 1

- 2 Wet the bristles and apply a small amount of toothpaste.
- 3 Place the bristles against the teeth at a slight angle towards the gum line.(Fig 2)
- 4 Press the mode button to select mode. If you know the mode used last time, you need not select mode again. There are 4 modes.


Fig 2

- a. **CLEAN:** Standard mode for superior teeth cleaning.
- b. **SOFT:** Gently caress, carefully clean sensitive mouth.
- c. **WHITEN:** Teeth cleaning with polishing and whitening.
- d. **MASSAGE:** Pulsates for gentle gum stimulation.

Notes:

During brushing, if press Mode button, toothbrush will work at next mode and start another 2 minutes cycle and the mode indicator will light up.

- 5 Apply light pressure and gently move the brush head slowly across the teeth in a small back and forth motion.
- 6 In order to maximize cleaning effectiveness, divide your mouth in 4 sections. Start with Section 1 and brush each section for 30 seconds. The brushing action will pause between each 30 second interval, indicating when to move on to the next section. The toothbrush will automatically stop after 2 minutes of use. (Fig 3)
- 7 Pull brush head out from the toothbrush.(Fig 4)
- 8 Rinse both the brush head and handle under running water wipe them dry. (Fig 5)
- 9 Put the handle back into the charger.


Fig 3


Fig 4



Fig 5

MAINTENANCE

Be sure to remove the plug from the electrical outlet before cleaning the product.


Toothbrush

- 1 Remove the brush head.
- 2 Wipe all surface areas with a damp cloth and place the handle back into the charger. Do not submerge the handle into water.
- 3 Rinse the brush head and bristles after each use.
- 4 Do not clean in the dishwasher.

Charger

- 1 Using a soft cloth, wipe off any dirt on the charger and the storage.
- 2 Clean the plug periodically to ensure it remains in good working order. The accumulation of dust or dirt on the plug can cause an electric shock, overheat or fire.

TROUBLESHOOTING


Problem	Possible cause	Solution
The toothbrush doesn't work.	You have just purchased the product or it has not been used for more than 3 months.	Charge the main unit for at least 24 hours.
Weak or diminished vibration.	The toothbrush is pressed hard against your teeth.	Place the brush head gently against your teeth.
The toothbrush operates only for a short period of time even immediately after being charged.	The battery has reached the end of its life cycle.	Contact service center for support.
	Charging time was insufficient.	Charge the main unit for at least 24 hours.

DISPOSAL

The power source of this product is Lithium-ion battery which can be recycled. At the end of the product's life, always remove the battery before you discard the product, and recycle or dispose of the battery according to local waste management requirements. If you have recycling questions, contact your local waste management office.




Removing and recycling the battery

Do not disassemble the product at any time other than when disposing of it. Doing so earlier will destroy the product. Observe basic safety precautions.

The appliance must be disconnected from the supply mains when removing the battery.

- 1 Before removing the battery, be sure to completely drain the unit of power.
- 2 Use a screwdriver to remove the screw cover and battery door.
- 3 Hold the handle upside down, press the toothbrush shaft, and remove the two hooks to release the internal components of the handle.
- 4 Insert the screwdriver next to the battery connections, and twist to break the connections. Remove the circuit board and pry the battery from the plastic carrier.
- 5 Recycle or dispose of the battery properly. Ask your local waste management office if they accept the remaining toothbrush components for recycling.

Be careful not to short-circuit the negative and positive poles of the battery. Do not connect the appliance to the mains again after you have removed the battery.

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

SI

Portable Water Flosser WT3500

Read this manual thoroughly before using and save it
for future reference

IMPORTANT SAFETY INSTRUCTIONS

DANGER

1. Never use this appliance unless understand the operation process and safety description to avoid being hurt by improper operation.
2. Only operate this appliance by following the operating way specified in this instruction manual.
3. Do not disassemble this appliance or change its components to avoid this appliance being damaged or causing electric shock.
4. Do not fill the water tank with water exceeding 40°C to do oral cleaning to avoid injury.
5. Please unplug and then remove the appliance from the water if it drops into water when being charged.
6. Please replace the charger if the power cord or plug is damaged.
7. Do not heat the appliance with fire and do not charge, use or dispose it at high-temperature places.

WARNING

- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.

CAUTION

1. Use this product only for its intended use as described in this


manual.

2. If you have had oral or gum surgery in the previous 2 months, consult your dentist before you use this appliance.
3. Do not press gums or teeth by the nozzle with excessive force to avoid the gums being hurt.
4. Do not immerse the appliance into water to avoid malfunction or causing fire, electric shock and explosion.
5. Do not plug or unplug with wet hands.
6. Make sure that your outlet voltage corresponds to the voltage stated on the rating label.
7. Do not operate any appliance with a damaged charger, cord or plug or after the appliance malfunctions, or has been damaged in any manner or dropped into water. Return appliance to an authorized service shop for examination, repair or adjustment.
8. If the supply cord is damaged, it must be replaced by the manufacturer, service shops or qualified persons in order to avoid a hazard.
9. Do not clean the nozzle, the handle, the charger, the water reservoir and its cover in the dishwasher or a microwave.
10. Do not use attachments other than the ones recommended by the manufacturer.
11. Unplug the appliance before perform maintenance.
12. Do not use the charger outdoors or near heated surfaces.

NOTE


1. Please do not use the appliance if without nozzle being attached.
2. Do not fill the reservoir with mouthwash which may corrode the appliance. Always clean the reservoir with clean water after using with mouthwash!
3. Please consult the doctor whether you can use the appliance if you have had dentistry surgical operation in past two months or are under being subject to dental treatment.
4. Save these instructions.

KNOW YOUR WIRELESS WATER FLOSSER


OPERATION INSTRUCTION

1. Check and make sure that your hands are dry before operate, then connect the appliance with power source.



2. Insert the body into the charger, the charging indicator blinks and it will keep blinking during charging. When the water flosser is charged well, the charging indicator will light solidly. Charge the appliance for 24 hours before first use.


NOTE: The water flosser needs to be charged if the charging indicator blinks before or after operating.


3. Remove the water flosser from the charger, then insert the nozzle into the body. And you can press the nozzle release button if need to remove the nozzle.

4. Load the water tank with water by following the two ways as below:

or

5. Select your required mode by the “Mode button”. There are three modes - “Normal”, “Soft” and “Pulse” - you can press the “Mode button” to make switch between modes during operating.


A. Normal: For the normal cleaning of teeth.

B. Soft: For the sensitive teeth and the user who uses the water flosser for the first time or only needs soft operation.

C. Pulse: For massaging gums to promote blood circulation and strengthen the teeth.

NOTE: If you use this water flosser for the first time, please choose the soft mode for a week, when you adapt to it, you can use the other mode you like.



6. Hold the water flosser vertically and then insert the nozzle into the mouth and aim at gums or teeth, at the same time, slightly open your mouth to let the water flow out smoothly.


7. Turn on the switch to activate the water flosser. The followings shall be paid attention to during operation.

A. The water shall flow out in a direction vertical to gums or teeth.


B. Slowly move the water flosser along the teeth and clean them in order. Do not spray water along the periodontal pocket directly.


8. After operation, press the switch to turn off the appliance.
The appliance will shut off automatically after operating for 2 minutes. Please press the switch again if need to continue operating.
9. After finish operating, you can follow the operation as below to remove the residual water:
- A. Open the water tank door and empty the water tank.

- B. Press the switch, the water flosser starts working to drain out the residual water inside automatically.

- C. Press the switch again to turn off the water flosser, and then close the water tank door.
- D. Wipe the water flosser with dry cloth.

NOTE: The water flosser shall be free of residual water to avoid bacterial infection. Therefore, if it is not used for a long time, please clean the water flosser well, then dry and store it.

CLEANING AND MAINTENANCE

Please clean the water flosser with clean water or neutral detergent, but never use corrosive or abrasive detergent as it will damage the water flosser, and do not wash it with water higher than 50°C.

1. Clean the body

- A. Please clean the body with dry or wet cloth.
- B. Do not immerse the body into water for cleaning.

2. Clean the water tank

- A. Wipe the water tank with dry or wet cloth, or detach it and rinse it with water.
- B. Please empty the water tank completely if it is not used for a long time.

3. Clean the nozzle

- A. Wipe the nozzle with dry or wet cloth, or detach it and rinse it with water.
- B. Do not bend, pull and warp the nozzle.

C. The dentist suggests replacing the nozzle every 6 months.

4. Clean the charger

Wipe the metal part of plug with dry cloth every 6 months to avoid dust accumulation which may result in poor contact.

TROUBLESHOOTING		
Problem	Possible cause	Solution
The water flosser can not work.	The water flosser is used for the first time or it has been not used for more than 3 months.	Please charge the water flosser for above 24 hours.
Water pressure is low when working.	There is something blocked in the nozzle. The nozzle is deformed.	Replace the nozzle
The water flosser can only be operated for short time after being charged fully.	The filter is blocked.	Clean the filter
	The life of battery ends.	Please contact service department
	The water flosser has not charged well.	Please charge the water flosser for more than 24 hours before operation.
The charging indicator does not illuminate when charging.	The power plug has not been connected well.	Please connect the power source well.
	Circuit malfunction	Please contact service department

WASTE DISPOSAL


The water flosser adopts recyclable rechargeable battery. When the life of water flosser ends, remove the battery and then recycle or dispose it basing on the local waste disposal

regulation. Please consult the local waste disposal authority if you have any questions.


DETACH RECHARGEABLE BATTERY

Do not dismantle the water flosser unless it needs to be discarded, or it will be damaged. Please unplug the appliance before remove the battery.


1. Make sure that the battery has been exhausted before remove it. To run out the battery, press the switch to activate the water flosser and make it run until it stops working automatically. Repeat operating for several times until the battery exhausts.
2. Remove the water tank and then unscrew the two screws at the bottom cover with a Phillips screwdriver.

3. Remove the two screws at the top cover with a Phillips screwdriver, then detach the top cover and remove the unit core.
4. Unscrew the screws at the bracket of unit core with a Phillips screwdriver.


5. Shear off the connecting wire of battery with scissors, then remove the battery.

6. Recycle or dispose the battery correctly. And please consult the local waste disposal authority for the recycle of other components of water flosser.

NOTE: Do not short-circuit the battery and never energize the water flosser again after remove the battery.


ENVIRONMENT FRIENDLY DISPOSAL

You can help protect the environment!
Please remember to respect the local regulations: hand in the non-working electrical equipments to an appropriate waste disposal center.

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