RED LIGHT THERAPY PANEL COMFORT

USER MANUAL

PROFESSIONAL TREATMENT AT HOME

DEAR CUSTOMER,

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EN

PALSAR7 - ABOUT THE BRAND

PALSAR7 IS A BRAND SPECIALIZING IN INNOVATIVE PRODUCTS FOR PROFESSIONAL SKIN CARE FROM THE COMFORT OF YOUR HOME. IT OFFERS MODERN TECHNOLOGIES THAT SUPPORT SKIN HEALTH AND BEAUTY WITHOUT THE NEED TO VISIT A BEAUTY SALON. THE BRAND EMPHASIZES QUALITY, SAFETY AND EASE OF USE OF ITS PRODUCTS. THANKS TO THIS, YOU CAN INDULGE IN PROFESSIONAL CARE ANYTIME, ANYWHERE.

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
EN

ABOUT THE PRODUCT

The red and near-infrared light panel is a modern home device designed for non-invasive care of the skin, scalp, and frequently stressed areas of the body. It uses red light with a wavelength of 660 nm and near-infrared light with a wavelength of 850 nm. Red light primarily affects the surface layers of the skin and can help give it a smoother and firmer appearance. Near-infrared light helps create a pleasant feeling of warmth and relaxation after physical activity or prolonged sitting or standing. Thanks to its larger coverage area, the panel is particularly suitable for the back, shoulders, arms, legs, or chest. However, you can also use it on your face, décolletage, or scalp. The long design with three adjustable joints allows you to easily rotate the panel and direct the light exactly where it's needed.

PRODUCT DESCRIPTION



INSTRUCTIONS FOR USE


1. Connect pole parts 1. and 2.. Make sure both parts are properly inserted and secured. Incorrect connection may make it difficult or impossible to separate them later.
2. Insert the assembled pole into the base and secure it by turning it clockwise.
3. Attach the telescopic arm to the top of the pole and push it down until it is firmly seated.
4. Attach the light fixture to the telescopic arm and secure it by tightening the fastening element.

RECOMMENDATIONS:

THE RECOMMENDED DURATION OF EACH USE IS 30 MINUTES, 3-5 TIMES A WEEK.

1. Insert the DC end of the cable into the panel.
2. Plug the adapter into the socket.
3. Press the button to turn on the light.

PRODUCT FEATURES

The panel combines red and infrared light. These technologies work together to enhance your at-home relaxation routine for joints and sore muscles.

RED LIGHT (660 nm)

Red light is a modern technology for body and skin care. It acts directly within cells, where it naturally helps trigger healing processes. At the same time, it has been scientifically proven to effectively reduce inflammation.

INFRARED LIGHT (850 nm)

Infrared light is invisible to the human eye, but it is highly popular in body care, relaxation, and thermal comfort. It deeply warms the treated area, aiding in the body's natural regeneration and relaxation. The heating function also pleasantly relieves muscle tension, contributes to greater comfort in the joints, and prepares the body for subsequent care.

RESULTS

MORE ENERGY FOR CELLS

These technologies support mitochondrial function and cellular energy.

POST-WORKOUT RECOVERY

Ideal after training, physical activity, or a long day.

COMFORT FOR MUSCLES, JOINTS, AND THE BACK

Infrared light targets deeper, strained areas of the body.

HEALTHIER SKIN APPEARANCE


Supports collagen, elasticity, and vitality, while helping to reduce the appearance of fine lines, wrinkles, and scars.

SCALP & HAIR REVIVAL

Also suitable as part of a regular beauty and self-care routine.

EVENING RELAXATION RITUAL

Soothing light for relaxation, calm, and at-home wellness care.

EFFECTIVE

104 diodes ensure highly effective therapy.

EASY TO USE

Easy to use for everyone.

CLEANING AND MAINTENANCE

- Always turn off the device and unplug it from the power outlet after use.
- Do not leave the device unattended while it is plugged in.
- The device is not waterproof. Do not allow water or other liquids to enter any part of the device.
- Do not use chemical cleaning agents, alcohol, or solvents.
- Store the device in a dry, cool, and shaded place, away from direct sunlight and high temperatures, and out of the reach of children.
- Do not expose the device to impacts, high temperatures, fire, or mechanical stress to prevent deformation or damage.
- Use the device according to the instructions, do not exceed the recommended usage time, and avoid excessive use.
- Do not disassemble the device or attempt to repair it yourself.
- In case of damage or a warranty claim, contact the manufacturer or the store retailer where you purchased the product.

ANTI SLIP SILICONE

PANEL CAN BE PLACED
VERTICALLY OR HORIZONTALLY

WARNING

THIS PRODUCT IS NOT INTENDED FOR:

- children and pregnant women
- people with pacemakers or other electronic medical implants
- people with reduced sensitivity to heat, touch, or pain
- individuals with reduced physical, sensory, or mental capabilities, unless they are under the supervision of a person responsible for their safety, people with skin conditions, allergic reactions, dermatitis, eczema, severe sunburn, or irritated, injured, swollen, inflamed, or significantly painful skin at the application site
- people with suspected fractures, sprains, muscle strains or tears, thrombosis, or other acute medical conditions
- people with epilepsy; certain lights or frequencies may trigger a seizure
- if you have a long-term medical condition, have suffered an injury, are recovering from surgery, are experiencing severe pain, or feel uncertain, consult a doctor before use
- direct viewing of the light source, especially for extended periods
- when using the product near the face, protective eyewear is recommended

NOTICE

- Read these instructions carefully before first use and keep them for future reference.
- Use the product only in accordance with the instructions provided in this manual.
- This product is not a medical device and is not intended to diagnose, treat, alleviate, or prevent diseases, injuries, or other medical conditions.
- The device is not suitable for people with increased sensitivity to light. If in doubt, consult your doctor before use.
- Do not use this device if you are taking medications that may cause your skin to become sensitive to light.
- Do not touch the device, adapter, or power cable with wet hands.
- Do not use the device while walking, driving, or engaging in sports activities.
- Stop using the device if you experience discomfort, excessive heat, burning, or if redness or skin discoloration persists for more than 24 hours. In such cases, consult a dermatologist or physician.
- Do not use the device if it is visibly damaged or if the power cable, adapter, or plug is damaged.
- Use the device only with the supplied adapter or an adapter recommended by the manufacturer.
- Remove metal jewelry and accessories from the treatment area before use. Do not use the device over clothing with metal zippers, snaps, buttons, or other metal elements.

TECHNICAL SPECIFICATIONS:

Product size: **38 × 22 cm**

Net weight: **0,5 kg**

Red LED Light: **660 nm**

Infrared LED Light: **850 nm**

Number of LED Lights: **104 pcs**

Power: **24 W**

THE PACKAGE CONTAINS:

- 1× LED Panel
- 1× Telescopic rod
- 1× Base
- 1× Adjustment bracket
- 1× Bracket
- 1× Adapter
- 1× Controller
- 1× Goggles

THE IMAGES ARE FOR ILLUSTRATIVE PURPOSES ONLY. MADE IN CHINA.

FOR MORE INFORMATION OR TO DOWNLOAD MORE DETAILED DATA, PLEASE CONTACT YOUR DEALER.

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Wear and tear of the product and its components caused by their common use (discharged batteries, loss of the original color of the product, decreased capacity of accumulator, wear and tear of mechanical or displaying components, wear and tear of supply cables, etc.) shall be subject neither to defect warranty nor quality warranty in the moment of take-over.

Furthermore, complaints about product parameters, not included in the operating instructions or in any commercial-technical documentation of Seller, related to the product, shall not be subject to quality warranty in the moment of take-over.

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