
RED LIGHT THERAPY PANEL CORE

USER MANUAL

PROFESSIONAL TREATMENT AT HOME

DEAR CUSTOMER,

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EN

PALSAR7 - ABOUT THE BRAND

PALSAR7 IS A BRAND SPECIALIZING IN INNOVATIVE PRODUCTS FOR PROFESSIONAL SKIN CARE FROM THE COMFORT OF YOUR HOME. IT OFFERS MODERN TECHNOLOGIES THAT SUPPORT SKIN HEALTH AND BEAUTY WITHOUT THE NEED TO VISIT A BEAUTY SALON. THE BRAND EMPHASIZES QUALITY, SAFETY AND EASE OF USE OF ITS PRODUCTS. THANKS TO THIS, YOU CAN INDULGE IN PROFESSIONAL CARE ANYTIME, ANYWHERE.

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
ABOUT THE PRODUCT

The red and near-infrared light therapy panel brings professional red light therapy into the comfort of your home. It combines 660 nm and 850 nm wavelengths, which have been extensively studied in relation to cellular activity and the body's natural regenerative processes. It supports the appearance of smoother, firmer, and more even-looking skin. At the same time, it helps maintain the vitality of the scalp and promotes a pleasant feeling of relaxation. Thanks to its high power output and large treatment area, it is suitable for use on the back, shoulders, legs, and smaller areas of the body. The intuitive controls allow you to quickly adjust the light mode according to your preferences. The panel can easily become part of your daily routine for beauty care, recovery, and overall well being.

PRODUCT DESCRIPTION


INSTRUCTIONS FOR USE


1. Connect the power cable to the power input on the back of the panel and plug the other end into an electrical outlet.
2. A **POWER ON/OFF** button is located on the back of the panel. Press it to turn on the device.
3. Use the **MODE** button to select either the 660 nm red light mode, the 850 nm infrared light mode, or the combined mode using both wavelengths simultaneously.
4. Use the **TIMER** button to set the treatment duration. Available time settings are 5, 10, 15, 20, 25, and 30 minutes.
5. The **PULSE** button allows you to activate pulse mode and select a frequency of 10 Hz or 40 Hz.
6. After selecting the desired mode and treatment time, the device will operate according to the chosen settings.
7. Once the selected time has elapsed, the therapy session will automatically end.

IRRADIANCE AND COVERAGE

RED LIGHT 660 NM

INFRARED LIGHT 850 NM

RECOMMENDATIONS:

THE RECOMMENDED DURATION OF EACH USE IS 30 MINUTES, 3-5 TIMES A WEEK.

PANEL CONTROLS

PULSE LIGHT

The advantage of pulsed light is that, thanks to the rapid alternation of light pulses, light energy can be utilized more efficiently and in a more targeted manner.

1. For best results, position the panel approximately 15 cm from the treatment area.
2. The recommended duration of a single therapy session is 10 to 30 minutes per body area.
3. The device can be used daily. However, do not exceed 30 minutes of treatment on the same area during a single session.

PRODUCT FEATURES

The panel combines red and infrared light. These technologies work together to enhance your at-home relaxation routine for joints and sore muscles.

RED LIGHT (660 nm)

Red light is a modern technology for body and skin care. It acts directly within cells, where it naturally helps trigger healing processes. At the same time, it has been scientifically proven to effectively reduce inflammation.

INFRARED LIGHT (850 nm)

Infrared light is invisible to the human eye, but it is highly popular in body care, relaxation, and thermal comfort. It deeply warms the treated area, aiding in the body's natural regeneration and relaxation. The heating function also pleasantly relieves muscle tension, contributes to greater comfort in the joints, and prepares the body for subsequent care.

RESULTS

MORE ENERGY FOR CELLS

These technologies support mitochondrial function and cellular energy.

POST-WORKOUT RECOVERY

Ideal after training, physical activity, or a long day.

COMFORT FOR MUSCLES, JOINTS, AND THE BACK

Infrared light targets deeper, strained areas of the body.

HEALTHIER SKIN APPEARANCE


Supports collagen, elasticity, and vitality, while helping to reduce the appearance of fine lines, wrinkles, and scars.

SCALP & HAIR REVIVAL

Also suitable as part of a regular beauty and self-care routine.

EVENING RELAXATION RITUAL

Soothing light for relaxation, calm, and at-home wellness care.

EFFECTIVE

60 diodes ensure highly effective therapy

EASY TO USE

Easy to use for everyone

CLEANING AND MAINTENANCE

- Always turn off the device and unplug it from the power outlet after use.
- Do not leave the device unattended while it is plugged in.
- The device is not waterproof. Do not allow water or other liquids to enter any part of the device.
- Do not use chemical cleaning agents, alcohol, or solvents.
- Store the device in a dry, cool, and shaded place, away from direct sunlight, high temperatures, and out of the reach of children.
- Do not expose the device to impacts, high temperatures, fire, or mechanical stress to prevent deformation or damage.
- Use the device according to the instructions, do not exceed the recommended usage time, and avoid excessive use.
- Do not disassemble the device or attempt to repair it yourself.
- In case of damage or a warranty claim, contact the manufacturer or the retailer where you purchased the product.

WARNING

THIS PRODUCT IS NOT INTENDED FOR:

- infants and pregnant women
- children over 12 years of age may use the device in the same manner as adults
- for children under 12 years of age, we recommend reducing the treatment time to approximately half the recommended adult dose
- people with pacemakers or other electronic medical implants
- people with reduced sensitivity to heat, touch, or pain
- people with reduced physical, sensory, or mental capabilities, unless they are under the supervision of a person responsible for their safety
- people with skin conditions, allergic reactions, dermatitis, eczema, severe sunburn, or irritated, injured, swollen, inflamed, or significantly painful skin at the application site
- people with suspected fractures, sprains, muscle strains or tears, thrombosis, or other acute medical conditions
- people with epilepsy; certain lights or frequencies may trigger a seizure
- if you have long-term health issues, an injury, are recovering from surgery, are experiencing severe pain, or are unsure, consult a doctor before use
- do not look directly into the light, especially for extended periods
- when using near the face, we recommend wearing protective eyewear

NOTICE

- Read these instructions carefully before first use and keep them for future reference.
- Use the product only in accordance with the instructions provided in this manual.
- This product is not a medical device and is not intended to diagnose, treat, alleviate, or prevent diseases, injuries, or other medical conditions.
- The device is not suitable for people with increased sensitivity to light. If in doubt, consult your doctor before use.
- Do not use this device if you are taking medications that may cause your skin to become sensitive to light.
- Do not touch the device, adapter, or power cable with wet hands.
- Do not use the device while walking, driving, or engaging in sports activities.
- Stop using the device if you experience discomfort, excessive heat, burning, or if redness or skin discoloration persists for more than 24 hours. In such cases, consult a dermatologist or physician.
- Do not use the device if it is visibly damaged or if the power cable, adapter, or plug is damaged.
- Use the device only with the supplied adapter or an adapter recommended by the manufacturer.
- Remove metal jewelry and accessories from the treatment area before use. Do not use the device over clothing with metal zippers, snaps, buttons, or other metal elements.

TECHNICAL SPECIFICATIONS:

Power: **120 W**
Input voltage: **100 to 240 V**
Lifespan: **20 000 hours**
LED quantity: **60 pcs**
Red LED light: **660 nm**
Infrared LED light: **850 nm**
Weight: **1,6 kg**

THE PACKAGE CONTAINS:

- 1× LED panel
- 1× Power cord
- 1× Protective goggles
- 1x User manual

THE IMAGES ARE FOR ILLUSTRATIVE PURPOSES ONLY. MADE IN CHINA.

FOR MORE INFORMATION OR TO DOWNLOAD MORE DETAILED DATA, PLEASE CONTACT YOUR DEALER.

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Wear and tear of the product and its components caused by their common use (discharged batteries, loss of the original color of the product, decreased capacity of accumulator, wear and tear of mechanical or displaying components, wear and tear of supply cables, etc.) shall be subject neither to defect warranty nor quality warranty in the moment of take-over.

Furthermore, complaints about product parameters, not included in the operating instructions or in any commercial-technical documentation of Seller, related to the product, shall not be subject to quality warranty in the moment of take-over.

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