
INSTRUCTIONS

Teeth whitening kit


English	2
Deutsch	7
Français	12
Čeština	17
Polski	22
Slovenčina	27

The SMILI teeth whitening kit combines unique light technology with a shaped mouthpiece to effectively remove stains from the tooth surface.

The mouthpiece is designed to fit the anatomy of the human mouth, making it easy to use and extremely comfortable. It has been designed to fit perfectly with the shape of your teeth, minimising irritation to the mouth and allowing maximum comfort.

Before starting the procedure, it is very important to familiarize yourself with the whitening process and understand it completely. To achieve optimal results, please read the instructions and safety guidelines.


Package contains


LED light device

Instructions for use

Whitening results guide

Pens

Whitening in 6 steps

Rinse your mouth with water and rinse the mouthpiece.

Apply the gel to the surface of the teeth using the pen with whitening gel.

Place the mouthpiece into your mouth and bite down. Then turn on the light.

15 minute treatment (blue mode)
10 minute treatment (blue-red mode)

Rinse your mouth with water.

Rinse the mouthpiece in water.

How to use

1. Before use, rinse your mouth briefly and rinse the mouthpiece.
2. Compare the shade of your teeth in front of the mirror with our whitening results guide and note the stage you are in.
3. Rotate the whitening gel pen until there is enough gel on the brush. Apply a thin, even layer to the surface of the teeth, avoid direct contact with your gums. If you feel slight tingling, don't worry. The gel has gotten on your gums, so rinse them carefully with water. In case of severe burning, stop the whitening process.
4. Before use, make sure the device with the LED light is fully charged.
5. Insert the mouthpiece into your mouth and bite down, then turn on the LED light. Take deep and regular breaths to prevent excessive saliva production. (Excessive saliva can dilute the gel concentration.)
6. The whitening device switches off automatically after 15 minutes. Remove the mouthpiece and rinse it. Continue by rinsing your mouth. Compare your teeth to the whitening results guide and note the gradual improvement.
7. We do not recommend eating or drinking anything except plain water for up to 30 minutes after whitening.

Whitening intensity

Blue mode - 15 min

Strong whitening treatment. If your teeth are badly discoloured and you need to get rid of coffee or cigarette stains, use the 15-minute whitening treatment with whitening gel.

Blue-red mode - 10 min

Delicate and gentle whitening treatment. For sensitive teeth and daily cleaning. We recommend a 10-minute whitening treatment without applying the whitening gel. The device with the LED light switches to this mode by double tapping.

Storage

Keep out of direct sunlight, store in a cool and dry place.

Durability

12 months

Warning

Do not use:

- if you are under 18 years of age,
- if you have very sensitive teeth,
- if you have rotten, badly damaged or loose teeth,
- during pregnancy and when breastfeeding,
- if you are allergic to the ingredients listed on the packaging.

The whitening product is not effective on tooth replacements such as fillings, crowns, bridges. The product should only be used to whiten teeth. Keep out of reach of children. Do not swallow.

- Avoid applying the gel directly to your gums. Wipe it off as soon as possible.
- If the gel is on your face or clothes, you should wash it off as soon as possible.
- Once the whitening process is complete, put the cap on the pen to prevent it from drying out.
- Keep the product in a safe place out of the reach of children.
- Not suitable for ceramic and artificial teeth.
- Not suitable for damaged enamel and damaged teeth.
- Do not apply the gel directly to your gums, as this may cause a burning sensation. If this situation occurs, keep calm, the burning sensation will disappear within 24 hours. Just remember to rinse your mouth out with water.
- This item is a cosmetic product.
- If you are unsure about using this product, consult your dentist.

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