
medivon

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

User manual EN

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Medivon Lumi 3D Pro

Model: G11P

A
B
C

M1
M2
M3
M4
M5
M6

1
2
3
4
5

PL

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EN

This manual applies to the Medivon Lumi product. Before starting to use the device, read this user manual and keep it in a safe place.

Technical Data: Power supply: 5 V / 2 A | Rated power: ≤ 16 W Battery capacity: 5000 mAh | Number of LEDs: 735 | Wavelengths: 415/590/633/850/1072 nm (only the face mask supports the 1072 nm band)

1. Intended Use

The LED mask for the face, neck, and décolleté is designed for light therapy and is a medical device intended to support facial skincare through controlled exposure to light at specific wavelengths.

Effects of individual types of light:

Yellow light – evens out skin structure, stimulates circulation, reduces redness, swelling, and fine lines.

Red light – exhibits anti-aging effects, supports natural collagen production, and improves skin firmness.

Near-infrared – supports skin regeneration processes, improves microcirculation and hydration levels.

Blue light – supports the reduction of bacteria responsible for acne formation, reduces inflammation, and the visibility of acne scars.

2. Contraindications and Warnings

Do not use the product in case of pregnancy, planning pregnancy, or breastfeeding, light hypersensitivity (photosensitivity), photosensitive epilepsy, photosensitive diseases (e.g., lupus erythematosus), eye diseases, active skin diseases, damage or irritation at the application site, or by persons under 18 years of age.

Do not use the device while taking medications or using cosmetics that may cause light hypersensitivity.

During use, apply protective eye shields, avoid looking directly into the light source, and use the product solely in accordance with this manual. The product is not waterproof; do not disassemble or modify it, and keep it out of reach of children.

In case of any adverse reactions, such as: severe redness, burning, pain, itching, swelling, dizziness, worsening of well-being, or other atypical symptoms, immediately stop using the product.

If symptoms persist or worsen, contact a doctor or dermatologist immediately.

The mask is intended exclusively for external use, in accordance with this manual.

3. Operating Instructions and Method of Use

Step A – Preparing the LED Mask for the Face, Neck, and Décolleté
Remove the protective film from the inner surface of the LED mask for the face, neck, and décolleté. Make sure that the surface of the mask is clean and dry.
Step B – Applying Eye Protection
Put on the silicone protective eye covers as intended in order to protect your eyesight from light emission.
Step C – Connecting the Controller
Connect the cable of the LED mask for the face, neck, and décolleté to the controller, and then ensure that the device is fully charged or connected to a power source.

Device operation is performed using the remote control marked with numbers 1-5:
(1) Power Port – used for charging the device and connecting the LED mask for the face, neck, and décolleté to the remote. Before use, make sure the device is fully charged.
(2) Display – shows the remaining treatment time, the current operating mode, and the device status.
(3) Power / Mode Selection Button – press and hold for about 3 seconds to turn on the device. Subsequent short presses allow you to select one of the available operating modes (M1–M6).
(4) Intensity Adjustment Button – a short press allows you to set the light intensity level according to the user's preference (50%, 75%, 100%).
(5) Pulsation Mode Button – enables or disables pulsation mode in any operating mode.

After startup and selection of appropriate settings, the device works automatically for the default time of 10 minutes, after which it turns off automatically. During the treatment, it is recommended to remain in a sitting or lying position and avoid moving the mask.

After finishing the session, remove the mask and turn off the device. After every use, wipe the mask for face, neck and décolleté with a soft, dry, or slightly damp cloth. Do not immerse in water. Do not use aggressive cleaning agents.

4. Light Sensitivity Test

Before the first startup of the product, it is mandatory to perform a test by applying the switched-on device to the arm for 10 minutes. If you feel any discomfort during the treatment, stop using it immediately. The appearance of itching or redness after the test indicates skin hypersensitivity and excludes the use of the mask, whereas slight redness that disappears spontaneously within 1 to 2 hours is a completely normal reaction.

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EN | Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution. Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment. Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m^2 devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market after 13 August 2005.

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Model: G11P
Basic UDI-DI: 697600643LEDMaskZH

Date of preparation of the user manual: March 2026

CE Declaration of Conformity

Manufacturer:	Shenzhen Sungrow LED Technology Co., LTD 3F, 2# Factory Bldg(Phase I) Bangkai Sci & Tech Ind Park, Guangming HighTech Ind Park, Shenzhen, GuongDong Prov, China 518000
SRN #:	CN-MF-000047979
European Representative:	Riomavix Sociedad Limitada Calle de Almansa 55, 1D, Madrid 28039 Spain Tel.: +34 658 396 230 SRN #:ES-AR-000001202
Trade Mark:	SGROW
Product name:	LED Light Therapy Mask
Modle:	G10/ G11/ G11P/ G12/ G13/ G14/ G15/ G15P/ G15K/ G16/ G17/ G18/G19/ G20 / G21 / G22/G25/G56 /G240/ G160/G120 /G2030/ G3040/ G3060 /G1980 /EM01 /EM02/EM03/ EM04 /EP01/ EP02 /EP03/MS01/MS02/MS03 /LP01/LP02/LP03 /BE-G122/ BE-G80
Basic UDI-DI:	697600641LEDMaskZH
Intended Use:	LED Light Therapy Mask is an over the counter device that is intended to use LED light for the alleviate wrinkles and mild to moderate acne.

MDR Classification: Class I, According to Rule 1, Annex VIII of Regulation.
Conformity Assessment Route: EU DECLARATION OF CONFORMITY following the Annex II + Annex III + Article 19 of MDR (EU) 2017/745.
We herewith under our sole responsibility declare that the above mentioned products meet the transposition into national law, the provisions of the following EU Regulation and Standards. All supporting documentations are retained under the premises of the manufacturer.
The manufacturer is exclusively responsible for the declaration of conformity.General applicable regulations, directives:Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.

Applied standards, common specification, guidance:
EN ISO 15223-1:2021,
EN ISO 20417:2021,
EN ISO 14971:2019,
IEC 62366-1:2015+A 1:2020,
EN ISO 10993-1:2010,
EN ISO 10993-5:2009,
ISO 10993-10:2010,
EN ISO 13485:2016,

Signature: Cong Fang
Date: MAy. 22, 2025
Name:Cong Fang
Title:Manager
Position:Shenzhen , CN

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Wytwórca:	Shenzhen Sungrow LED Technology Co., LTD 3F 2# Factory Bldg (Phase I) Bangkai Sci & Tech Ind Park, Guangming HighTech Ind Park, Shenzhen, Guangdong Prov Chiny 518000
SRN #:	CN-MF-000047979
Autoryzowany przedstawiciel:	Riomavix Sociedad Limitada Calle de Almansa 55, 1D, Madrid 28039 Hiszpania Tel.: +34 658 396 230 SRN #:ES-AR-000001202
Znak towarowy:	**SGROW**
Nazwa wyrobu:	Maska do terapii światłem LED
Modele:	G1O/ G11/ G11P/ G12/ G13/ G14/ G15/ G15P/ G15K/ G16/ G17/ G18/ G19/ G20 / G21/ G22/ G25/ G56 / G240/ G160/ G120/ G2030/ G3040/ G3060/ G1980/ EMO1/ EMO2/ EMO3/ EMO4/ EPO1/ EPO2 / EPO3/ MS01/ MS02/ MS03/ LPO1/ LPO2/ LP0O3/ BE-G122 /BE-G80
Basic UDI-DI:	1697600641LEDMaskZH
Zamierzone zastosowanie:	Maska do terapii światłem LED to dostępne bez recepty urządzenie, które wykorzystuje światło LED w celu łagodzenia zmarszczek oraz trądziku o nasileniu od łagodnego do umiarkowanego.

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Shortened warranty terms
This Guarantee is provided by Armare S.A., a company governed by Polish law, with its registered office in Warsaw, ul. Grzybowska 87, 00-844, Poland. You can find us in the official Polish register of companies under the KRS number: 0000730590, NIP: 5213833614 and REGON: 380721860. We guarantee that the equipment you have purchased from us (hereinafter referred to as the "Product") will work properly. If something goes wrong for a reason inherent in the Product, this warranty protects you. The warranty starts from the day you receive the Product and lasts:
24 months, if you are a consumer or a natural person running a business and the purchase of the Product is not of a professional nature for you. Otherwise, the warranty period is 12 months.
If the law in your area requires a longer warranty, the warranty period will be extended to meet local law requirements.

In addition, if the Product does not comply with the contract, you have the right to legal protection from and at the expense of the seller. This warranty does not affect these rights - it does not limit, suspend or exclude them.

Detailed warranty conditions, how to submit complaints, what happens when you find a defect, what are exceptions to the warranty, when the warranty expires or is extended, and what are the rules regarding the processing of personal data, you can find on our website:
http://medivon.pl/warranty

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Tłumaczenie sporządzone przez importera Armare S.A., ul. Grzybowska 87, 00-844 Warszawa, Polska