PL | EE | EN

Momini BubbleBreath

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HUNAN ACCURATE BIO-MEDICAL TECHNOLOGY CO.,LTD. 湖南准生物医疗科技有限公司 430105013731G

1 Instructions for Use

Thank you for purchasing this nasal irrigator manufactured by Hunan Accurate Bio-Medical Technology Co., Ltd. ("Accurate"). To ensure proper use of the nasal irrigator, please read this manual carefully and fully understand it before first use. After reading, keep it in an easily accessible place for future reference.

Intellectual Property

The Importer owns the copyright to this unpublished user manual and may treat it as confidential information. This manual is intended as a reference for the operation, maintenance, and repair of this product. No one may use, disclose, or otherwise obtain the information contained in this manual, in whole or in part, in any way without the Importer's written consent.

The Importer reserves the right to final interpretation of this manual.
The Importer reserves the right to change the contents of this manual without prior notice.
The Importer reserves the right to make technical changes without prior notice.
The Importer reserves the right to change the product specifications without prior notice.

Statement

The Importer disclaims any warranty of any kind with respect to this manual, including (but not limited to) the implied warranties of merchantability and fitness for a particular purpose.

The Importer is responsible for the safety, reliability, and operation of the device only under the following conditions:
• The device is repaired by persons authorized by the Importer;
• The power supply used meets the requirements specified in the manual;
• The user operates the device in accordance with the manual.

The Importer is not responsible for the safety, reliability, or operating condition of the device in the following situations:
• The device is repaired or modified by persons not authorized by the Importer;
• The device is not used properly in accordance with the User Manual.

Repair Service

The customer is entitled to a one-year free repair service if the device complies with the Importer's warranty service regulations.

Scope of Paid Services

If the device does not comply with the Importer's warranty service regulations, the Importer will provide paid services. Even during the warranty period, repair of the device is subject to a fee if the cause is: damage caused by the user; Force majeure events, such as natural disasters. The Importer is not liable for any direct, indirect, or consequential losses, damages, or delays caused by the following situations (including, but not limited to): improper use; replacement with parts not approved by the Importer; or repair of the device by persons not authorized by the Importer.

Return

If returning the product to the Importer is absolutely necessary, please follow the steps below: Obtain a Return Merchandise Authorization (RMA). Contact the Importer's Customer Service Department or your local reseller and provide the product's serial number (located on the device). Provide the device model and serial number and briefly describe the reason for the return. Shipping costs: The cost of shipping the device to the Importer for repair is the user's responsibility.

The following circumstances are not covered by the warranty:

Failure due to failure to use the product in accordance with the instructions.
Failure due to impact or dropping the device due to carelessness during use or handling.
Failure due to lack of proper maintenance (see "Precautions during use") or operation in an unsuitable environment.
Failure caused by improper connection to other equipment.
Failure caused by disassembly, modification, or repair without the Importer's consent.
Failure caused by the use of parts that do not meet the Importer's technical requirements without its consent. Failure caused by accidents.
Transportation Costs: Transportation costs associated with repairing or replacing the product.
The product's wiring diagram, component list, and other materials necessary for repair or maintenance are provided only to technical service personnel trained and approved by the Importer.
If warranty or repair service is required, contact your supplier directly by phone, telegram, mail, fax, or other means of communication.

Device Shipping

If it is necessary to return the device to the manufacturer or retailer to prevent damage during transport, the device must be properly packaged—if possible in the original packaging—and the services of a qualified transport company must be used. If the original packaging is lost, select appropriate packaging and take appropriate protective measures to prevent damage to the device during transport. For long-distance shipments, please pack the device with particular care. Accurate is not responsible for any damage to the device incurred during shipment.

Product Information

Product Name: Nasal Irrigator
Product Model: NS1A3, NS2A3, NS3A3, NS1B3, NS2B3, NS3B3, NS1C3, NS2C3, NS3C3, NS1D3, NS2D3, NS3D3, NS1E3, NS2E3, NS3E3
Manufacturer: Hunan Accurate Bio-Medical Technology Co., Ltd.
Address: Accurate Industrial Park, No.108, Zhixian Road, Xuelian Community, Xueshi Street of Yuelu District, 410208 Changsha, Hunan Province, PEOPLE'S REPUBLIC OF CHINA
Tel: +86-731-85598539

Authorized representative: Shanghai International Holding Corp. GmbH (Europe), Eiffestraße 80, 20537 Hamburg, Germany

Revision History

This manual has a revision number. This revision number changes whenever

the manual is updated due to software or technical specification change.
Contents of this manual are subject to change without prior notice.
- Revision number: 1.1
- Release time: 2024-04

© Copyright Hunan Accurate Bio-Medical Technology Co., Ltd.
All rights reserved.

1.1 Basic Information of Product
Product Name: Nasal Irrigator
Model/Specification: NS1A3, NS2A3, NS3A3, NS1B3, NS2B3, NS3B3, NS1C3, NS2C3, NS3C3, NS1D3, NS2D3, NS3D3, NS1E3, NS2E3, NS3E3
Manufacturing Date: See the label at the bottom of the product
Service Life: 3 years
Configuration and Components: The nasal irrigator consists of main unit, washing head, internal lithium battery, power cord and mobile terminal software. The washing head consists of spray nozzle, adapter, cover, suction tube, water tank, etc.

NOTE:
Please purchase certified nasal wash solution/nasal salt for use in conjunction with the product.

1.2 Technical Performance Indicators
1.2.1 Key Technical Indicators

① Spray rate
Model: NS1A3, NS2A3, NS3A3, NS1B3, NS2B3, NS3B3, NS1C3, NS2C3, NS3C3, NS1D3, NS2D3, NS3D3, NS1E3, NS2E3, NS3E3
When the main unit is set to full-rate spraying, i.e., 12mL/min (±3mL/min), the spray is generated continuously;
When the main unit is set to low-rate spraying, i.e., 10mL/min (±3mL/min), the spray is generated continuously;
When the main unit is set to low-frequency pulse, i.e., 7mL/min (±3mL/min), the pulse spraying cycle is 2.0s (error ± 25%), spraying 1.0s (error ± 25%), pause 1.0s (error ± 25%);
Spray Height >200mm.
The maximum capacity of water tank is 25mL (±5ml).
Residual liquid volume: 2ml.

Operating principle:
Based on the Venturi ejector principle, the compressed air produced by the pump inside the main unit forms high-speed gas flow through a small tube; the negative pressure generated brings the nasal wash solution to eject toward the hole. Under the impact of high-speed gas flow, the liquid flow becomes fine droplets and mist particles which come out of the spray nozzle and enter the patient's nasal cavity for rinsing.

Contraindications: The product is made of such materials as ABS, PC and medical silicone; it is forbidden to use the product for those allergic to any of these materials.

NOTE:
It is strictly forbidden to use the product in an environment with any explosive, flammable or corrosive gases.

WARNING
• Where the storage temperature is below 5°C, the device should be placed in an environment under normal operating temperature for over 4 h prior to use. If the device works beyond normal conditions, its service life will be shortened.
• Please read the manual carefully before use. The device should be used strictly according to the scope of application and operation procedure described in this manual. If you have any questions, please contact your supplier or the manufacturer.

1.2.2 Power Parameters
Specific power input: 5V (Voltage), 0.5A
Lithium battery is used in the main unit, which can be charged repeatedly for approximately two years. At the end of the battery life, please contact our After-sales Service Department for replacement.

1.2.3 Type of Nasal Irrigator
• Type of power supply: Internal electrical power source
• Degree of protection against electric shock: Type BF applied parts
• Ingress protection rating: The ingress protection rating of main unit is IP22
• Not for use in the presence of flammable gases, and not category AP or APG equipment
• Working mode: Continuous operation
• Without applied parts protected from the defibrillator discharge effect
• Without signal output part
• Non-permanently installed equipment

1.2.4 Working Environment
Ambient temperature: 5°C–40°C
Relative humidity: ≤90%
Barometric pressure: 70106kpa
Temperature of the spray nozzle when contacting human body: ≤48,3°C

1.2.5 Transportation and Storage Conditions
Ambient temperature: –20°C–55°C
Relative humidity: ≤90%
Barometric pressure: 70106kpa
1.3 Features and Scope of Application of Nasal Irrigator
1.3.1 Features of Nasal Irrigator

• Multiple spray modes;
• Three-level adjustment;
• User-friendly design, convenient, time-saving and comfortable operation;
• Facilitating use by different family members.

1.3.2 Scope of Application
It is intended for nasal rinsing and washing at home or in clinics, hospitals or other places.


1.3.3 Intended Users
For users who need nasal washing or therapy, the product can be used in conjunction with nasal wash solution to wash the nasal cavity (children should use the product under adult supervision).

1.4 Composition of the Nasal Irrigator
1.4.1 Nasal Irrigator
The nasal irrigator consists of main unit, washing head, internal lithium battery, power cord and mobile terminal software. The washing head consists of spray nozzle, adapter, cover, suction tube, water tank, etc.


Packing List: ① Main unit × 1; ② Washing head × 1; ③ User Manual × 1; Power cord × 1.

1.4.2 Information of Mobile Terminal Software
Version of the software launched: V1.0.0
Operating Environment Requirement
The device supports devices with IOS 10.0 or Android 6.0 and above, as well as devices with Bluetooth 4.0 and above.
Security Software
Installation of security software such as antivirus or firewall is not required.
Data and Device Interface
Bluetooth 4.0 and above.
User Access Control Mechanism
User identity authentication is done via an account and password.
Software Environment and Related Security Software Update Requirements
Users can manually update the device software from the app market; no security software updates are required.
Installation and Maintenance Updates
Users can manually download, install, and update the software from the App Store or other app markets.
Uninstallation
For software uninstallation procedures, please refer to the relevant user manual of your smart device.

1.4.3 Diagram of Nasal Irrigator
Model: NS1
Description of Detachable Parts and Consumables
• Removable components of the product include cover and washing head (including spray nozzle, adapter, suction tube and water tank).

Model: NS2, NS2+

• User must purchase detachable parts from the original manufacturer for replacement; otherwise, the effect of the product cannot be achieved.
• Washing heads are consumable, and it is suggested to use one washing head for one person only. For replacement, user must purchase washing heads from the original manufacturer; otherwise, bacterial infection could be caused.

1.5 How to Operate the Nasal Irrigator
1.5.1 Cleaning and Disinfection
Cleaning and disinfection should be carried out before initial use. The operation method is detailed as follows:
1.5.1.1 Cleaning
① If stain exists on the surface of main unit shell, or there is residual spray solution after each use, use a clean soft cloth dipped with 75% medical alcohol to wipe the surface gently after wringing out, and then wipe with medical gauze to remove any moisture.
Remove the washing head from the main unit; use purified water (the temperature of which should not be higher than 50 ☒) to wash for 3~5min. After washing, wipe with medical gauze to remove any moisture and ensure full drying.

NOTE:
① DO NOT wash the device body in water to avoid water ingress which will damage the main unit; take care to prevent water ingress into the main unit.
② DO NOT use any volatile liquid (e.g., benzene) other than medical alcohol to wipe the main unit.
1.5.1.2 Disinfection
The device should be thoroughly cleaned before disinfection. Disinfection method: Use 75% medical alcohol to wipe the surface

of parts of the irrigator, and ensure that all contaminated surfaces come into contact with the alcohol. The device can be powered on only after all surfaces become fully dry.

WARNING
The product shall not be sterilized under low/high temperature or by γ-rays.

1.5.2 Charging
For charging with the USB cable, connect one end of the USB cable to the charging port at the bottom of the main unit and the other end to the power adapter, and connect power supply to start charging. The device can work normally after charging.


NOTE:
To ensure safe use, please purchase a power adapter conforming to the medical electrical safety standard GB9706.1-2020 or IEC60601-1: 2005.
For models NS1A3, NS1B3, NS1C3, NS1D3 and NS1E3, in case of low battery during working, the battery indicator will flash continuously, indicating charging is necessary, and the nasal irrigator will stop working.
For models NS2A3, NS3A3, NS2B3, NS3B3, NS2C3, NS3C3, NS2D3, NS3D3, NS2E3 and NS3E3, in case of low battery during working, the battery indicator on the display screen will flash continuously, indicating charging is necessary, and the nasal irrigator will stop working.

1.5.3 Continuous Operation Time
The NS nasal irrigator fully charged can work continuously for at least 30min.

WARNING
Before using the device, first insert the rubber plug to the charging port to prevent failure of the main unit caused by liquid ingress.
Without permission, DO NOT attempt to remove or repair any component of the device or disassemble the battery. If replacement is needed, please contact Accurate's authorized service center or return the device to Accurate's factory for replacement with the original part or component.

1.5.4 Installation
• Installation Method

① As shown in the installation diagram, install the suction tube, cover and adapter (already assembled before delivery) in

sequence;
② Fix the main unit, and simply tighten the washing head clockwise.

NOTE:
① Prior to use, please ensure that the washing head is free of foreign matters;
② After spray solution is loaded, avoid inverting the device, or randomly opening the cover, or shaking the device violently, to avoid liquid overflow.

• Add Solution to the Water Tank
Solution (for reference: 10ml, 15ml, 25ml or other spray wash solution and seawater wash solution) should be added to the sink on the suction tube side;

NOTE:
When adding solution, confirm that the sink is the right one, and close the cover and adapter.
Normally, 10–25 ml of wash solution is needed for each nostril. You can adjust the dosage based on your own condition or according to your doctor's advice.
For models with display screen, after pressing the ON/OFF button, the screen will turn on; after pressing the ON/OFF button again, the device will start working.
For models without display screen, after pressing the ON/OFF button, the mode indicator will turn on; after pressing the ON/OFF button again, the device will start working.
Remove the washing head for timely cleaning, disinfection and maintenance.

Attention for people using the product for the first time!!
To prevent cross-infection, it is suggested to use one washing head for one person only.
For first use after purchasing, and before and after each use thereafter, each component of the washing head should be cleaned, disinfected, and air-dried for storage.
1.5.5 Description of Functions
1.5.5.1 Mode Selection

(for models NS1A3, NS1B3, NS1C3, NS1D3, NS1E3)
The default mode is low-rate spraying. Press the ON/OFF button to start or stop; press the MODE button to switch among full-rate spraying, low-rate spraying and low-frequency pulse. Mode indicator or working status is provided during spraying.

1.5.5.2 Mode Selection
(for models NS2A3, NS3A3, NS2B3, NS3B3, NS2C3, NS3C3, NS2D3, NS3D3, NS2E3, NS3E3)
The default mode is low-rate spraying. Press the ON/OFF button to start or pause; press the MODE button to switch among full-rate spraying, low-rate spraying and low-frequency pulse. During spraying, short press the ON/OFF button to paused or start working, or long press the ON/OFF button to exit operation. During spraying, the mode indication or working status is displayed on the screen.

1.5.5.3 Alarm Clock Setting
(for models NS2A3, NS3A3, NS2B3, NS3B3, NS2C3, NS3C3, NS2D3, NS3D3, NS2E3, NS3E3)
Alarm Clock Setting: In the default interface, long press the MODE button for 4s and release it to enter the setting state. Short press the ON/OFF button repeatedly, Hour, Minute and Alarm Clock (24H) will flash respectively. Press the MODE button to set the time. After setting, press the ON/OFF button; the device will automatically enter the Alarm Clock Setting interface. Press the SET button to respectively set the time and the working status of alarm clock ("ON/OFF"). Up to 3 alarm clocks can be set. After setting, press the ON/OFF button to go back to the default interface. In the setting process, you can also long press the ON/OFF button to exit setting.

1.5.5.4 Volume Setting
(for models NS2A3, NS3A3, NS2B3, NS3B3, NS2C3, NS3C3, NS2D3, NS3D3, NS2E3, NS3E3)
In the default interface, short press the SET button to increase the volume (by 1 level each time this button is pressed); there are totally 5 volume levels for adjustment. Repeatedly press the SET button to adjust the volume. You can turn up/down the volume or mute the sound.

1.5.5.5 Alarm Ring OFF Setting
(for models NS2A3, NS3A3, NS2B3, NS3B3, NS2C3, NS3C3, NS2D3, NS3D3, NS2E3, NS3E3)
When an alarm clock rings, press the ON/OFF button to stop the alarm.

1.5.6 Storage on the Device
Model: NS3A3, NS3B3, NS3C3, NS3D3 and NS3E3 have the function of storing the number of washing times.

1.5.7 Data Transmission
Model: NS3A3, NS3B3, NS3C3, NS3D3, NS3E3
After turning on the device and launching the app on your mobile phone, your mobile phone will automatically search Bluetooth devices. After searching out the device name "Nasal Irrigator", click this name to connect. After the screen of the device displays the Bluetooth connection icon, data can be transmitted to the mobile terminal software.

1.5.8 Nasal Washing Operation
• Align the spray nozzle of washing head with a nostril.
• Gently press the button on the main unit to turn on the device, and use the device to start washing according to the instructions.

NOTE:
During washing, moderately bend and lean forward your head to achieve a better washing effect.

1.6 Contraindications
(1)Patients after inner ear or eardrum surgery, or whose eardrum is injured, or whose mouth and pharynx are injured or have just gone through a surgery, or who cannot cooperate due to any nerve disease, or with dysgnosia, epilepsy, blood disease, or serious infectious disease.
(2)For neonates and pregnant women, the product should be used with caution.
Precautions and Reminders
(1) To ensure correct use of the product, please use it under your doctor's guidance, and carefully read the User Manual before use and keep the manual properly. If the manual is lost, please contact our After-sales Service Department.
(2) To avoid the risk of apnea, please keep the product out of the reach of infants and children. In case of any accident, please seek medical advice promptly.
(3) To prevent any component from being sucked in or swallowed, please keep the product out of the reach of infants and children. In case of any accident, please seek medical advice promptly.
(4) Without permission, DO NOT attempt to disassemble the device or use accessories other than the specified ones; otherwise, damage to the device could be caused.
(5) DO NOT connect the device to other equipment; otherwise, unpredictable damage could be caused to the device. If such damage occurs, please contact our After-sales Service Department.
(6) DO NOT attempt to modify the device; otherwise, unpredictable damage could be caused to the device. If such damage occurs, please contact our After-sales Service Department.
(7) It is forbidden to use the device in an environment with inflammable anesthetic gas mixed with air, oxygen or nitrous oxide; otherwise, the explosion hazard could be caused. If such hazard occurs, please contact your local fire department.
(8) DO NOT use the device in an environment with strong electromagnetic interference; otherwise, the device may fail. If such failure occurs, please contact our After-sales Service Department.
(9) The washing head can be used for single person only; otherwise, cross-infection could be caused. If cross-infection occurs, please seek medical advice promptly.
(10) After use, please timely wash the washing head. Remember that the main unit should not be washed in water; otherwise, unpredictable damage could be caused to the device.
(11) For nasal wash solution, 0.9% normal saline is recommended, or follow your doctor's advice.

WARNING:
Always pay attention to the pressure in your nasal cavity and middle ears. For some patients with nasal obstruction, water flow

from the nasal irrigator could cause discomfort due to excessive pressure on nasal cavity or middle ears. In such case, irrigate the nostrils one by one, that is, letting water flow back from the nostril without blocking the nostril with the spray nozzle. If serious discomfort occurs to your middle ears during or after use of the product, please stop using it immediately and seek medical advice. DO NOT blow your nose within 5 min after using the product but use tissue paper to gently squeeze away water in the nasal cavity and let it flow out; otherwise, the pressure of middle ears will be increased, causing discomfort. The product should not be used in one nostril for over 5 min each time.

1.6.1 Use of the Product by Children

Children and those special people must use the product under correct guidance and supervision. The nasal irrigator cannot be used for any purposes other than nasal washing. Each time the device is used, or for storage of the device, components of the washing head must be cleaned and disinfected. If it is found that any part is damaged or drops into water due to carelessness, DO NOT continue using the device.

1.6.2 Use of the Product by Surgical Patients

For patients with organic lesion in nasal cavity, first receive treatment and recover from the organic disease, get the certificate allowing use of the product issued by the local medical organization and obtain doctor's guidance before using the nasal irrigator. This will help the patients to greatly relieve their symptoms and control their conditions. When using the product, surgical patients must be guided by a doctor or seek medical advice.

Take and add drugs: If you are taking drugs for treatment of rhinitis, please use the nasal irrigator under the guidance of a doctor. You are not suggested to directly add drugs to salt solution; if it is needed to do so, doctor's guidance is required. DO NOT irrigate your nasal cavity using salt solution added with drugs for a long period of time.

1.6.3 Operations after Use

After each use, timely clean, disinfect and thoroughly air-dry the nasal irrigator; then install the washing head back to the main unit, and put the main unit in the packaging box along with other parts. Store the device in an environment meeting the requirements, and pay attention to avoiding impact.

2 General Information

2.1 Precautions for Safe Use of the Nasal Irrigator

• This product is a rechargeable device; for any questions, please consult our after-sales personnel.
• To minimize injury, close observation is required when children use the electric irrigator;
• DO NOT unplug the charging cable by wet hand;
• DO NOT use the product during bathing;
• DO NOT place or store the product in an area where it could fall or be pulled into a bathtub or shower cubicle;
No matter whether it is the first time that the product is used, always check whether the charging cable is damaged as long as it is within the service life of the product;
• DO NOT use the product when the charging cable is damaged or when the product does not work normally, gets damaged due to drop or falls into any liquid. Please timely contact our company or authorized service center;
• DO NOT insert any foreign object into the device hole or power socket;
• DO NOT use the product in an environment where O2 or aerosol is used;
• DO NOT wear any nose ornament when using the product;
• DO NOT place the product close to any hot gas, electric heater or heating oven to avoid damage to the product due to high temperature;
• DO NOT use parts (e.g., spray nozzle) not supplied by our company;
• DO NOT use the product for other purposes.
• DO NOT directly inject water into your eyes, ears or other sensitive parts since improper use could result in unexpected injury;
• To avoid electromagnetic interference, the place where the nasal irrigator is used should be absent of medical high-frequency equipment or other high-frequency RF equipment;
• The product should be used in a clean, well-ventilated, heat-insulated, damp-proof and dust-proof environment without toxic gases and corrosives;
• Upon expiry of the service life of the product, operation abnormality or other potential safety hazards could occur. You are not suggested to continue using the product beyond its service life. Please contact your local dealer or our company to purchase a new one;
• The nasal irrigator and its replacement parts after expiry of their respective service life should be disposed of according to applicable national measures for product disposal, thus to achieve the goal of environmental protection;
• To avoid infection, one spray nozzle should be used for one person only;
• The internal battery of the device must be replaced by our authorized maintenance personnel; user shall not attempt to replace it. DO NOT disassemble the battery, throw it into fire, or short-circuit it. Combustion, explosion or leakage of the battery could cause injury;
• If the device is stored under the minimum or maximum temperature, before use of the device, please place it under normal temperature for at least 4h.

WARNING:

If the device cannot work normally due to expiry of service life, or when the battery is replaced and removed from the device, please properly dispose of the device/battery in accordance with applicable national requirements, and scrap it if the conditions for scrapping of medical devices are satisfied. To avoid explosion, DO NOT throw the battery into the fire.

2.2 Maintenance and Cleaning of the Nasal Irrigator

2.2.1 Routine Maintenance and Cleaning after Each Use

NOTE

① After use, timely carry out cleaning and disinfection according to the above method to avoid residual spray solution which will affect next use.
② To avoid deformation caused by heating, the spray nozzle, adapter, cover, suction tube, water tank and other components of the washing head should not be washed by boiling or in boiled water.
③ During cleaning and disinfection, use clean and sterile medical gauze and medical cotton swabs. Residual disinfectant on any component should be washed clean to guarantee the safety during next use.
④ DO NOT use any sharp hard object to wipe the surface.
⑤ When the washing head is broken or aged, it must be replaced. If the device is used continuously, it is suggested to replace the washing head at the interval of about 3 months, depending on the actual condition of use. You are suggested to

directly contact our company to purchase components of the product.
⑥ The product does not need special maintenance during its service life, and it should be charged at least once a month; otherwise, the product could be damaged due to electrolyte leakage.

2.2.2 Preventive Checks
No matter whether the product is used frequently, it should be checked at least once a month;
Whether the surface of the product is damaged;
Whether the charging cable is damaged;
Whether the main unit, washing head, cover and suction tube are broken or obstructed;
Whether any detachable part is missing and whether the socket can be inserted smoothly.

2.3 EMC
Medical electrical equipment requires specialized instructions on electromagnetic compatibility, as well as installation and use based on the electromagnetic compatibility information provided in the accompanying documents. Portable and mobile RF communication equipment may affect the use of the equipment.
WARNING
Except for those sold as spare parts for internal components by manufacturers of equipment or systems, the use of accessories or cables outside of regulations may result in an increase in equipment or system emissions or a decrease in immunity.
WARNING
The device or system should not be used in close proximity or stacked with other devices. If it is necessary to use in close proximity or stacked, it should be observed and verified that it can operate normally in the configuration in which it is used.

The device meets the requirements of IEC 60601-1-2.

Guidance and manufacturer´s declaration - Electromagnetic Emissions		
The device is suitable for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emission tests	Compliance	Electromagnetic environment - guidance
Radio frequency (RF) emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Radio frequency (RF) emissions CISPR 11	Class A	The device is suitable for use in all establishments other than domestic and those indirectly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer´s declaration - Electromagnetic Immunity

The device is suitable for use in the electromagnetic environment specified below.
The customer or the user of the device should assure that it is used in such an environment.

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V/m 150 kHz to 80 MHz	3 V/m, 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the Pulse Oximeter, including cables, than the recommended separation distance calculated from the
Radiated RF IEC 61000-4-3	10 V/m, 80 MHz to 2,7 GHz, 80% AM at 1K Hz	10 V/m	equation applicable to the frequency of the transmitter. Recommended separation distance: $d = \left[\frac{3,5}{V_1}\right]\sqrt{P}$ $d = \left[\frac{3,5}{E_1}\right]\sqrt{P}$ dla 80 MHz–800 MHz $d = \left[\frac{7}{E_1}\right]\sqrt{P}$ dla 800 MHz–2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,[1] should be less than the compliance level in each frequency range.[2] Interference may occur in the vicinity of equipment marked with the following symbol:

Note 1: At 80 MHz to 800 MHz, the higher frequency range applies.
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer´s declaration – electromagnetic immunity

The device is suitable for use in the electromagnetic environment specified below.
The customer or the user of the device should assure that it is used in such an environment.

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge IEC 61000-4-2	±8 kV contact, ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact; ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst (EFT) IEC 61000-4-4	±2 kV power supply lines; ±1 kV signal input/output, 100 kHz repetition frequency	±1 kV signal input/ output, 100 kHz repetition frequency	N/A
Surge – IEC 61000-4-5	±0.5 kV, ±1 kV differential mode; ±0.5 kV, ±1 kV, ±2 kV common mode	N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	N/A	N/A	N/A
Power frequency (50/60 HZ) magnetic field IEC 61000-4-8	30 A/m, 50/60 Hz	30 A/m, 50/60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT is the a.c. mains voltage prior to application of the test level.

Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the Device

The device is suitable for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communication equipment.

Rated Maximum Output Power of Transmitter (W)	150 kHz–80 MHz $d=[\frac{3.5}{V_1}]\sqrt{P}$	80 MHz–800 MHz $d=[\frac{3.5}{E_1}]\sqrt{P}$	800 MHz–2,5 GHz $d=[\frac{7}{E_1}]\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz to 800 MHz, the separation distance for the higher frequency range applies.
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

2.4 Troubleshooting

The following simple checks should be carried out before requesting repair service:

Common Problems of the Device and Solutions

No.	Problem	Possible Cause	Rozwiązanie
1	The device cannot work normally after it is turned on	Dead battery	Charge the device
		No spray solution is loaded, or spray solution is inadequate	Check the condition of spray solution in the washing head
		The ON/OFF button gets stuck ON/OFF	Check whether the button is stuck
2	The spray rate is too low	The spray solution is inadequate	Add spray solution
		The device is not properly assembled	Reassemble and restart the device
3	The device shuts down automatically during use	Dead battery	Charge the device
		The device vibrates or shakes violently during use	Hold the device steadily
4	Spray solution is used up, but the device cannot shut down automatically	Some spray solutions could produce foam which will be left in the washing head	Press the ON/OFF button, and clean the foam in the washing head
5	Spray solution overflows from the washing head	Breakage or plastic aging	Replace the washing head

If the problem cannot be solved, please contact our After-sales Service Department.

2.5 Description of Symbols and Marks Used on the Device

Symbol	Describe	Symbol	Describe
	The device belongs to Type BF equipment.		Notice! Refer to accompanying documents.
	The transport package cannot be directly exposed to sunlight.		This side up during transport or placement.
	The packaging box contains fragile items. Handle with care.		Protect the package from rain.

Symbol	Describe	Symbol	Describe
	The device belongs to Class II equipment.	IP22	The main unit is designed with Level 2 protection against the ingress of water drops.
	DO NOT throw it away		Manufacturer
	ON/OFF button		Observe the instructions for use
— max —	Maximum water level		

Declaration of Conformity		
Legal Manufacturer:	Name: Hunan Accurate Bio-Medical Technology Co., Ltd. Address: Accurate Industrial Park, No.108, Zhixian Road, Xuelian Community, Xueshi Street of Yuelu District, 410208 Changsha, Hunan Province, People's Republic of China	
Authorized Representative:	Name: Shanghai International Holding Corp. GmbH (Europe) Address: Eiffestraße 80, 20537 Hamburg, Germany SRN: DE-AR-000000001	
Product name:	Nasal irrigator	
Trade name:	//	
Model:	NS1A3, NS2A3, NS1B3, NS2B3, NS1C3, NS2C3, NS1D3, NS2D3, NS1E3, NS2E3	
EMDN code:	Q030103: Nasal irrigation solutions	
Basic UDI-DI:	69273945NI02UY	
Intended purpose:	Intended for rinsing and washing the nose at home or in clinics, hospitals, or other locations.	
Risk class:	Class I	
Classification rule:	Rule 13 in Chapter III of Annex VIII to Regulation (EU) 2017/745	
Classification rule:	Rule 13 in Chapter III of Annex VIII to Regulation (EU) 2017/745	
We hereby declare that the above-mentioned products comply with Regulation (EU) 2017/745 of the European Parliament and of the Council. All supporting documentation is kept at the manufacturer's premises. We, as the manufacturer, bear sole responsibility for this Declaration of Conformity.		
Applied standards:	MDR (EU) 2017/745	ISO 10993-5:2009
	ISO 13485:2016	ISO 10993-10:2010
	ISO 14971:2019	ISO 10993-23:2021
	IEC60601-1-2:2014+A1:2020	Directive 2011/65/EU
	IEC60601-1-11:2015.AMD1:2020	IEC60601-1-6:2010+A1:2013+A2:2020
	ISO 10993-1:2020	MEDDEV.2.7.1 Rev4
	IEC 62304-2006+A1:2015	IEC 62366-1:2015+A1:2020
	ISO 15223-1:2021	AAMI TIR 12-2020
	AAMI TIR30:2011 (R2016)	
Place, date of issue:	Hunan, Changsha 2024-12-03	
Signature:	Name: ZHANG LI Position: DYREKTOR GENERALNY	

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Translation prepared by the importer Beauty Factor Sp. z o.o. (Ltd.)