Büg

EN Insect bite healer

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
AOJ HEALTH Technology Co., Ltd.

Room 302, Block A, Building A, Jingfa Intelligent Manufacturing Park,
Xiaweiyuan, Gushu Community, Xixiang Street, Bao'an District
Shenzhen 518126 CHINA

	Caution - Prudence - Precaución - Attenzione - Vorsicht - Voorzichtig - cautela
	Refer to instructions manual - Se référer au manuel d'instructions - Consulte el manual de instrucciones - Fare riferimento al manuale di istruzioni - In der Bedienungsanleitung nachlesen- Raadpleeg de handleiding - Consulte o manual de instruções
	Type BF Applied part - Type BF Partie appliquée - Parte aplicada tipo BF - Parte applicata di tipo BF - Angewandter Teil Typ BF - Toegepaste deel type BF - Tipo BF Peça aplicada
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive - Eliminación de acuerdo con la Directiva CE de Residuos de Aparatos Eléctricos y Electrónicos- Smaltimento conforme alla Direttiva CE sui rifiuti di apparecchiature elettriche ed elettroniche- Entsorgung gemäß der Richtlinie über Elektro- und Elektronik-Altgeräte (WEEE) - Verwijdering volgens de EG-richtlijn betreffende afgedankte elektrische en elektronische apparatuur - Eliminação de acordo com a Diretiva CE relativa aos resíduos de equipamentos eléctricos e electrónicos
	Manufacturer - Fabricant - Fabricante - Produttore - Hersteller - Fabrikant - fabricante
	Date of Manufacture - Date de fabricantion - Fecha de fabricación - Data di fabbricazione - Herstellungsdatum - Fabricagedatum - data de fabrico
	Batch code - Code de lot - Código de lote - Codice lotto - Chargencode - Partijcode - código de lote
	MR Unsafe - Non compatible IRM - Inseguro para RM - Non compatibile con RM - Nicht geeignet für MRT - Niet compatibel met RM - MR Inseguro

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WORKLING PRINCIPLE

The principle of Insect bite healer is to use heat energy to act on the surface of the skin, so as to reduce or eliminate the itching and discomfort caused by insect bites. At the same time, the heat effect can also promote local blood circulation, accelerate the metabolism and excretion of toxins, reduce itching and discomfort.

WORKLING PRINCIPLE

- Do not use the device if you suffer from diabetes. Sensitivity to pain could, in certain circumstances, be reduced and diabetics may burn themselves.
- Do not use the device on people with disabilities, children under 3 years old or people with reduced sensitivity to heat. This includes people with skin alterations due to illness or scarred tissue in the application area, after taking pain relief medication, drinking alcohol or taking drugs.
- Do not use the device on acutely or chronically diseased (injured or irritated) skin (e.g. inflamed skin, whether painful or not, reddened skin, rashes, (e.g. allergies), burns, frostbite, bruises, swellings (without signs of an insect bite/sting), both open and healing wounds, and surgical scars where the healing process could be affected).
- Do not use the device if you are suffering from a fever.
- Do not use the device on sensitive areas of skin.

• Do not use the device if you have any sensory impairment that reduces the feeling of pain (e.g. metabolic disorders).

• Do not use the device if you are using lotions, creams and gels locally at the same time.

• Do not use the device if you suffer from persistently irritated skin due to long-term heat application on the same area of skin.

GENERAL WARNINGS

• Over 3 years of age under supervision, can be used for self-treatment by anyone over 12 years of age.

• Before using the device, consult your doctor if you have any skin conditions, serious illness (in particular if you have a propensity to thromboembolic conditions or recurrent malignant growths), unexplained chronic pain in any part of the body.

•The device is not a substitute for medical consultation and treatment. If the symptoms persist despite using the Insect bite healer and stop using it. Seek medical advice immediately if the symptoms worsen whilst or after using the Insect bite healer and stop using it.

• Seek medical advice immediately if it is a tick bite. The bite may result in transmission of pathogens such as tick-borne encephalitis (TBE) or borreliosis.

• Seek medical advice immediately if the first signs of increased itching, skin redness, swelling, dizziness, nausea, shortness of breath, formation of wheals or a drop or increase in blood pressure occur.

• Do not use near to the eyes, on the eyelids, mouth or on mucous membranes. Furthermore, the device must not be used on moist skin.

• Do not allow children to play with the device.

• Store the device out of the reach of children and pets.

• Keep packaging material away from children (risk of suffocation).

• Please note the general hygiene measures when using the insect bite healer.

• If the device does not work properly, or if you feel unwell or experience pain, stop using it immediately.

•If users or patients have occurred any serious incident that relation to the device, please report to manufacturer and the competent authority of the Member State in which you are established.

INSERTING THE BATTERIES

1. Push open the battery compartment lid (5) .

2. Insert the batteries included in delivery into the battery compartment as shown inside the battery compartment. When inserting the batteries, ensure that the correct polarity (+/-) is observed.

3. Slide the battery compartment lid (5) closed until you feel and hear it click into place. Install the battery as above.

Note : Remove the batteries if you do not intend to use your Büg for a long time (over 3 months.)

USAGE

Warning : Before use, make sure that there is no insect sting remaining in the area that was stung. Remove the sting carefully before use. Use immediately after the insect bite or sting achieves the best results. If you wait too long to treat the area that has been stung or bitten and itching and swelling have already started, the Büg can only have a limited effect. Nevertheless, treat the area with your Büg. Usually, itching can be relieved, and healing can besped up.

1. Press the power button (2) to start your Büg. The device status LED (1) lights up green, a buzzer sounds and indicates to you that thedevice is ready to use.

2. Mode 1: This short application time is recommended for first use or use on sensitive skin, on a young child for example.

After being stung or bitten, apply the heating plate (4) of the Büg to the affected area and press the activation button (3) once. The device status LED (1) will illuminate blue, indicating that the application has started. The Büg maintains a treatment temperature of 47°C ±2°C (116.6°F ±3.6°F) for 3 seconds.

3. Mode 2: This application time is recommended for regular use on normal skin, on an adult for example. After being stung or bitten, apply the heating plate (4) of the Büg to the affected area and press the activation button (3) twice in quick succession. The device status LED (1) will illuminate yellow, indicating that the application has started. The Büg maintains a treatment temperature of 50°C ±2°C (122°F ±3.6°F) for 6 seconds.

4. After the end of the application, a buzzer sounds, heating of the heating. Plate (4) stops automatically and the device status LED (1) lights up green again. Remove the heating plate (4) from the area of skin once the signal sounds.

Important : Depending on the sensitivity of the person using the Büg and of the area of skin, the treatment temperature may feel unpleasant and cause slight skin redness and, in rare cases, irritated skin. To achieve this itch-relieving effect on the skin nerves, a short period of heat pain is necessary. Only mild heat pain is required for treatment success. If the heat function feels too hot, stop treatment immediately.

5. If required, you can use the heat function of the Büg on the sting/bite again after a break of 2 minutes or also use it immediately on another sting/bite in a different area. The

maximum number of 5 applications per hour on the same treatment area must not be exceeded.

6. If the device is no longer required, press the power button (2) to turn off the Büg. The green device status LED (1) turns off. If there is no further operation within 5 minutes, the Büg will automatically shut down.

7. Before storing, make sure that the Büg is switched off and the heating plate (4) has cooled down.

Note : Before the battery level becomes too low for other applications, the device status LED starts to flash 3x and an acoustic signal sounds 3x. Replace both batteries soon.

SAFETY NOTES FOR BATTERY

• If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.

• Choking hazard! Small children may swallow and choke on batteries. Therefore, store batteries out of the reach of small children.

• Risk of explosion! Do not throw batteries into a fire.

• If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.

• Do not disassemble, open or crush the batteries.

• Observe the plus (+) and minus (-) polarity signs.

- Protect the batteries from excessive heat.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a relatively long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries! Warning on electromagnetic compatibility
- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated above, this device as well as the other devices must be monitored to ensure they are working properly.
- Failure to comply with the above can impair the performance of the device.

CLEANING AND CARE

• The device is intended for multiple uses.

• It is recommended that the device is cleaned before first use and each additional use.

• The device must be switched off, allowed to cool down and the batteries removed each time before cleaning.

• Use only soft, dry cloths for cleaning the housing.

• Never use any substances for cleaning that could potentially be toxic if they came into contact with the skin or mucous membranes, or if they were swallowed or inhaled. Do not use any abrasive cleaning products and never submerge the device in water.

• Do not clean the device or accessories in the dishwasher!

• Avoid mechanical cleaning with brushes or similar, as this could result in irreparable damage.

• When using the device on multiple areas that have been stung and bitten, it is recommended that you clean the heating plate each time. This reduces the possibility of germs and bacteria being transferred from one area of skin to the next. Use a commercially available alcohol disinfection wipe to disinfect the heating plate.

This device is manufactured and registered in compliance with FDA regulations.

PROBLEM	CAUSE	SOLUTION
The Büg does not switch on.	The battery installation is incorrect	Install the battery over again.
	The batteries are flat	Replace both batteries as described in these instructions for use.
The Büg no longer gets hot enough.	The battery level is too low.	Replace both batteries as described in these instructions for use.
The Büg outputs 3 acoustic signals/the device status LED (1) flashes – even after switching the device OFF and ON again.	Battery level is too low, meaning that the Büg may no longer be able to be operated.	Replace both batteries as described in these instructions for use.
Application cannot be started, even though the device status LED (1) lights up green.	The Büg needs to be restarted.	Press the power button (2) again.
Device outputs 5 acoustic signals/the device status LED (1) flashes – even after switching the device OFF and ON again.	The Büg may be faulty.	Please contact the Bblüv Customer Services (info@bbluvgroup.com)

GENERAL PRECAUTIONS

• Do not touch the device with wet hands. Do not allow water to spray onto the device. Only operate the device if it is completely dry.

• Do not use any detergents or solvents to clean and maintain the device. If liquid penetrates the device this could cause damage to the electronics or other device parts and lead to a malfunction.

• The device should only be used by one person; use by multiple people is not recommended.

• Protect the device against direct sunlight and store in a dry place.

• Protect the device from impacts, moisture, dust, chemicals, marked temperature fluctuations and nearby sources of heat (ovens, heaters).

• Ensure that no water enters the device. If this should happen, remove the batteries and only use the device again when it is completely dry.

• If you store the device, no heavy objects should be placed on top of it.

• Before use, ensure that there is no visible damage to the device or accessories.In case of any doubt, do not use the

device, and contact your retailer or the specified Customer Service address.

• Store the device in a location protected against climatic influences. The device must be stored in the environmental conditions specified.

• When ambient temperature is less than -10 or higher than 50 .The time from minimum storage temperature to use is at least 4 hour; The time from maximum storage temperature is at least 4 hour.

DISPOSAL

This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Remove the batteries before disposing of the device/unit. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop.

Consult your local authority or your supplier for information about disposal.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

—Reorient or relocate the receiving antenna.

—Increase the separation between the equipment and receiver.

—Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

—Consult the dealer or an experienced radio/TV technician for help.

TECHNICALS SPECIFICATIONS	
Model	BÜG / IBR30
Size	109.5mm x 34mm x 24.2mm
Weight	25g
Power supply	2 x batteries, type AAA(LR03 1.5V)
Max. temperature	50°C ±2°C (122°F ±3.6°F)
Functions	sensitive application: 3s,47°C±2°C(116.6°F±3.6°F)
Regular application	6s,50°C±2°C (122°F±3.6°F)
Expected operating duration	Approx. 300 applications (per batteries)

PERMISSIBLE STORAGE AND TRANSPORT ENVIRONMENT	
Temperature	-10°C~+50°C
Relative humidity	30% to 85%
Pressure	700hPa to 1060hPa

PERMISSIBLE OPERATING ENVIRONMENT	
Temperature	+10°C~+40°C
Relative humidity	30% to 80%
Pressure	700hPa to1060hPa

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